

Pro Lab
User Manual



Pro Lab User Manual

v 1.181

November 2021

All rights reserved.

Copyright © Tobii Pro AB (publ)

The information contained in this document is proprietary to Tobii Pro AB. Any reproduction in part or whole without prior written authorization by Tobii Pro AB is prohibited.

Products that are referred to in this document may be either trademarks and/or registered trademarks of the respective owners. The publisher and the author make no claim to these trademarks.

While every precaution has been taken in the preparation of this document, the publisher and the author assume no responsibility for errors or omissions, or for damages resulting from the use of information contained in this document or from the use of programs and source code that may accompany it. In no event shall the publisher and the author be liable for any loss of profit or any other commercial damage caused or alleged to have been caused directly or indirectly by this document.

Content subject to change without notice.

Please visit [Tobii Pro Connect](#) for updated versions of this document.

Table of Contents

1 General Information	8
1.1 Welcome to Tobii Pro Lab	8
1.2 Symbols used in this document	8
1.3 How to collect diagnostic data for Pro Lab	8
1.4 How to cite Tobii Pro Lab	8
1.4.1 APA example for Tobii Pro Lab (software)	9
1.4.2 APA example for Tobii Pro (User Manual)	9
1.4.3 APA's own tex files (6th edition)	9
1.4.4 BIBTEX entries for Pro Lab (software and User Manual)	9
2 Introduction	11
2.1 Pro Lab editions	11
2.2 Project types	11
2.3 How Pro Lab is structured	12
2.3.1 Project Overview	12
2.3.2 Design	12
2.3.3 Record	13
2.3.4 Analyze	14
2.4 Pro Lab licenses	15
2.5 System requirements	15
3 Install or update Pro Lab	16
3.1 Windows versions	16
3.2 Pro Lab license key	16
3.2.1 Activate or deactivate a Pro Lab license	16
3.3 Install Pro Lab	17
3.4 Upgrade contracts	17
3.5 Update Pro Lab	18
4 Start a project	19
4.1 Create a new project	20
4.2 Project Overview	20
5 Manage projects	22
5.1 Open an existing project	22
5.2 Rename or remove a project	22
5.3 Export and import a project	22
5.3.1 Distributed data collection	23
5.3.2 Distributed manual coding	25
5.3.3 Back up projects	25
5.3.4 Import data from SD card	26
5.4 Manage recordings	26
5.4.1 Replay	27
5.4.2 Delete	27
5.4.3 Import and Export recordings	27
5.5 Manage Event Types	27
5.5.1 Create an Event Type	28
5.5.2 Edit an Event Type	28
5.5.3 Delete an Event Type	28

5.6 Analysis tools	28
6 Design an experiment	30
6.1 Timelines	30
6.1.1 Create a new timeline	30
6.1.2 Delete a timeline	32
6.1.3 Duplicate a timeline	32
6.1.4 Change the name of a Timeline	32
6.2 Stimulus properties	32
6.2.1 Add a stimulus to the timeline	32
6.2.2 Advance on	33
6.2.3 Calibration stimulus	34
6.2.4 Image stimulus	36
6.2.5 Video stimulus	37
6.2.6 Stimulus editor (images and video)	38
6.2.7 Web stimulus	38
6.2.8 Screen Recording stimulus	39
6.2.9 Text stimulus (Reading studies and instruction)	39
6.3 Groups	41
6.3.1 Edit a group's contents	41
6.3.2 Duplicate a group	42
6.3.3 Delete a group	42
6.4 Stimulus/group variables	42
6.4.1 Create variables	42
6.4.2 Assign variables	43
6.5 Use TTL markers to indicate stimulus onset for external listeners	43
6.6 Change the appearance of the Timeline stimulus cards	44
6.6.1 Change the size of the stimulus cards	44
6.6.2 Emulate the screen resolution of the presentation screen	44
6.7 Correct resolution and scaling of the presentation screen	45
6.7.1 Set the screen resolution of the presentation screen	45
6.8 Record in 100% scale	45
6.9 Timeline zoom	45
7 Make recordings	46
7.1 Mouse Tracker	46
7.2 Screen-based recording	46
7.2.1 Start a screen-based recording	46
7.2.2 Prepare timeline	47
7.2.3 Moderate a screen-based recording	47
7.2.4 Perform a screen-based calibration	48
7.2.5 End a screen-based recording	50
7.2.6 Lab browser	51
7.3 Scene camera recording	52
7.3.1 Start a scene camera recording	52
7.3.2 Moderate a scene camera recording	53
7.3.3 Perform a scene camera-based calibration	53
7.3.4 End a scene camera recording	54
7.3.5 Participant Variables	54
7.4 External Presenter recording	56

7.4.1	Connect the computer and Tobii Pro eye tracker	56
7.4.2	Configure the computer and install the software	56
7.4.3	Connect Tobii Pro Lab to the eye tracker	57
7.4.4	Set up an experiment in E-Prime Extension and Tobii Pro Lab	57
7.5	Rename a recording	58
7.6	Export calibration results	59
7.7	Change device settings	59
7.7.1	Add a device	59
7.7.2	Remove a device	59
7.7.3	Change participant camera and microphone settings	60
7.7.4	Change the eye tracker settings	60
7.7.5	Change TTL parallel port settings	61
7.8	Stimuli presentation timing data	62
7.8.1	Modify screen latency settings	62
8	Biometrics	63
8.1	Set up a Shimmer GSR sensor	63
8.2	Use the GSR Data Chart	63
8.3	GSR data filtering and analysis	64
8.3.1	GSR data filtering	64
8.3.2	SCR detection analysis	64
8.3.3	GSR metrics	65
8.4	Export GSR Data	66
9	Analyze data	67
9.1	Replay a recording	67
9.1.1	Replay controls	67
9.1.2	Replay settings	67
9.1.3	Timeline tools	68
9.1.4	Replay tools	71
9.2	Multi-select recordings	76
9.3	Export recordings	76
9.3.1	Export a video clip	76
9.3.2	Export a whole recording	77
9.4	Map eye tracking data	78
9.4.1	Assisted and manual mapping	78
9.5	Visualize eye tracking data	82
9.5.1	Create and customize a Heat map	82
9.5.2	Create a Heat Map	82
9.5.3	Customize the appearance of a Heat Map	83
9.5.4	Save a Heat Map as an image	84
9.5.5	Create and customize a Gaze Plot	84
9.5.6	Create a Gaze Plot	84
9.5.7	Customize the appearance of a Gaze Plot	85
9.5.8	Save a Gaze Plot as an image	86
9.5.9	Modify the Gaze Filter setting	87
9.5.10	Change the Gaze Trail settings	88
9.6	Metrics visualization	88
9.6.1	Create a Metrics Visualization	89
9.6.2	View a Metrics Visualization	89

9.7 Times of Interest (TOI)	90
9.7.1 Automatic Times of Interest	91
9.7.2 Custom Times of Interest	92
9.7.3 Add offsets to Custom Times of Interest	93
9.8 Areas of Interest (AOI)	94
9.8.1 Draw AOIs on Snapshots, image stimuli, or frames saved as media	94
9.8.2 Dynamic AOIs	95
9.8.3 Edit the shape of an AOI	95
9.8.4 Change the size of an AOI	95
9.8.5 Move an AOI	96
9.8.6 Cut, copy or paste AOIs	96
9.8.7 Change the name, appearance, or visibility of an AOI	97
9.8.8 Undo or redo AOI actions	97
9.8.9 Pan and zoom in the AOI editor	97
9.8.10 Activate or deactivate an AOI	98
9.8.11 Export an image of the AOIs on the image or Snapshot	98
9.9 AOI tags	98
9.9.1 Create AOI Tags	99
9.9.2 Aggregated AOIs	101
9.9.3 Export and import AOIs	104
9.9.4 Automatic AOIs (Text stimulus)	105
9.10 Analyze a Screen Recording	106
10 Metrics	107
10.1 Understand the metrics	107
10.2 Export metrics data to a file	109
10.3 Binned metrics	111
10.4 Reading metrics	111
10.5 Metric export formats	113
10.5.1 Interval-based TSV file	114
10.5.2 AOI-based TSV file	124
10.5.3 Event-based TSV file	131
10.5.4 Excel Report	133
11 Data export	138
11.1 Export eye tracking, GSR and recording data in a *.tsv file	138
11.1.1 How to export data	138
11.2 Data export formats	138
11.2.1 Computer timestamp (screen-based)	139
11.2.2 Computer timestamp (wearable)	139
11.3 Data export information	139
11.3.1 General	139
11.3.2 Eye tracking data	141
11.3.3 Media	142
11.3.4 Gaze events	143
11.3.5 Web data	144
11.3.6 Other sensor data	145
11.4 Plof file format	145
Appendix A Keyboard shortcuts	146
A1 Dashboard	146

A2 Design	146
A3 Recording	146
A4 Web recording	146
A5 Replay	146
A6 Areas of interest (AOI)	147
Appendix B Coordinate systems and data mapping	148
B1 Display Area Coordinate System (DACS)	148
B2 Scene Camera Projects and the Media Coordinate System (MCS)	149
B3 Head Unit Coordinate System (HUCS)	151
B4 Map data from Glasses recordings onto Snapshots	152
Appendix C Gaze Filter functions and effects	154
C1 Raw Gaze Filter preset type	154
C2 Tobii I-VT Gaze Filter definition	154
C2.1 Gap fill-in (interpolation)	155
C2.2 Eye selection	155
C2.3 Noise reduction (Moving Median)	155
C2.4 Noise reduction (Moving Average)	156
C2.5 Velocity calculator	156
C2.6 I-VT fixation classifier	157
C2.7 Merge adjacent fixations	157
C2.8 Discard short fixations	157
C2.9 Tobii I-VT Gaze Filter default parameters	158
C3 Effect of mapping data onto Snapshots with Tobii I-VT or Raw Gaze Filter	159
Appendix D Calculate heat maps	160
D1 Basic concept	160
D2 The Heat Map	161
Appendix E Pro Lab project and media limitations	164
Appendix F Analysis of an External Presenter project recorded by E-Prime	165
Appendix G Support, Training, and Warranty	167
G1 Customer Support	167
G1.1 Get help online	167
G2 Training and Education Services	167
G3 Warranty information	167

1 General Information

1.1 Welcome to Tobii Pro Lab

Congratulations on your purchase of Tobii Pro Lab!

Pro Lab adds efficiency to every step of your research process with a visual and functional user interface and specialized software features. You can easily create complex experiments, collect eye tracking data, observe and qualitatively analyze individual recordings, and aggregate data for quantitative analysis and visualization - all in one solution.

Pro Lab's analysis and visualization tools allow you to easily process and prepare large volumes of eye tracking data for useful comparison, interpretation, and presentation. You can calculate a variety of eye tracking metrics and create visual representations of your data thus obtaining a top-down overview of your data and report findings.

Pro Lab is built on a platform that ensures precise and consistent timing accuracy - down to the millisecond. Thanks to transparent tools for controlling your collected data, you can rely on your analysis and findings with confidence. This includes customizable gaze filters, full access to all raw gaze data, as well as insights into metric calculations.

1.2 Symbols used in this document

Three different symbols are used in this document:



The Information symbol means something is important or needs special attention.



The Tip symbol denotes additional information that can make a process or function easier.



The Warning symbol means there is a possible risk of harm if the warning is ignored.

1.3 How to collect diagnostic data for Pro Lab

1. Open Pro Lab.
2. In the main menu on the left, select **Help and Learn**.
3. Select **Collect diagnostics data for Pro Lab**.
4. Choose a location to save the diagnostic file.
5. Click **OK**.
This generates a .tar file that you can provide to Tobii Pro Support.
6. A confirmation message displays when the file has been completed.

1.4 How to cite Tobii Pro Lab

Here are two examples (APA and BIBTEX) of how to cite Tobii Pro Lab software or this manual.

Please make sure to update the version number in the example in case it does not match the version you are citing. The current version number of your Pro Lab software can be found in the *About and Updates* section.

1.4.1 APA example for Tobii Pro Lab (software)

Tobii Pro AB (2014). Tobii Pro Lab (Version x.xx) [Computer software]. Danderyd, Sweden: Tobii Pro AB.

1.4.2 APA example for Tobii Pro (User Manual)

Tobii Pro AB (2014). Tobii Pro Lab User Manual (Version x.xx). Tobii Pro AB, Danderyd, Sweden.

1.4.3 APA's own tex files (6th edition)

```
@SOFTWARE{TobiiProLab,  
ENTRYSUBTYPE = \{Computer software},  
AUTHOR = {Tobii Pro AB},  
TITLE = {Tobii Pro Lab},  
VERSION = {1.145},  
PUBLISHER = {Tobii Pro AB},  
URL = {https://www.tobii.com/},  
DATE = {2014}  
}
```

1.4.4 BIBTEX entries for Pro Lab (software and User Manual)

```
@software{TobiiProLab,  
author = Tobii Pro AB,  
title = {Tobii Pro Lab},  
url = {http://www.tobii.com/},  
howpublished = {Computer software},  
address = {Danderyd, Stockholm},  
version = {1.145},  
year = {2014},  
}
```

```
@misc{TobiiProLab,  
author = Tobii Pro AB,  
title = {Tobii Pro Lab},  
url = {http://www.tobii.com/},  
howpublished = {Computer software},  
address = {Danderyd, Stockholm},  
type = {bibcomputersoftware},  
version = {1.145},  
year = {2014},  
}
```

```
}  
@manual{TobiiProLabManual, author = Tobii Pro AB,  
title = {Tobii Pro Lab User Manual},  
url = {http://www.tobiiipro.com/},  
howpublished = {Computer software},  
address = {Danderyd, Stockholm},  
version = {1.145},  
year = {2014},  
}
```

2 Introduction

Tobii Pro Lab provides a comprehensive platform for the recording and analysis of eye gaze data, which helps in the interpretation of human behavior, consumer responses, and psychology. Combining simple preparation for testing procedures and advanced tools for visualization and analysis, eye-tracking data is easily processed for useful comparison, interpretation, and presentation. A broad range of studies are supported, from usability testing and market research, to psychology and oculomotor physiological experiments. Pro Lab's intuitive workflow, along with its advanced analysis tools, enables large and small studies in a timely and cost-efficient way without the need for extensive training.

In addition to offering powerful analysis tools, Pro Lab is also designed to work with other software commonly used for recording and analyzing data. This is done by synchronizing with recording software using TTL, as well as by enabling data exports in standardized formats, for example for Microsoft Excel, Matlab, and SPSS.

2.1 Pro Lab editions

Pro Lab is available in three editions: **Full Edition**, **Screen-Based Edition**, and **Analyzer Edition**.

The editions have different combinations of enabled modules (Design, Record, Analysis) and compatible hardware (screen-based and wearable).

- The **Screen-Based Edition** contains all three modules and works with Tobii Pro Screen-based hardware. The featured project types are Screen-based and Scene Camera. In addition, External Presenter projects created with the Full Edition can be opened.
- The **Analyzer Edition** contains only the Analyzer module and features the Glasses project type to create new projects and import data recorded with Tobii Pro Glasses 2 & Tobii Pro Glasses 3. This edition opens Screen-based, Scene Camera, and External Presenter projects recorded with either the Full or Screen-based Edition.
- The **Full Edition** contains all modules (Design, Record, Analyze) and works with all Tobii Pro hardware; screen-based and wearable. It opens and creates all project types: Screen-based, Glasses, Scene Camera, and External Presenter.



Tobii Pro Lab is continuously being developed and refined. Please visit tobiipro.com for the most recent specifications for the software and for the latest version of this document.



VR 360 was discontinued in releases after version 1.162. This was the last version to support the VR 360 project type. For more information, read [System requirements for Tobii Pro and VR-capable setups](#) on [Tobii Pro Connect](#).

2.2 Project types

Data about study layout, stimuli, participants, and recordings is stored in Pro Lab projects. There are four kinds of projects:

- **Screen project:** when stimuli is presented on a screen by Pro Lab and gaze data is captured using a screen-based eye tracker
- **Glasses projects:** when data is collected by Pro Glasses 2 or Pro Glasses 3
- **Scene Camera projects:** when you use an external video camera to record events in the real world
- **External Presenter projects:** when you use third-party software such as E-Prime together with Tobii Pro Lab



In Glasses projects, only the Project Overview and Analyze module are available.



Each project can contain many recordings, participants, timelines, snapshots, mapped data, and events, etc.

2.3 How Pro Lab is structured

Pro Lab software is comprised of three modules: the **Design** module, the **Record** module, and the **Analyze** module, all of which are accompanied by a **Project Overview** tab. Depending on which edition (license) of Pro Lab you use, available modules and functions vary. Access to the modules also varies depending on what kind of project you are working on.

The following tables show you which features are available with each edition.

2.3.1 Project Overview

Regardless of whether you use a screen-based or a wearable eye tracker, the Project Overview section provides information about the elements of your project, such as what recordings are in it and what Events are associated with the recordings. It also provides quick access to some analysis tools.

Project Overview is included in all licences.

2.3.2 Design

You can create experiments in the Design module based on Timelines containing different stimuli. You can edit stimuli presentation settings like display position, background color, presentation time and stimulus advancement methods, (i.e. end on a mouse click or key press to adapt your experiment). In this module, you also get a preview of what the stimuli will look like on the screen. The table below shows the features available for each edition of Pro Lab.



The Design module works with selected screen-based eye trackers from Tobii Pro. It does *not* work with Tobii Pro Glasses 2 and Tobii Pro Glasses 3.

The table below shows the Design module features available for each edition of Pro Lab.

Feature	Screen-based	Analyzer	Full
Design experiments with multiple timelines or use hierarchical structures with randomized presentation (shuffled order, randomized sampling), and repetitions with	•		•

Feature	Screen-based	Analyzer	Full
stimuli			
Add text stimuli with automatic Area of Interest definitions	•		•
Batch editing of stimuli settings	•		•
Use multiple stimuli advance options, either alone or in combination (advance on time, key press, mouse click)	•		•
Configure stimulus onset markers (TTL) for synchronization purposes	•		•
Designate a gaze trigger zone to advance to next stimulus when viewed.	•		•

2.3.3 Record

The Record module lets you configure eye trackers from Tobii Pro and present different stimuli, with high timing accuracy. You can read more about this in the Tobii Pro Learning article [Stimulus presentation timing in Tobii Pro Lab](#). You can validate a calibration, record eye tracking data, mouse clicks and key presses, as well as Galvanic Skin Response (GSR) data from Shimmer3 devices. The participant camera with audio lets you record the participant. The Record module turns into a Moderator view during live viewing of the track status, stimuli displayed, and gaze data.



The Record module works with selected screen-based eye trackers from Tobii Pro. It does *not* work with Tobii Pro Glasses 2 and Tobii Pro Glasses 3.

The table below shows the Record module features available for each edition of Pro Lab.

Feature	Screen-based	Analyzer	Full
Scene camera project (support for real world experiments using screen based eye trackers)	•		•
External Presenter project			•
Configure eye tracker settings	•		•
Define experiment participants	•		•
Calibrate eye tracker (regular and infant calibration)	•		•
Numeric calibration results (accuracy and precision values)	•		•
Present image and video stimuli	•		•
Record eye tracking, mouse, and	•		•

Feature	Screen-based	Analyzer	Full
keyboard data			
Recording of galvanic skin response data from Shimmer3 GSR+ sensors	•		•
Moderator view: track status, stimuli displayed and gaze data live	•		•
Send stimulus onset markers (TTL) for synchronization purposes	•		•
Receive TTL-in markers and the value for synchronization (available for Pro Spectrum and Tobii Pro TX300 eye trackers only)	•		•
Participant camera	•		•
Present webpages and make screen recordings	•		•

2.3.4 Analyze

The Analyze module enables you to replay, visualize and analyze your recorded data. It provides data-filtering features, visualizations and the ability to export data for presentations and for further processing in third-party software. In addition, it also provides assisted and manual mapping.

The list below shows the Analyze module features. They are available for *all* editions of Pro Lab unless otherwise noted.

Analyze module features:

- Replay of recordings
- Import Tobii Pro Glasses recordings*
- Manual mapping onto Snapshot images
- Assisted mapping onto Snapshot images
- Create and edit static and dynamic Areas of Interest (AOIs) on images and videos
- Add Areas of Interest on text stimulus (character, word, sentence)
- AOI Tags and Grouping by tags
- Log Events for behavioral coding
- Times of Interest: define time intervals based on recording and logged Events
- Selecting a frame as background and pairing it with Time of Interest (Screen and Scene camera projects only)
- Plot gaze x and y coordinates as well as eye movement velocity over time
- Metrics Visualizations let you plot how metrics are affected by experiment conditions (AOI tag groups, stimulus variables, or participant variables)

- Plot and visualize galvanic skin response (GSR) data over time (together with gaze video replay and eye movements)**
- GSR data analysis: noise reduction filters and detection of Skin Conductance Responses (SCRs) and Event Related SCRs**
- Static Heat Map Visualizations on images
- Static Gaze Plot Visualizations on images
- Video export of recordings and recording segments
- Export eye tracking metrics
- Export Event and time interval based metrics
- Export GSR metrics**
- Export binned metrics
- Export visualizations as images (.png and .jpg)
- Export numeric calibration results (accuracy and precision values)
- Export calibration results as images (.png format)
- Recording data to text file (.tsv)

*not available for Screen, Scene Camera, and External presenter projects

**not available for Glasses projects

2.4 Pro Lab licenses

In order to use Pro Lab, you need to have a license for it. Pro Lab has two different licenses: a perpetual license and a subscription-based license.

A perpetual license grants you one year of free upgrades. One- to four-year upgrade contracts are available for perpetual licenses.

A subscription license provides you with access to the latest software versions as soon as they become available.

Each license is associated with a specific edition of the software: Full Edition, Screen-Based Edition, or Analyzer Edition. The different editions provide access to different modules and features in the software. Read more about [Pro Lab editions](#).



A license can only be active on one computer at a time. If you attempt to use a license that is already active on another computer, you will be asked to first deactivate the license on the other computer.

2.5 System requirements

For the most up-to-date information about Pro Lab's software system requirements, please read [Minimum System Requirements for Tobii Pro](#) on tobiipro.com.

3 Install or update Pro Lab

This chapter describes how you install and update Pro Lab on your Windows computer. For optimal performance, Pro Lab should be installed on a computer that adheres to the Pro Lab system requirements. For more information, read [System requirements](#). You can also review the latest system requirements online at tobiipro.com or ask your sales representative.

3.1 Windows versions



Tobii Pro Lab is designed and tested to run on Windows 10 Pro.

3.2 Pro Lab license key

You will receive a license key(s) and a link to the latest version of Pro Lab via e-mail when you purchase Tobii Pro software.



You need the license key information during your first installation on any new computer, so please keep it in a convenient place for future reference.

3.2.1 Activate or deactivate a Pro Lab license

How to activate a license:

1. Ensure your computer is connected to the internet.
2. Start Pro Lab.
3. Select **License** from the left-hand menu under "Tobii Pro Lab."
4. Enter the license number into the License key field.
5. Click **Activate**.

How to deactivate a license:



Before deactivating the license, check that you have saved the license key in a safe place. When the deactivation is completed, the software has no record of the license key. This means that if you want to use the key again on the computer, you must enter it manually.

1. Ensure your computer is connected to the internet.
2. Start Pro Lab.
3. Select **License** from the left-hand menu under "Tobii Pro Lab."



Warning! The Deactivation button is final. After pressing Deactivate, you will need to manually reenter the license key.

4. Click the **Deactivate** button.

3.3 Install Pro Lab

Pro Lab uses several Windows programming components and drivers. If they are not already on your computer, the Pro Lab Installer automatically installs them.

How to install Pro Lab:

1. Download the installation file for the latest version of Pro Lab, go to tobiipro.com
2. Locate and run the installation file on your computer and follow the instructions on the screen.
3. When prompted to accept the license agreement, read the license agreement and click **Accept**
4. When the installation is finished, start the application.
5. Enter the license key, provided in the product purchase e-mail, in the license key input fields and click **Activate**.



You must have local admin rights when installing Tobii Pro Lab.

If more than one person uses Pro Lab on the same computer, you must enter the license key for each user.

One license can only be used on one computer at a time. However, you can deactivate the license on one computer and activate it on another. For more information, read [Pro Lab licenses](#).

You cannot utilize roaming user profiles on the computer running Tobii Pro Lab.



If nothing happens when you start Tobii Pro Lab and there are no error messages, check which edition of Windows 10 you have on your computer. If it is one of the Windows 10 N editions, you need to install a Media Feature Pack from Microsoft. For more information, read the Connect article, [Tobii Pro Lab does not start](#).

3.4 Upgrade contracts

Tobii Pro Lab is license-based software.

When you purchase a software license key, your purchase includes a one-year upgrade contract. The upgrade contract gives you access to updates for your software without additional charges in this one year period.

When your upgrade contract expires, your software keeps working as usual, but you will not be able to access any updates.

You can purchase upgrade contracts for additional years. The [Tobii Pro Sales Team](#) can help you find a plan that works for you. If you are using a subscription-based or temporary license, the upgrade contract is included in your purchase. When your subscription or temporary license expires, the upgrade contract expires as well, and your license will not be usable anymore.

If you would like to renew your subscription or purchase a temporary license, our [Tobii Pro Sales Team](#) can help you find a plan that works for you.

3.5 Update Pro Lab

If the computer running Pro Lab is connected to the internet, it automatically checks for application updates. If a newer version is available, a red symbol appears at “About and Updates” in the left-hand navigation panel when you start the software.

How to update Pro Lab:


1. Ensure your computer is connected to the internet.
2. Start Pro Lab.
3. Click **About and Updates** in the left-hand navigation panel.
4. Click **Install** and wait for the download to finish.
5. Click **Install now**. The installation is completed automatically, and the new version of the application starts as soon as the installation is finished.



If the computer running Pro Lab doesn't have an internet connection, you can download the installer manually from [Tobii Pro Connect](#).

4 Start a project

When you start Pro Lab, you have the option to create a new project, open an existing project, enter or modify your license details, check for software updates, and get help with the software. You can also access these options while working on a project by clicking the Pro Lab logo and text in the top left corner.




Read more about [Project types](#).

4.1 Create a new project

How to create a project:

1. Click **Create new Project**.
2. Select one of the project types.



3. Enter the project's name in the Name field. If you don't enter a name, Pro Lab suggests one, e.g., "Project1".
4. Select the folder where the project should be saved by clicking **Browse**. If no selection is made, the project is saved in a subfolder named "Tobii Pro Lab" in the Documents folder of the user currently logged on to the computer.



Pro Lab doesn't support projects stored on network disks, USB-disks, or folders synchronized to cloud storage services (such as Dropbox, OneDrive, etc).

5. Click **Create**.

When the project has been created, the Project Overview tab displays.

4.2 Project Overview

The Project Overview section is the starting point for working with a project in Pro Lab. It lets you easily jump to project specific content. You can export and import projects and more, as well as access many more tools at the top of Project Overview window. The tabs on the left-side toggle on and off the panels—when the icon is orange, the panel is displayed.

The amount and type of panels available on the left varies depending on the type of project.

Recordings: This lists all the recordings in the project. For more information, read [Manage recordings](#).

Participant Variables: Participant Variables are used to filter data when generating visualizations, calculating eye tracking metrics and comparing the behavior of different participant groups. Click the plus button to add a new Participant Variable. For more information, read [Participant Variables](#).

Participants: This is a list of the participants in your study, including which variables have been specified. Click the plus button to add a new Participant. For more information, read [Participant Variables](#).

Snapshots: Snapshots are still images of environments and objects of interest. Click the plus button to add a new Snapshot. For more information, read [Map eye tracking data](#).

Events: An Event Type is a definition of an Event that you want to mark in your recording. Click the plus button to add a new Event Type. For more information, read [Manage Event Types](#).

Stimulus Variables: A Stimulus Variable is a specific feature of your stimulus that is a part of the context in which a behavior occurs. Click the plus button to add a new Stimulus Variable. For more information, read [Stimulus/group variables](#).

5 Manage projects

5.1 Open an existing project

When you want to continue working with an existing project, select the *Open existing Project* option.

How to open an existing project:

1. Click on the *Open existing Project* option. A list of recently-opened projects appears. Information displayed about each project includes the name of the project, when it was created, and where it is stored.
2. Double-click the desired project.
3. If the desired project isn't in the list, click the **Browse** button at the top right. Locate the .project file associated with the desired project and click **Open**.

If the project is already open, the Project Overview window will display.



If you want to open a project that was created in an older version of Pro Lab, the project will be migrated to the version of Pro Lab you currently have. When you attempt to open the project from the list, a pop-up message explains that after migration, your project won't open in earlier versions. It is important to back up your project if you plan to open it later in an older version of Pro Lab. Make sure the backup checkbox is checked in the pop-up message.

5.2 Rename or remove a project

In the **Open existing project** window, right-click a project name and select “Rename Project” or “Remove project from list.” The project *will not be deleted*, only removed from the list of recently opened projects.

5.3 Export and import a project

There are three different reasons you may need to export and import a project.

- Your project may have more than one person collecting data (distributed data collection) and you need to share the recordings. For more information, read [Distributed data collection](#).
or...
- Your project may have more than one person coding data (distributed manual coding) and you need to share the recordings. For more information, read [Distributed manual coding](#).
or...
- You want to back up your project using the export/import functionality. For more information, read [Export and import a project](#).

5.3.1 Distributed data collection

Researchers can save time and money by doing multiple data collection sessions across different locations and eye tracking setups. For example, you can reduce the time needed to reach the goal sample size by running simultaneous data collection sessions with two or more eye tracking setups, or collaborate with other research teams to perform large scale and/ or geographically separated studies. In this section you will learn how to conduct a distributed data collection in Pro Lab.

In this example, Researcher A designs the project and send it to a data collector, Researcher B. Researcher A might even have more than one data collector for their project, in this case Researcher C and D.

If independent variables are used in this “best practice example,” they are related to participant characteristics so the project should include predefined variables and values. These need to be assigned to the listed participants. Typically researchers create a protocol with information on recruitment and data collection procedure. For more information, go to the [Tobii Pro Learning Center](#).

Researcher A:

1. Create a project. For more information, read [Create a new project](#).
2. In Project Overview, create the participant variables and values.
3. Define the participant variables, values, and participants:
 - a. **Independent variables** (these affect the recruitment related to a predefined sampling order):
Create all the participants that should be recorded and specify their variables and variable values
For more information, read the online Learning article [Participant Variables and Groups](#) in the Learning Center.
 - b. **No independent variables:**
Specify interesting variables for the Researcher B to provide for each recorded participant.
For more information, read the online Learning article [What are experimental variables?](#) in the Learning Center.
4. On the **Design** tab, create the timeline(s) that should be recorded, including all stimuli, groups, timings, etc.
5. On the **Project Overview** tab, select **Export >Project** at the top.
6. In the **Project Export** dialog, in the **Export type** list, select "Project for data collection" in the drop-down and divide the participants into different child projects by selecting/deselecting the participants. For example, if there are 10 participants, you could check five of them to send to Researcher B and later choose the remaining five to send to Researcher C.
7. Click the **Export** button after choosing where to save the file locally. The file will be saved in Zip format.
8. Repeat the process to select and send other Participants to Researcher C, etc.



You don't have to wait for one project to be finished before queuing the next one. When all the child projects are created, you can close the dialog.

9. Send the generated Zip files to Researcher B and C.



After a Participant has been exported from the list, a double check mark displays next to the name in the list until the Project Export window is closed.

Researchers B and C:

1. Unzip the project.
2. Open the unzipped project in Pro Lab.
3. Recruit participants (typically completed earlier according to the recruitment protocol):

- a. **Independent variables:**

Match the specified variable values specified by Researcher A.

- b. **No independent variables:**

No match required.



If additional participants need to be added, make sure to specify the variables and values that have been applied to previous participants.



Do not change anything in the Design module because this can make the analysis phase more difficult.

4. Make participant recordings.
5. On the **Project Overview** tab, select **Export >Project** at the top.
6. In the **Project export** dialog, select the export type "Project with recordings" in the drop-down list to view the list of recordings. Deselect any recordings you do not want to include.
7. Click the **Export** button after choosing where to save the file locally. The file will be saved in Zip format.
8. Send the generated Zip file back to Researcher A.

Researcher A:

1. Receive the zipped project file from Researcher B.
Pro Lab can import a zipped file and you do not need to unzip it elsewhere.
2. Open the parent project and in Project Overview, select **Import > Project**.
3. Select the Zip file from the stored location.
4. In the Import dialog, select which recordings to import (typically all) and how to handle any merge conflicts (typically different participant names).
5. After import, review the confirmation dialog and the potential list of conflicts that were generated during import. The list displays so that you know about any duplicates that were created. Duplicates can be identified by the number in parenthesis after the name. You can manually adjust or delete any duplicated participants, variables, and/or values afterwards

in Project Overview.

6. Repeat steps 1-6 for all the files you receive from external data collectors.
7. Start analyzing the aggregated data.

5.3.2 Distributed manual coding

When there is a lot of coding to do, it may help to do the work simultaneously on several computers.



Conflicting data during import can be overwritten so we advise that you make a backup of the target project before doing an import. For more information, read [Distributed manual coding](#). This safeguard enables data recovery in case the data import causes the unplanned deletion of data that should have been preserved.

How to export coding data from another project:

1. In **Project Overview**, select **Export >Project** at the top.
2. In the **Project Export** dialog, select the Export type "Send events and mapped data" in the drop-down list at the top.
3. Click the **Export** button after choosing where to save the file locally. The file will be saved in the .plex format and will contain all the coding-related data.
4. Send the generated .plex file back to the researcher doing the analysis.

How to import coding data from another project:

1. Open the project into which you will import the coding data.
2. In **Project Overview**, select the **Import** drop-down menu.
3. Select "Coding data."
4. Use the file browser to locate the project file (*.plex) associated with the parent project.
5. Click **Open** and select either:
 - a. *Keep data from existing Project*: If there is conflicting data, you want to keep the data in the currently-open project and discard the corresponding imported data; or
 - b. *Replace existing data with imported data*: If you want to overwrite the data in the currently-open project with imported data.
6. Click **Start import**.
7. When the import is finished, click **Done**.

5.3.3 Back up projects

We strongly recommend that you back up projects regularly and store the backups on multiple external drives.

How to export a full project backup:

1. In Project Overview, select **Export >Project** from the menu at the top.
2. In the Project Export dialog, select *Full project* from the Export type drop-down list. Everything in the project will be included.

3. Browse to where you want to save the file.
4. Click the **Export** button.

5.3.4 Import data from SD card

When you use either Pro Glasses 2 or Pro Glasses 3 to record data, you save the data on a SD card. You then need to import that data into a Glasses project in Pro Lab in order to be able to analyze it. Importing data from the SD card differs slightly between Pro Glasses 2 and Pro Glasses 3 as explained below.

Pro Glasses 2:

How to import recordings from Pro Glasses into a Pro Lab Glasses project:

1. Open the Glasses project into which you will import the recordings.
2. Go to Project Overview and select the **Import** drop-down menu.
3. Select *Glasses 2 Recording*. A file browser opens.
4. Locate the data file (*.ttgp) on the SD card using the file browser.
5. Import the data by clicking the **Open** button.
6. Select the recordings you want to import.

Pro Glasses 3:

1. Open the Glasses project into which you will import the recordings.
2. Go to Project Overview and select the **Import** drop-down menu.
3. Select *Glasses 3 Recording*. A folder browser opens.
4. Locate the SD card root folder using the folder browser.
5. Import the data by clicking the **OK** button.
6. Select the recordings you want to import.

5.4 Manage recordings

The Recordings list contains information about the recordings in the current project. Details displayed include:

- **Recording:** The name of the recording
- **Participant:** The name of the participant
- **Timeline:** The name of the recording's timeline
- **Duration:** The duration of the recording
- **Date:** The time and date when the recording was performed
- **Resolution:** The resolution of the recording
- **Gaze Samples:** The percentage is calculated by dividing the number of correctly identified eye-tracking samples by the theoretical maximum. An eye tracker with a 50 Hz sampling frequency generates 50 samples per second. If the software can use all samples to calculate gaze points, the value in the Gaze Samples column would be 100%. However, this

percentage is rare, because some samples are always lost due to the participant blinking, or looking away from the monitor in the case of a screen-based eye tracker. Blinking usually causes around 5-10% data loss during a recording.

5.4.1 Replay

- Open and replay a recording in the **Recordings** list by double-clicking it or by right-clicking and selecting *Open*. The replay tab displays in the **Analyze** module.
- If you want to open recordings *without* getting redirected to the **Analyze** module, right-click the recording you want to replay and select the option *Open in background*. Pro Lab then opens the recording in the replay tab in the Analyzer module, but you won't see it until you go to that module by clicking the Analyze option in the top menu.



You can start and stop the recording replay using the space bar.

For more information, read [Replay a recording](#).

5.4.2 Delete

If you no longer want to keep a recording in a project, you can delete it.



A deleted recording is permanently removed from the system. It cannot be restored after deletion.

How to delete a recording:

1. Right-click on the recording and select *Delete*. A warning message appears.
2. Confirm your choice by clicking the *Yes, delete this Recording* option.

5.4.3 Import and Export recordings

For more information about distributed data collection, read [Export a whole recording](#). For more information about exporting and importing a recording, read [Export and import a project](#).

5.5 Manage Event Types

The Event Types list in the Project Overview lists all Event Types that are associated with a project. An Event Type is a definition of an Event that you want to mark in your recording. The Event Types are used to highlight either single Events or the start and end of a sequence, i.e. a Time of Interest. In the Replay tab, you can mark a specific point in time with an Event Type to indicate something of interest. This creates an *Event* (an instance of that Event Type).

Details displayed in the Event Types list include:

- **Color of the Event marker:** The colored symbol to the left of the Event Type name represent the marker that will be displayed on the recordings Timeline when an Event of that type has been created. Each Event Type has a unique color that is assigned by Pro Lab.
- **Name:** The name of the Event Type, assigned when the Event Type was created.
- **Shortcut:** A keyboard shortcut can be used to add an Event of that type to a recording Timeline in the Replay tab in the Analyze module.

5.5.1 Create an Event Type

How to create a new Event Type:

1. Click the plus symbol (+) in the top right corner of the Event Type list.
2. Enter a name of the Event Type in the left input field of the dialogue window. If you don't enter a name, Pro Lab suggests one, e.g. My Event001.
3. In the drop-down menu to the right, select which keyboard shortcut creates an Event of that Event Type on the Timeline.
4. Click *OK*.
5. The named Event Type appears in the Event Types list.

5.5.2 Edit an Event Type

How to edit an Event Type:

1. Put the mouse pointer on the row of the Event Type you want to edit. A pen icon appears
2. Click the pen icon. A menu appears.
3. Edit the Event Type details you want to modify.
 - Edit the Event Type's name by selecting the current name and entering a new.
 - Edit the associated keyboard shortcut by opening the drop-down list and making a new selection.
4. Click *OK*.

5.5.3 Delete an Event Type



Deleting an Event Type also deletes all instances of it in every recording in the project, along with all Times of Interest associated with it.

How to delete an Event type:

1. Put the mouse pointer on the row of Event Type you want to edit. A trash can icon appears
2. Click on the trash can.
3. Click *Delete Event*.
4. Click *OK*.

5.6 Analysis tools

The Project Overview contains links to the analysis tools in the Analyze module. When you click a tool, the software switches to the Analyze module and opens a tab for the indicated tool. You can also access these tools by clicking the drop-down symbol next to the Analyze option in the top navigation.

Analysis tools provide the following options:

- **Visualizations:** Here you create Visualizations (Heat Maps or Gaze Plots) based on the gaze data on top of your stimuli or Snapshots. Visualizations only appear on a Snapshot

image if data has been mapped onto the Snapshot.

- **AOI editor:** Here you draw areas of interest (AOIs) on your stimuli or Snapshot images. AOIs enable numerical/ statistical analysis based on regions. When an AOI has been created, you can export eye-tracking metrics for the stimulus or Snapshot on which it is created.
- **Metrics Export:** Here you export your eye tracking metrics based on AOIs or Events for further analysis in third party software, such as SPSS, Microsoft Excel, and MATLAB.
- **Data Export:** Here you export your eye-tracking data for further analysis in third-party software, such as SPSS, Microsoft Excel, and MATLAB. Unlike in Metrics Export, the data available here is not tied to AOIs. Instead, you access raw gaze data, such as gaze points in different coordinate systems, pupil diameters for each eye, the eye position, and information about the recordings in general. Information about whether gaze data fell within AOIs is also included, but metrics regarding AOIs are not provided.



You can choose between two file formats, depending on which the intended third-party application. Go to the Settings pane and choose select either *Standard*, (a .tsv file suitable for any application, e.g. MATLAB), or *Excel compatible* (suitable for use in Microsoft Excel). Both formats have the same data content.

6 Design an experiment

Before you do a Recording in Pro Lab, you need to design an experiment in the Design module. An experiment can consist of one or more Timelines. A Timeline outlines the order in which stimuli are displayed to the study participant.



When you work with a Glasses project, the Design module is not visible because Pro Lab doesn't handle stimuli presentation or recording in Pro Glasses 2 or Pro Glasses 3. Hence, you don't need access to the Design module for Glasses projects.



When you work with a Scene Camera project, the Timeline is not visible because the stimulus is the video feed from the external camera.

6.1 Timelines

A Timeline can contain multiple stimuli of different types. Stimuli types that are currently supported include

- images (.bmp, .gif, .jpg, and .png)
- videos (.avi and .mp4)
- screen recording stimuli
- webpages
- calibration stimuli

6.1.1 Create a new timeline

How to create a new Timeline:

1. Click **Design** in the top navigation.
 - If no Timeline is available, a Timeline with one stimulus gets created automatically.
 - If there already are Timelines, you create a new one containing one stimulus by clicking **New Timeline** in the **Add** section of the Design module interface.
 - The automatically included stimulus is a calibration object. When it is reached during a recording, the participant undergoes a procedure that calibrates the eye tracker to their eyes. Data collected before the calibration most likely has poor accuracy and precision and should therefore not be used for analysis. Read more about [Calibration stimulus](#).
2. Add more stimuli in one of the following ways. (When stimuli get added to the Timeline, visual representations of their contents appear on the Timeline. For more information, read [Stimulus properties](#).)
 - Click the empty stimulus container on the Timeline. A file browser appears; use it to locate the file(s) you want to use as stimuli. Make your selection and click *Open*.
 - Click **Add stimuli** in the **Add** section of the Design module interface. A file browser appears; use it to locate the file(s) you want to use as stimuli. Make your selection and click **Open**.

- Locate the file(s) you want to add as stimuli by using Windows File Explorer and drag them onto the Timeline.
 - Change the order of the stimuli by clicking on those you want to move and drag them it to the desired position (s).
3. Add a Group by clicking *Add Group Element* in the *Add* section of the Design module interface. A new group appears in the timeline. You can drag and drop it to the desired place. For more information, read [Groups](#).
 4. Go to the **Properties** panel, enter the group's name in the *Name* field, and set its desired **Action**, that is, determine in what way you want the content of the group to appear:
 - a. **Show all**: All stimuli in the group are displayed sequentially.
 - *Shuffle*: All stimuli in the group are displayed in a random order. Each stimulus appears only once.
 - b. **Repeat**: The group's content is displayed sequentially several times over. Enter the number of times the sequence shall be repeated. Max value is 100.
 - c. **Sample**: One, and only one, stimulus in the group is displayed. Determine how it will be selected:
 - *Sequential*: Each time the Sample Action is executed, it displays the next stimulus in strict order.
 - *Random with replacement*: Each time the Sample Action is executed, a random stimulus from the group is displayed. It can reappear later in the action.
 - *Random without replacement*: Each time the Sample Action is executed, a random stimulus from the group is displayed. It will not reappear later in the action, unless there are more repetitions than stimuli. For example, if you set 20 repetitions and have 10 stimuli, the cycle restarts when all stimuli have been displayed once and all stimuli get displayed twice.
 5. Add a Screen Recording stimulus by clicking **Add Screen Recording stimulus** in the **Add** section of the Design module. A new Screen Recording stimulus appears in the timeline. (For more information on Screen Recording stimulus properties, read [Stimulus properties](#) and [Areas of Interest \(AOI\)](#) on how to analyze a screen recording).
 6. Add a Text Stimulus by clicking **Add text** in the **Add** section of the Design module. A new Text stimulus appears in the timeline. Double-click the Text stimulus to open the text editor. Enter and edit your text in the text editor.
 7. Add an Instruction element by clicking **Add instructions** in the **Add** section of the Design module. A new Instruction element appears in the timeline. Double-click the Instruction element to open the text editor. Enter and edit your text in the text editor.
 8. The additions are saved automatically. You can continue working with Pro Lab.



Displaying stimuli before the calibration stimulus may make the participant relax or catch their attention before the data collection begins. However, data collected before calibration is not supposed to be used during analysis because its accuracy and precision are highly likely to be poor.

6.1.2 Delete a timeline

How to delete a timeline:

1. Go to Timelines to the right and hover the mouse cursor over the desired Timeline. A trash can icon appears.
2. Delete the Timeline by clicking on the trash can.



Deleting a Timeline doesn't affect the recordings made using it. This deletion means that you cannot make any more recordings with that Timeline, whereas the associated stimuli remain available for analysis.

6.1.3 Duplicate a timeline

How to duplicate a Timeline:

1. Go to Timelines to the right and hover the mouse cursor over the desired Timeline. A duplication icon appears.
2. Duplicate the Timeline by clicking on the icon. A copy of the Timeline appears.

6.1.4 Change the name of a Timeline

How to change the name of a Timeline:

1. Click on the current Timeline name. (When a new Timeline is created, it gets a name automatically, e.g. Timeline1.)
2. Delete the current name and enter the new name.
3. Press "Enter" on the keyboard.

6.2 Stimulus properties

6.2.1 Add a stimulus to the timeline

Click a type of element button in the Add menu or drag and drop your images and video elements from a folder on your computer to the Timeline.



Figure 1. The Add menu on the Timeline in the Design module.

The elements that can be added from the **Add** menu are:

	<p>Add stimuli from your computer:</p> <p>.bmp, .jpg, or .png image files. Pro Lab will scale images according to the screen size.</p> <p>.gif, .avi, or .mp4 movie files - Select video files to be presented to your participants. Pro Lab will scale images according to the screen size.</p>
	<p>Add a group element (adds rules and actions to your presentation)</p>

	Add web stimulus
	Add screen recording stimulus
	Add text (or instructions)
	Delete selected Timeline elements



Pro Lab defines a stimulus, for example an image file, video, or calibration stimulus, as one entity that can be used many times. As a consequence, properties set for that stimulus get applied to all of its instances on all Timelines in which it appears. For example, if the viewing time is changed for a particular image, that also changes the viewing time for all occasions it is shown during all recordings, regardless of what Timelines it appears on. If you want to use the same stimulus with different properties, each image or video of it must get a unique file name.



On the **Design** tab, the **Properties** panel on the right-side of the screen changes depending on which kind of stimuli or element you are working with.

6.2.2 Advance on

The Advance on options you select in the Properties panel on the right are displayed as icons below the stimulus on the Timeline.




Figure 2. Advance on options in the Properties panel




Figure 3. The selected Advance on options are displayed as icons below the stimulus on the Timeline.

6.2.3 Calibration stimulus

Add a calibration stimulus by toggling on **Eye tracking** in the Advanced settings section in the right-hand side panel on the **Design** tab.

Click the Eye tracker Calibration card in the Timeline to display its properties in the Properties panel on the right-hand side.

The following calibration stimulus properties can be changed:

- **Validate calibration:** Determines if a calibration validation will happen after the calibration or not. When toggled on, the calibration will be validated.
- **Randomize target point order:** Determines if the targets are shown in a random order or not. When toggled on, the targets are shown in a random order.
- **Background color:** Determines what color the screen background gets when the calibration is done. The object that the participant will look at is displayed on top of this color.
 - ✔ Tobii Pro recommends using a calibration background color as close as possible to the stimuli background color.
 1. Select background color by opening the drop-down and the color picker.
 2. Click on the desired color or enter a hexadecimal color code in the number entry field.
 3. The background color is updated and appears in the Timeline as the background of the stimuli cards displayed there.
- **Number of targets:** Specifies how many calibration targets (points) are displayed during the calibration procedure. The default value is 5.
- **Calibration Type:** Determines how the advancement from one calibration target to the next is handled.
 - *Timed:* If this option is selected, the next calibration point gets displayed when the eye tracker has signaled that it has sufficient data to proceed or when sufficient time has passed for the collection of a reasonable amount of data.
 - *Manual:* When this option is selected, you proceed to the next stimuli point by pressing the spacebar on the keyboard.
- **Calibration Target:** This property specifies the appearance of the object the participant is asked to look at during calibration. The Calibration Target object gets displayed at as many

different locations on the presentation screen as is specified by the “Number of points” property. (The Calibration Target can only be set if the Calibration Type is set to Manual.)

- **Point:** is selected, the displayed calibration target is a simple filled-in circle.
 - *Point color:* This property specifies the color of the calibration Point.
 1. Select the Point color by opening the drop-down and the color picker.
 2. Click on the desired color or enter a hexadecimal color code in the number entry field.
 3. The Point color is updated and appears in the Timeline as the Point color on the calibration stimuli cards.
 - **Video:** If this option is selected, the study participant is shown a small video clip with accompanying sound in the location of the calibration target. This is particularly useful when the study participant is a child or someone else who has difficulty paying attention for long periods of time.
 - *Select video clip:* Select among seven thumbnails of video clips that you can use as the Calibration Target.
 1. Start a preview by clicking the *Play* button on a thumbnail.
 2. If you click elsewhere on the thumbnail, you set that clip as the Calibration Target for the calibration stimuli on all Timelines.
 3. A small visual representation of it appears on the calibration stimuli cards on the Timeline in the locations where the Calibration Target is displayed during the calibration procedure.
 - If you want to add your own video as a Calibration Target:
 1. Click the plus sign (+) next to the **Add video** option. A selection dialogue displays.
 2. Select your video file.



The video clips can be of many different formats. If Pro Lab can handle a file as video stimuli, then it should work as a Calibration Target as well. Note that the file is not scaled, so it should have a suitable size for your screen.

- **Resizing Video:** This calibration target is designed to maintain the attention of an infant. The study participant is shown a dynamic video that plays continually throughout the calibration process. The video begins in full screen and then scales down to the specific calibration target locations to capture data for the calibration. Between calibration points, the video expands back to full screen.

During calibration, if the participant becomes uninterested, it is possible to switch to a different target video clip ad hoc from the Moderator view without stopping the recording (requires a dual screen setup). Manual calibration type is required.

- *Animation Duration:* Specifies how long it takes the video to transition between the full screen size and the specified target minimum size. Use the drop-down

menu to choose a default speed or enter a custom duration in seconds.

- **Target minimum size:** Specifies the minimum size of the video at the calibration points. This value assumes that the participant is positioned perfectly centered and 650 mm from the center of the presentation screen. Targets at the corners or edges of the screen will display larger than the center target to ensure that the angular approximation is respected on the overall screen.
- **Target video clip:** This video is presented to the participant throughout the calibration. To add a custom target video clip from your local computer (.avi or .mp4), click the plus sign (+) next to **Add video**. The video will be shown full screen so ensure proper resolution.
To download an example video clip (caterpillar), click the **Download video** icon (requires an internet connection). The suggested video will open in your default browser. Download it to your computer and then browse to it using the **Add video** button as above.

6.2.4 Image stimulus

The following properties of an image stimulus can be changed:

- **Name:** By default, the stimulus name is identical to the file name of the image file. You can change it to any desired name. However, each stimulus must have a unique name.
- **Scaling:** This determines how the stimulus appears on the screen.
 - If you select *Original*, the size (in pixels) of the image file is retained and displayed as is. This means that, if the resolution of the screen is lower than the image size, only parts of the image are shown.
 - If you select *Fit to screen*, the image is either enlarged or shrunk (keeping its aspect ratio) to fit on the screen during stimuli presentation.
- **Position:** This governs where on the screen the stimulus appears (e.g. centered or upper-left corner). Change the display position by clicking on one of the position boxes. The active position is represented by a turquoise square.
- **Background color:** Most stimulus images don't cover the entire screen. This setting determines the color of the unused part of the screen.
 1. Select the background color by opening the drop-down and the color picker.
 2. Click on the desired color or enter a hexadecimal color code in the number entry field.
 3. The background color is updated and appears in the Timeline as the background of the stimulus card for the selected stimulus.
- **Show mouse cursor:** Click to record and show the mouse cursor.
- **Advance on** determines for how long a stimulus is displayed. This property has three settings that can be used one by one or in combination. For example, if both *Time* and *Key press* are selected, the next stimuli gets displayed when the time has run out or if any key has been pressed before the expiration of the time set in the Time setting. The available settings are:

- **Time:** You specify for how long a stimulus is displayed, using seconds and milliseconds. The actual viewing time can vary dependent on the update frequency of your screen and the capabilities of your computer. The setting has a minimum value of 50ms; This value is selected based on tests showing that even high-end computers cannot guarantee consistent exposure times at lower values.
- **Key press:** If this option is selected, you advance from the one stimulus to the next on the Timeline by pressing any key on the keyboard.
- **Mouse click:** If this option is selected, you advance from the one stimulus to the next on the Timeline by clicking the mouse.
- **Gaze:** The stimulus is displayed until the participant fixated a certain amount of time within a defined area (Trigger zone) on the stimulus. The time is specified in seconds and milliseconds. If there are gaze samples outside the Trigger zone, or there is >34ms of data loss before the trigger time is reached, the Gaze trigger time will start over from zero.
- **Data loss reset:** The default value for data loss is 34 ms but you can adjust the time using the up and down arrows.

6.2.5 Video stimulus

The following properties of a video stimulus can be changed:

- **Name:** By default, the stimulus name is identical to the file name of the video file. You can change it to any desired name. However, each stimulus must have a unique name.
- **Scaling:** Determines how the stimulus appears on the screen.
 - If you select *Original*, the size (in pixels) of the video file is retained and displayed as is. This means that, if the resolution of the screen is lower than the video size, only parts of the video are shown.
 - If you select *Fit to screen*, the video is either enlarged or shrunk (keeping its aspect ratio) to fit on the screen during stimuli presentation.
- **Position:** Governs where on the screen the stimulus appears (e.g. centered or upper-left corner). Change the display position by clicking on one of the position boxes. The active position is represented by a turquoise square.
- **Background color:** Most stimulus videos don't cover the entire screen. This setting determines the color of the unused part of the screen.
 1. Select the background color by opening the drop-down and open the color picker.
 2. Click on the desired color or enter a hexadecimal color code in the number entry field.
 3. The background color is updated and appears in the Timeline as the background of the stimulus card for the selected stimulus.
- **Advance on** determines for how long a stimulus is displayed. This property has two settings that can be used one by one or in combination. For example, if both *Key press* and *Mouse click* are selected, the next stimuli gets displayed when any key has been pressed or if any mouse button has been clicked. However, when the clip has reached its end, the

next stimuli on the Timeline appears automatically. The available settings are:

- *Key press*: If this option is selected, you advance from one stimulus to the next on the Timeline by pressing any key on the keyboard.
- *Mouse click*: If this option is selected, you advance from the one stimulus to the next on the Timeline by clicking the mouse.
- *Gaze*: The stimulus is displayed until the participant fixated a certain amount of time within a defined area (Trigger zone) on the stimulus. The time is specified in seconds and milliseconds. If there are gaze samples outside the Trigger zone, or there is >34ms of data loss, before the trigger time is reached the Gaze trigger time will start over from zero.
- *Data loss reset*: The default value for data loss is 34 ms but you can adjust the time using the up and down arrows.

6.2.6 Stimulus editor (images and video)

To view and edit an image or video stimulus in a larger format, hover the mouse over the stimulus card on the timeline and click the **Edit** button (pencil) in the lower right corner. The selected image is then displayed in a larger format. If **Advance on Gaze** is selected you can edit the size and placement of the trigger zone in this view. For example, you may want to advance when the participant looks at one part of the screen or at a particular object. The default zone is in the middle of the stimulus.

To navigate back to the Timeline view, use the **Back to Timeline** button in the upper left corner.



Sometimes stray coordinates end up outside the trigger zone and these can then restart the trigger time calculation. Tobii Pro recommends making the trigger zone as large as possible so that any stray coordinates are more likely to be inside the gaze trigger zone.



Note that any modifications to the Trigger zone will be reset if the option is deselected and then reselected.

6.2.7 Web stimulus

The Web stimulus is used to display web pages to participants during a recording. In order to use this element during a recording, your computer must be connected to the Internet or the websites need to be stored locally on your computer. Tobii Pro Lab uses the website's URL or a local address to open the site in the built-in Lab browser. This browser automatically records all mouse clicks, keystrokes, and web pages visited during the experiment. The following properties of a Web stimulus recording can be changed when you select a web element on the Timeline on the **Design** tab:

Name: By default, the stimulus name will be Web page (n). You can change it to the name of the web page, however, each stimulus must have a unique name.

URL: Enter a URL or a local computer address (you do not need to type http://). If you leave it empty, the Lab browser will launch without a URL. For example, this might be useful if the researcher wants the participant to select their own starting web page. The Lab browser launches in full screen mode.



The local files need to be available on the computer you are using. This is because the web stimulus is used to display webpages to participants during a recording. In order to use this element during a recording, your computer needs to be connected to the Internet or the websites need to be stored locally on your computer.

Description: A description of the web stimulus.

Frame rate: Available values are 5, 10, 15, 25 (default), and 30fps (frames per second). In most cases, the default rate of 25 fps is good enough.

Video quality: High means that the screen is recorded at native resolution. Medium is 75% of native and Low is 50% of native resolution. Use high resolution if possible. If you need to save computer CPU or disk space and do not require a detailed recording, select a lower video quality.

Advance on: If the box “Time” is checked, you can then set the amount of time before the recording moves to the next stimulus on the Timeline.

Web screenshots: Toggle on and off automatic screenshots. Alternatively, you can capture screenshots manually instead using the F8 key or clicking at the **Capture web screenshot** button at the bottom of the Moderator screen.

For more information about recording, read [Screen-based recording](#) and for more information about analyzing recordings with web stimuli, read [Replay a recording](#).

6.2.8 Screen Recording stimulus

The following properties of a Screen Recording stimulus can be changed:

- **Name:** Enter an appropriate name. If left as-is, Pro Lab sets the name to Screen Recording with a unique number if there is more than one stimulus with this exact name, e.g. *Screen Recording (1)*.
- **Frame rate:** Available values are 5, 10, 15, 25 (default), and 30 fps (frames per second). In most cases, 25 fps (default) is good enough.
- **Advance on:** If your experiment requires that you advance to the next stimulus after a set time period, check the Time checkbox and enter the length of the period.
- **Video quality:** *High* means that the screen is recorded at native resolution. *Medium* is 75% of native resolution and *Low* is 50% of native resolution. Use native resolution if possible. However, if you need to save computer CPU or disk space and you are not interested in details in the recording, select a lower video quality.

For more information about analyzing a Screen Recording, read [Replay a recording](#).

6.2.9 Text stimulus (Reading studies and instruction)

Use the Text stimulus to write text you want the participant to read during the recording. This type of stimulus is intended for reading research but you can also use it to give instructions or anything else conveyed with text.

Tobii Pro Lab will compute AOI reading metrics on Automatic AOIs, generated on three levels (depending on what you choose): Characters, Words and Sentences.




Figure 4. Sample instructions using the Text editor

Click the **Edit text** button (or double-click the element on the Timeline) to write the text and to see the full editing interface with additional functionality including a font selector, font size, and margin control. You can also choose to visualize boundaries to see the Automatic AOIs' size.

The following Text stimulus properties can be changed:

Name: Enter an appropriate name. If left empty, Pro Lab sets the name to Texts with a unique number if there is more than one element with this name, e.g. Text (1).

Background color: Sets the background color of the entire Text element.

Text color: Sets the font color

Margins (%): Controls the amount of space between the edge of the element and the text.

Display margins: Click the Display Margins button in the full editor to see the margin's dotted lines as you adjust the percentage. These dotted lines won't show in the element.

Edit text Click the button to see the full editing interface with additional functionality including a font selector, font size, margin control and boundaries. Note that the line height will affect the vertical size of AOIs.

Show mouse cursor: Select this to enable the participant to see their mouse cursor during the recording.

Advance on: Determines what triggers the next stimuli to display. The available settings are:

- **Time:** Time: The minutes and seconds clock appears on the Properties panel when this is checked. Participant advances to the next stimulus after the set amount of time.

- **Key press:** Participant advances to the next stimulus by pressing any key on the keyboard.
- **Mouse click:** Participant advances to the next stimulus by clicking the mouse.

Punctuation

Pro Lab supports the symbols using a question mark (?); exclamation mark (!); and period (.) as an end to a sentence.



The punctuation marks must be followed by a space so that they indicate the end of a sentence.

Text Editor

- The text editor supports Unicode characters in the Basic Multilingual Plane (BMP). Among other things, in the BMP, there are the Latin characters and symbols, transcriptions, other European characters and writing systems. For more information, search for unicode information online.
- The text editor does not support Input Method Editors (IME). Using an IME might cause instability and cause Pro Lab to stop working.
- The text editor supports copy and pasting from external text editors as long as Unicode characters in the Basic Multilingual Plane (BMP) are used. Pay special attention to the font properties when copying and pasting because these are not always copied correctly.

6.3 Groups

Groups are containers of timeline elements and enable more compact timeline designs. The behavior of the elements inside a group is determined by the action property set for the group. You can nest groups within groups to create more complex stimuli presentation behaviors.

6.3.1 Edit a group's contents

How to edit the group properties and content:

1. There are 3 ways to open a group:
 - Click the triangle at the bottom right-hand corner of group's tile; or
 - Double-click the group; or
 - Highlight the group and press Enter on your keyboard.
2. The group's content appears below the group element and displays its contents.
3. Add a new stimulus or group:
 - by dragging-and-dropping; or
 - by clicking-and-browsing in the empty tile in the Timeline; or
 - by dropping onto the unexpanded group icon.
 - If you wish to add a stimulus or to group them with another group, you must select the target group (click in the expanded area) and open it before doing a drag-and-drop.
4. Change the order of the group elements on the Timeline by dragging-and-dropping its Stim-

uli and Groups in the desired order.

5. The addition is saved automatically. You can go on working with Pro Lab.



If you delete a stimulus in a group, your action only deletes that copy of it, whereas all other copies elsewhere among groups and timelines remain unaffected. However, if you change the properties of a stimulus in a group, your action affects that stimulus everywhere among groups and timelines, also in such instances that have been copy-pasted.

6.3.2 Duplicate a group

How to duplicate a group:

1. Duplicate a group by selecting it and copying it (CTRL + C).
 - All stimuli and groups nested inside it are duplicated, too.
 - This is an appropriate method when you want to nest a group inside another group.
2. Open the target group and select the item just to the left of the desired location.
3. Paste (CTRL + V) the copy. It appears in the desired location.



You can change the Timeline inside a duplicated group without affecting the original version, and vice versa. However, if you change the properties of a stimulus in a group, your action affects that stimulus everywhere among groups and timelines, also in such instances that have been copy-pasted.

6.3.3 Delete a group

How to delete a group:

1. Delete a group by selecting it and clicking on the garbage can icon, or by pressing the delete button on the keyboard.
2. This action also deletes all groups and stimuli nested inside it.

6.4 Stimulus/group variables

A stimulus/group variable is a specific feature of your stimulus or group of stimuli that are part of the context in which the behavior occurs. It is often an expression of or a subset of your independent variable and covariates. Examples include number of items, item category, stimuli density, color, brightness, and contrast.

In Tobii Pro Lab, stimulus/group variables are a way of adding metadata related to the stimuli that you can use later during your analysis, e.g. independent variables or co-variate that you wish to add as a column in your data export.

6.4.1 Create variables

How to create a stimulus/group variable

1. In Pro Lab, open your project in the *Project Overview*.
2. Open the *Stimulus* panel by clicking on “Stimulus Variables” in the left column.

3. Create a new variable by clicking on the plus sign in the top right corner of the *Stimulus* panel. A new variable appears with a default value.
4. Enter the variable's name in its name field to the right of the chevron (single "arrowhead").
5. Press Enter after writing the name. The cursor moves to the first value field.
6. Enter the value's name in its name field. Add a new value by pressing Enter or by clicking on the "Add value" plus sign.

6.4.2 Assign variables

How to assign a stimulus/group variable:

1. Open the *Design* module.
2. Add stimuli by using the **Add** buttons located in the upper right and then browsing and selecting them on your computer.
3. Select the desired stimulus or group. Multiselect is enabled.
4. Assign a value by clicking the checkbox next to it in Stimulus variables on the Properties panel. Use the **Edit variables** button to add or edit variables.
5. Select the value you want to assign to the stimulus or group.

6.5 Use TTL markers to indicate stimulus onset for external listeners

TTL signals are the current standard for synchronizing data between stimuli presentation software, biometric data sources (e.g., EEG, GSR), and data collection software. In Pro Lab, you can send TTL signals through a parallel port on the computer running the software.

Pro Lab sends a TTL marker only for an *image* or *video* stimulus onset.



If this setting is activated, TTL markers are sent for image and video stimuli on the Timeline, unless manually specified otherwise. In addition, settings for marker value type and marker bit depth are also set for the stimuli collectively. Only the TTL marker value can be specified per stimulus, provided that you set this property to manual.



Pro Lab sends the TTL markers when it estimates that the image or video stimulus has appeared on the screen. This estimation uses, unless otherwise specified, a zero lag display and a 10ms screen latency. If you need to have high timing precision when sending TTL markers at stimulus onsets, read the [Timing guide for stimulus display in Pro Lab](#) to learn how to provide your setup with a precise stimulus timing.

How to activate the sending of TTL markers on stimulus onset:

1. Select an image or video stimulus on the Timeline.
2. Click the "Send Stimulus onset markers (TTL)" toggle switch under the Advanced settings in the right-hand panel.
3. Select Automatic or Manual as the marker value. The default setting is Automatic. If a stimulus appears several times on the same Timeline, all occasions get the same marker value.



If you select “Automatic,” Pro Lab automatically assigns a unique value to each stimulus on the Timeline, starting with 1 for the first and going up incrementally by 1 for each following stimulus. If you select “Manual,” Pro Lab still assigns a value to each stimulus, but you can modify that value in the properties of each stimulus.

How to change the TTL marker value for a stimulus:

1. Click the image or video stimulus on the Timeline.
2. Enter the marker value you want to use for that stimulus in the TTL marker value input field.
3. If you don't want to send a TTL marker for a specific stimulus, set its TTL marker value to 0.
4. Select the desired marker bit depth. The default is 8 bits. This value governs how many marker values you can use and it depends on how much the listener to the TTL signal can handle. An 8-bit marker depth enables the use of up to 255 different marker values, whereas a bit depth of 3 or 1 enables 7 or 1 different markers.



The marker signal is set to be active for approximately 34ms after the moment when Pro Lab assumes the stimuli starts to when it becomes visible on the presentation screen. This value was selected to allow for sampling frequencies of 30 Hz in devices listening to the signal while still being shorter than the shortest stimulus display time (i.e. 50ms).

6.6 Change the appearance of the Timeline stimulus cards

Each stimulus instance is represented by a Timeline card. It contains information about the stimulus name, what it looks like, and what its background color is. These settings are governed by the stimulus properties. However, a card's shape and size are determined by global settings that affect all cards at once.

6.6.1 Change the size of the stimulus cards

In the bottom left corner of the Design module interface, there is a toggle switch that enables you to display large or small representations of the stimulus cards. Change card size by clicking on the left part of the switch for small, or on the right part for large. The default setting is small cards.

6.6.2 Emulate the screen resolution of the presentation screen

You use the drop-down list “Preview screen resolution” to select what resolution is used on your presentation screen. The shape of the stimuli cards gets the same aspect ratio as the screen on which the stimuli are displayed. In addition, the stimuli are scaled on the cards to indicate how much screen area they cover during the presentation.

All Tobii Pro eye trackers delivered with screens get their recommended screen resolutions listed in the drop-down. So if you know that you will use, for example, a Tobii Pro Spectrum, you can select it in the list and see on the stimuli cards how the stimuli appear on the Pro Spectrum screen.

6.7 Correct resolution and scaling of the presentation screen

Each stimulus instance is represented by a Timeline card. It contains information about the stimulus name, what it looks like, and what its background color is. These settings are governed by the stimulus properties. However, a stimulus shape and size is determined by global settings for resolution and scaling that will also determine the shape and size of the timeline element.

6.7.1 Set the screen resolution of the presentation screen

Use the drop-down list “Presentation screen resolution” to select what resolution should be used on your intended presentation screen. The shape of the stimuli cards gets the same aspect ratio as the screen on which the stimuli will be displayed. In addition, the stimuli are scaled on the cards to indicate how much screen area they cover during the presentation. When the recording starts the chosen resolution will be applied on the presentation screen and then reverted back if the screen had a different resolution.

All Tobii Pro eye trackers delivered with screens get their recommended screen resolutions listed in the drop-down. So if you know that you will use, for example, a Tobii Pro Spectrum, you can select it in the list and see on the stimuli cards how the stimuli appear on the Pro Spectrum screen.



Always use the native resolution of the presentation screen to guarantee correct timing and correct aspect ratio.

6.8 Record in 100% scale

In order to start the recording, the scaling of the presentation screen must be 100% in order to ensure correct format of the stimuli. To change the scaling, adjust your computer’s Display Settings.

6.9 Timeline zoom

Sometimes you need to see the cards in your Timeline to be smaller or larger. In the top right corner of the Design module (next to the Recorder link button) you can switch between viewing large or small representations of the stimulus cards.

Change card size by clicking on the two symbols. The default setting is small cards.

7 Make recordings

The Recorder module in Pro Lab provides you with all controls required for making recordings for later analysis in the Analyze module.



If you work with a Glasses project, the Record module is not available because there is no need to access the Record module for Glasses projects.

7.1 Mouse Tracker

If you want to build and pilot a data project but you don't currently have access to an eye tracker, you can use the Mouse Tracker as a temporary substitute. This sensor produces gaze data at the location of the mouse cursor. Note that filters and analysis may not work like they do for real eye-tracking data therefore this feature is intended for piloting and demonstration and teaching purposes, not bona fide research. The data is sampled at 60 Hz or 250 Hz.

How to use the Mouse Tracker:

1. In the Record module, select "Mouse Tracker" from the list of eye trackers in the Eye Tracker card above the timeline.
2. Use the radio buttons to select between 60 Hz and 250 Hz.

7.2 Screen-based recording

If you only have one Timeline and you approve of your selections in Pro Lab, start a recording by selecting the **Start recording** button in the lower-right corner of the **Record** module. However, if you want to verify all selections, use the instructions below to start a recording.

A recording ends after the appearance of the last stimulus or by an action by the moderator.



Pro Lab does not support presentation screens that have the Windows screen duplication setting enabled

7.2.1 Start a screen-based recording

How to make a recording:

1. Ensure that the desired devices are enabled and selected. They (eg. the eye tracker, participant camera, microphone, GSR, and TTL port) are represented by cards above the Timeline(s). For more information, read [Change device settings](#).
2. Select the desired Timeline by clicking it. The active Timeline is highlighted in purple. The number before the name of the Timeline indicates how many recordings have previously been made using it.
3. Select the participant you will record from the list of participant names. If that person has no previous recordings, select the **+New** button and enter the participant's name in the name field. Note that the **Start Recording** button is disabled until a participant is selected.



You can define or edit already defined participant variables. For more information, read [Participant Variables](#).

4. Give the recording a name by entering it in the “Recording name” input field in the Recording section.
5. Specify which of your connected screens will display the stimuli. The active screen is highlighted in purple in the Target Screen section of the interface. If you are unsure which screen is which, click on the “Identify” button. A number corresponding to the representations in Pro Lab then appears for a moment on all connected screens. You can also modify screen latency. For more information, read [Modify screen latency settings](#).
6. Start the recording by clicking the **Start recording** button. If the Start recording button is grayed out after doing the above steps, you may need to prepare the timeline. For more information, read [Prepare timeline](#).

How to stop the recording:

The recording can be stopped using the ESC key (default method) or Shift+ESC.



The ESC key is sometimes custom-programmed on individual computers to close the browser or another application so it is recommended to use Shift+ESC for Web recordings.

7.2.2 Prepare timeline

The first time you create a timeline with a text element and want to record it, or if a previously created timeline with a text element has changed, then you will see the **Prepare timeline** button in the lower right instead of the **Start recording** button.

How to prepare the timeline:

1. Click the **Prepare timeline** button.
2. Pro Lab will prepare the timeline for recording. When it is finished, select **OK**.
3. Follow the steps in [Start a screen-based recording](#).

7.2.3 Moderate a screen-based recording

When a recording is running, the Pro Lab application window switches to the moderator view. It displays those stimuli that the participant sees on their screen and it additionally informs the recording moderator about whether the eye tracker is detecting a gaze, whether the participant is positioned in a way that enables data collection, for how long the recording has been running, and what participant name and recording name are used. It also provides the moderator with the opportunity to terminate the recording at any time.



The moderator view is only available if more than one screen is connected to the Pro Lab computer and the Target Screen is not the same as the screen on which the application window is shown.

7.2.4 Perform a screen-based calibration

Calibration

In human populations, there is a natural variation in the shape and geometry of the eyes. For example, the exact location of the fovea varies from individual to individual. Tobii Pro Eye Trackers, use a subject calibration procedure to optimize its gaze estimation algorithms (i.e. the 3D eye model) and account for this variation, resulting in a fully customized and accurate gaze point calculation. Additionally, some eye trackers use the calibration procedure to select the detection mode (dark or bright pupil) that provides the best data, and then lock, or use that mode predominantly during the recording.

During the calibration the participant is instructed to look at calibration targets that appear at multiple locations on the plane (commonly the surface of the display monitor) where the stimulus is located. The data collected during this period is then mapped to those locations using, a standard configuration of the 3D eye model or the last configuration saved to the eye tracker, e.g. the last previous test participant. The 3D eye model in Tobii Pro Eye Trackers compensates for drift, so the calibration is done only once before starting to collect data and does not need to be adjusted during recording. However, keep in mind that fluctuations in pupil size may cause shifts in the indicated gaze position, and the experiment should be designed accordingly. It is good practice to use a calibration background color that matches the stimuli.

Add a calibration stimulus to the Timeline in the **Design** module. For more information, read [Timelines](#). When that stimulus is reached on the Timeline, the calibration procedure is initiated.

Validation

In order to trust the quality of your recorded data, it is important to verify and quantify the eye-tracking data quality during the recording. Consequently, to evaluate the performance of this new algorithm configuration it is best practice to test it with a new data set. This provides a more realistic estimate of the data quality of your recording.

This process of collecting new data to quantify and evaluate the data quality after calibration is called a calibration validation.

Calibration validation is on by default. This means that it occurs directly and seamlessly after the calibration is done. During the validation, the participant sees four additional calibration targets on the screen.

The results of the validation are displayed after the calibration validation procedure. They can be reviewed and exported later, in the **Project Overview** section, or by using the Data Export tool. They are highly dependent on the calibration used at the time you do the validation. Any changes made to the calibration will have to be validated again.

Disable validation

Calibration validation is a part of the normal calibration process. But sometimes it is necessary to disable the validation when it is hard to hold the subject's attention, for example, when the participant is an infant. You can turn off the validation by selecting the Calibration element in your Timeline and then toggling the **Validate Calibration** switch in the Properties Panel, under Calibration.

The controls for the calibration can either be displayed on the target screen in front of the participant or in the moderator view. This setting is made on the Record section, under the

Calibration controls list box. If you display the controls in the moderator view, you get better control of the calibration.

Perform a calibration validation:

1. Select a Timeline that contains a **Calibration** stimulus, and start a **Start a screen-based recording**.
2. The participant should sit comfortably in front of the eye tracker in a way that provides the eye tracker with an unrestricted view of the participant's eyes. Do this by using the graphical representation of the track box (i.e. the location of the eyes in relation to the eye tracker) in the Moderator view or the Presentation screen. If the eye tracker can see the eyes, it depicts them as two white dots that ideally appear in the center of the track box.
3. The participant should let their gaze follow the target object that appears on the screen once the calibration starts.
4. Click the **Start calibration** button on the presentation screen. The calibration procedure starts. If requested, the validation starts immediately after the calibration and presents four additional targets.
5. When the calibration data has been collected and validated, a visual representation of the calibration results displays on the presentation or moderator screen.
6. **Inspect the calibration and validation results**. When the calibration and validation results are satisfactory, click **Continue**. The rest of the Timeline is presented to the participant.

Inspect the calibration and validation results

The visual results can be displayed either as averaged values for both eyes combined or separately for each eye. Select this using the two radio buttons at the top right of the Calibration window: Average Results or Left and Right Eyes Results.

Each colored dot represents a gaze sample collected during the calibration or validation, green for the left eye and blue for the right eye. Orange dots represent the average gaze sample between left eye and right eye. White dots represent the average gaze sample across entire targets. Colored lines are drawn between targets and white dots to represent accuracy. Colored circles are drawn around white dots to represent precision.

The crosses indicate where the validation target objects appeared on the screen and the pluses where the calibration target objects appeared on the screen. When your cursor hovers over a target, a pop-up window appears with the target's accuracy and precision displayed in degrees and pixels as well as the number of related gaze sample.



Tobii Pro eye trackers display different numbers of gaze samples for calibration targets. Some eye trackers will display only one gaze sample per calibration target.

To the right of the window, the **Calibration average results** table helps you interpret the results. This is the average total results across all calibration targets for all the calibrations done with the current participant. These results represent the estimated error in the gaze estimation algorithm (3D eye model). They do not represent the real accuracy and precision you can expect during the recording. Instead, if the calibration you just did was validated, it will show validation average

results below it. Validation average results are the estimate of the average total results across all validation targets for the relevant calibration. These results represent an estimate of the accuracy and precision you can expect during the recording.

If data is missing from some calibration targets, if data is mapped too far away from the targets or if the calibration average results are too high for your experiment, you should recalibrate. When recalibrating, you can choose to either redo the entire calibration or to recalibrate only selected points. To redo the whole calibration, click the **Recalibrate all points** button. If you want to redo the calibration for selected points, select them by clicking on them and then click **Recalibrate**.

If data for the calibration is good enough, you should inspect the data for the validation. Ultimately, the validation data is what matters if you want good accuracy and precision during the recording. If data is missing from some validation targets, if data is mapped too far away from the targets, or if the validation average results are too high for your experiment, you should validate the calibration again. If you want to validate the calibration again, you cannot select individual points but instead must run the validation again. Select **Revalidate**.

More information about accuracy and precision definitions can be found in our Metrics report on tobii.com. Please note that the precision values are provided without applying a filter.

Reuse calibration

In some circumstances you can reuse previously made calibrations (e.g. you are running multiple recordings in a row with the same participant). It is recommended to reuse calibrations only if conditions are similar from one recording to the next one: same participant, same experimental environment, same lighting conditions, similar stimuli background, etc. For each Participant, Tobii Pro Lab saves the last five calibrations.

How to reuse a calibration:

1. Start a **Screen-based recording** with an existing Participant that was previously calibrated.
2. When the recording starts, you see the track box and click on Reuse calibration.
3. The calibration results window opens. Select the calibration to use from the Calibration average results table and click Continue or Validate (if validation is active).



Reused calibrations will have to be re-validated if validation is active.

7.2.5 End a screen-based recording

A recording can either end when the last stimulus has been displayed or by an action of the moderator. To stop an ongoing recording, the moderator can:

- press F10 on the keyboard or click the **Next** button in Moderator view
- press the ESC button or press Shift + ESC on the keyboard
- click the **Stop recording** button in the Moderator view

When the recording finishes, you can Save or Discard the recording. To save it, you must confirm that you have obtained participant consent. Alternately, you can select *Discard recording* to delete the recording. The participant will also see the dialog about their consent to store/transfer

their eye tracking data on the target screen. For more information, go to [Tobii's Data Transparency policy](#).

A recording either ends when the last stimulus has been displayed or by the moderator's interruption.

7.2.6 Lab browser

Dynamic web content

Dynamic web content is a conceptually and technically difficult aspect of today's web pages.

Concept-wise, it is difficult to compare data from several recordings when the dynamic content makes each interaction with a web page different for each participant.

Technology-wise, it is challenging to capture content that changes over time in a useful way to be used in visualizations of eye-tracking data, especially if the content only changes on some parts of a page. Pro Lab captures a screenshot of each page the participant visits and whatever is shown on the page at that time will be what is visible in the Visualizations tab.



You can capture screenshots manually using the F8 key; or by clicking the Capture web screenshot button at the bottom of the Moderator tool; or adding your own screenshot in the Custom Web TOI dialog.

Screenshots: The Lab browser will capture an automatic screenshot for every viewport (the part of the page which is visible on the screen) and then stitch them together into a full screenshot. The screenshot is captured before leaving the page. Any fixed elements on the page will be displayed only once and in the place where they first appeared. Fixations on a fixed element will be transferred to the right place even if the page was scrolled.



Pro Lab might not be able to take a screenshot of Single-page applications or pages. (A Single-page application (SPA) is a web application or web site that dynamically rewrites the current page rather than loading completely new pages from a server.)

On the Design tab, screenshot capturing can be turned off to save on performance during recording. Screenshots can then be captured manually using the F8 key or clicking the Capture web screenshot button at the bottom of the Moderator tool.

You can also import a screenshot in the Custom TOI dialog.

Move to Next and Stop: When pressing F10 for next stimuli or Esc/Stop for ending the recording; it will take a few moments before the browser has reloaded or closed down. Make sure to not press the button again.

Tabs: The Lab browser does not display tabs. It will not display web pages prompted to be opened in a new tab or window. Participants can use Ctrl + click to open in the same window.

Full screen mode: The web page can only be displayed in full screen mode.

Download: Downloading is disabled although participants can still click the download link.

Bookmark: It is not possible to add bookmarks. This is because there is a new session of the web browser every time so customizations such as bookmarking are not possible.

Video content: Flash player is not, and cannot be installed. this may cause some video content to not be displayed.

Hamburger menu: The hamburger menu in the upper right corner is visible but has no function.

7.3 Scene camera recording

When you have created a new scene camera project, or opened an existing one, the Project Overview window appears on your screen. Open the Record module by clicking the Record option in the main menu bar.

A recording is terminated by an action of the moderator.

7.3.1 Start a scene camera recording

How to make a recording:

1. Check that the devices used for the recording are enabled and selected. They (e.g. the eye tracker, scene camera and TTL out port) are represented by cards at the top of the window. Note that there is no Timeline(s) in Scene Camera projects. For more information, read [Change device settings](#).



The Scene Camera card is only displayed when a Scene Camera is connected to the system.

2. Deactivate the scene camera's Autofocus function. See its user manual for more information. If the Autofocus function is active during a recording, focus may shift and cause incorrect data mapping.
3. Sound is recorded by a connected microphone by default. If you don't want to record sound, set this on the Scene Camera card by clicking on the card and disabling the "Use audio" switch. If more than one microphone are connected to your system, select which one to use in your recording.
4. Select the name that you will record from the list of Participant names. If that person has no previous recordings, select the **+ New** button and enter the participant's name in the name field. If no name is entered, Pro Lab suggests a name, such as "Participant1." Note that the **Start Recording** button is disabled until a participant is selected.



You can define or edit already-defined participant variables. For more information, read [Participant Variables](#).

5. Enter the recording's name in the "Recording name" input field in the Recording section. If no name is entered, Pro Lab suggests one, such as "Recording1."
6. Start the recording by clicking the "Start recording" button.
7. When the recording finishes, you can Save or Discard the recording. To save it, you must confirm that you have obtained participant consent. Alternately, you can select *Discard recording* to delete the recording. The participant will also see the dialog about their consent to store/transfer their eye tracking data on the moderator screen. For more information, read the [Tobii Eye Tracking Data Transparency Policy](#).



When you use an eye tracker that is not attached to a monitor, it must be configured so that the application registers where the eye tracker is located in relation to the stimulus and what area is supposed to be measured. The eye tracker should be configured for stand-alone use using Tobii's free Eye Tracker Manager software. There is also configuration information included in your eye tracker's user manual.



If you disregard configuring the eye tracker, you will get incorrectly mapped data.

7.3.2 Moderate a scene camera recording

When a recording starts, the Pro Lab application window switches to the moderator view. It displays the video of the scene camera and informs the moderator whether the eye tracker is detecting gaze, whether the participant is positioned in a way that enables data collection, for how long the recording has been running, and what participant name and recording name are used. It also provides the moderator with a way of terminating the recording at any time.

7.3.3 Perform a scene camera-based calibration

You ensure good eye-tracking data by performing a calibration before starting data collection for analysis. When you use a scene camera, you can calibrate in two ways: either with a pre-made calibration board, such as the one included with Mobile Device Stand accessory (or with one you make by yourself) *or* by using any physical object, such as a mobile phone or table.




Figure 5. An example of a calibration target

The calibration controls are displayed in the Record module.

How to perform a calibration:

1. Put the calibration board or a physical object in front of the scene camera at the required distance.
 - a. Set the active area that you want to measure by moving the four purple markers to the corners of your requested area.
 - b. Add the desired calibration points by dragging them from the small box in the top right corner of the window. You must use at least two calibration points, but your calibration result improves with more points. We recommend five points: four about 25% inwards from the four corners and one in the middle.

- c. If you use a physical object, select points that can be easily defined as targets for the participant.
 - d. When you are done, lock the calibration plane by clicking the pad-lock symbol on the right side of the window.
2. Start the calibration by clicking the **Calibrate** button. Two buttons appear: **Calibrate point** and **Abort calibration**.
 3. The first calibration point is indicated. Instruct the participant to look at it. When they do so, click **Calibrate point**. Repeat this step for each calibration point.
 4. You can stop the calibration by clicking **Abort calibration**.

7.3.4 End a scene camera recording

A recording can either end when the last stimulus has been displayed or by an action of the moderator. You stop a recording before the last stimulus has been completed by clicking the “Stop recording” button in the moderator view or by pressing the “Esc” button on the keyboard.

When the recording finishes, you may Save or Discard the recording. To save it, you must confirm that you have obtained participant consent. Alternately, you can select *Discard recording* to delete the recording.

The participant will also see the dialog about their consent to store/transfer their eye tracking data on the Moderator screen. For more information, read [Tobii's Data Transparency policy](#).



A Scene Camera recording can only be stopped by the moderator, as the recording doesn't follow a timeline.


7.3.5 Participant Variables

Participant Variables are used to filter data when generating visualizations, calculating eye tracking metrics and comparing the behavior of different participant groups. Variable and variable values are also included in the data export files for further statistical processing. Once created, you can select values for the participant variables when you edit your participants.

The participant variables are handled from the Participant card in the Record module or on the Project Overview. You can define participant variables before or after recordings.

Add a new participant variable

You can add, edit, or remove participant variables from either the Project Overview or Record panels:

1. In the Record module, click the **Edit variables**  button. (If you are on the **Project Overview** tab, the Participant Variables button is on the left-hand side.) Click the button to display the Participant variables panel on the right-hand side of the screen. Use the panel to add, edit, or remove participant variables.
2. If this is the first time you are adding participant variables, use the Plus “+” button to add a new participant variable and a generic participant variable. If you already have participant variables defined for this participant, they are displayed in the **Participant variables** panel.

For more information, read [Edit a participant variable](#).

3. Enter the name of the first participant variable in the text field and press Enter.
4. Add a new value for the participant variable by clicking the “+” button. A new participant variable value is added.
 - a. Enter the name of the value.
 - b. Add the wanted values by using the Add Value button.
 - c. If the participant variable is supposed to accept multiple values, select the “Allow multi selection” check box.
 - d. You can also add a new value by pressing Enter in a value field.
5. When you have added the desired values, you can add another participant variable by repeating steps 2-4d.

Add a new participant

How to add a new participant:

1. Click the Plus “+” button and enter the participant’s name.
2. If there are already defined participant variables, select the values for the participants.

Edit a participant variable

How to edit an existing participant variable in Project Overview:

1. If you need to expand the desired participant variable, click on the down arrow below it or click on the “Expand all” button to display all participant variables.
2. Edit the desired the data. When you are done, press Enter.
3. If you want to create new participant variables and values, see [Add a new participant variable](#).

Options of deletion:

- If you want to delete a value for a participant variable, click the trash can symbol to the right of the value.
- If you click the trash can symbol to the right of the participant variable name, you delete that participant variable and all of its values. A dialog box lists the affected participants.

Delete a participant

How to delete a participant:

1. Find the desired participant and click the trash can symbol to the right of the participant name.
2. If there are existing recordings for this participant, Pro Lab lists those that are about to get deleted. **If you want to keep these recordings, DON’T delete the participant.**

7.4 External Presenter recording

An external presenter is a piece of software, for example E-Prime or other third-party software, that communicates with your Tobii Pro Lab software. This section describes how you install the relevant software and connect the eye tracker to the test computer, to Tobii Pro Lab, and to E-Prime.



For more information, read [Analysis of an External Presenter project recorded by E-Prime](#).

The E-Prime Eye Tracking Extensions for Tobii Pro enable you to combine E-Prime's experiment design and stimulus presentation with Tobii Pro Lab's visualizations and AOI analysis. The extensions enable E-Prime to communicate with the Tobii Eye Tracker (TET package) and the Tobii Pro Lab software through a set of instructions (Routines). These instructions are organized in two E-Prime packages:

- The TET Package contains instructions that interact directly with the eye tracker.
- The Tobii Pro Lab Package (TPL) contains a collection of instructions that interact with Tobii Pro Lab. For example, it instructs Tobii Pro Lab to start recording.

7.4.1 Connect the computer and Tobii Pro eye tracker

How to connect your computer and the Tobii Pro eye tracker:

1. Connect peripherals (keyboard, mouse and other peripherals that you are supposed to use during your eye-tracking experiments) and power cable to the test computer.
2. Assemble the eye tracker (if necessary) and connect its power cable if there is one. For more information, check the user manual for your Tobii Pro eye tracker.
3. Connect the LAN cable and video cable to the eye tracker and the computer.
4. *Optional*: Connect the secondary screen to the computer. This requires a dual output graphic card or a splitter.
5. Switch on your computer and your eye tracker.
 - If no image appears on the eye tracker's screen, click the Source button located in the front of your eye tracker.

7.4.2 Configure the computer and install the software

How to configure and install the software:

1. Install Tobii Pro Lab, E-Prime and E-Prime Extensions for Tobii Pro on the test computer. For more information, read [Setting up an External presenter recording with E-Prime](#) in [Tobii Pro Learn & Support](#).
2. Connect the E-Prime dongle to one of the USB ports of the computer and install the dongle driver. Windows detects and installs the correct driver automatically. If Windows doesn't to find the correct driver, check this support page for advice (login is required): support.pstnet.com/hc/en-us
3. Reboot your computer.

7.4.3 Connect Tobii Pro Lab to the eye tracker

How to connect Tobii Pro Lab to the eye tracker:

1. Launch Tobii Pro Lab on your computer.
2. Create a new *External Presenter* project and a new test.
3. Go to the *Record* tab and check that the eye tracker is up and running.
4. Switch to E-Prime.

7.4.4 Set up an experiment in E-Prime Extension and Tobii Pro Lab

How to set up a simple experiment using various package calls to control the eye tracker and Tobii Pro Lab. (For more information about available routines and parameters, check the manual for E-Prime Extensions for Tobii Pro.) Before you start, make sure that E-Prime Extensions for Tobii Pro has been installed.

1. Open E-Studio and create a new *Basic experiment*.
2. Go to the *Experiment Object*, select the *Packages* tab, and add the Tobii Pro Lab package and the TET package
3. Go to the *Devices* tab and add the Tobii Pro Lab device and the Tobii Eye Tracker (TET) device.
 - Leave the default parameters for these two devices as they are, unless you have problems connecting to your eye tracker (for example, if you use several eye trackers on the same network).
4. Start Tobii Pro Eye Tracker Manager and select your eye tracker. Pro Eye Tracker Manager is free Tobii Pro software available on the [Tobii Pro Products webpage](#).
5. Open the *SessionProc*, which is the topmost procedure in the experiment.
6. At the beginning of the *SessionProc*, add the following routines in order:
 - a. *TETOpen* (*comment: opens a connection to the eye tracker*)
 - b. *TETCalibRegular* (*comment: performs a calibration of the user*)
 - c. *TPLOpen* (*comment: opens a connection to Tobii Pro Lab*)
 - d. *TPLStartRecording* (*comment: instructs Tobii Pro Lab to start the recording*)
7. At the end of the *SessionProc*, add the following routines in order:
 - a. *TPLStopRecording* (*comment: instructs Tobii Pro Lab to stop the recording*)
 - b. *TPLClose* (*comment: closes the connection to Tobii Pro Lab*)
 - c. *TETClose* (*comment: close the connection to the eye tracker*)
8. Open the trial-level procedure in the experiment, called *TrialProc*.
9. In the beginning of the *TrialProc*, add the *TPLSetDisplayEvent*. This entails setting a few parameters:
 - *<MediaName>*: Enter the name of the media that you want to use for Tobii Pro Lab. The media name is also the main channel for send information about the conditions applied to this stimulus.

- In this case: we write "Stim_[Text]" as the parameter. This appears in Tobii Pro Lab as "Stim_Cat" if the Text variable in that trial was "Cat". All relevant experiment factors are supposed to be included in the name. Example: If you have Position (left/right) and Animal (cat/dog) as factors and they exist as E-Prime variables, your parameter may be "[Position]_[Animal]_[Image]".
 - <ObjectName>: This is the name of the E-Prime object that you create as a stimulus image and send to Tobii Pro Lab. In this experiment, use the name of the object *Stimulus*, so add "Stimulus" as the parameter.
10. If you have AOIs that correspond to the positions and sizes of E-Prime objects, define the AOIs and send them to Tobii Pro Lab. Add a *TPLSetAOI* routine after the previous routine but before your experiment stimuli. By default this routine has *c* (context) as a parameter, but you are supposed to add more parameters:
 - The media name argument from the *TPLSetDisplayEvent*, to make sure that you put the AOIs on the right image.
 - The name of the AOI, e.g. "AOI1". If the name isn't unique, Tobii Pro Lab changes it to make it unique.
 - The E-Prime object that defines the AOI. Example: An E-Prime Slide object "Stimulus" has a *SlideImage* subobject called "AOIleft", and therefore you are supposed to pass "Stimulus.AOIleft".
 - You may pass more AOI arguments, as long they come in pairs of the name of the AOI and the of corresponding E-Prime object.
 - A full parameter list can look like: *c*, "Stim_[Text]", "AOI1", "Stimulus.AOIleft", "AOI2", "Stimulus.AOIright"
 11. If you intend to have AOI tags, insert the *TPLSetTag* routine and add the following arguments:
 - a. The media name from the two previous routines.
 - b. The pairs of AOI names and tags, for example:
 - "AOI1"
 - "left" or "Animal.cat" with the former being an ungrouped tag, and the latter being a tag of the tag group Animal.
 12. At the end of the *TrialProc*, add a *TPLCompleteLabEvents* routine.
 13. Your simple experiment is now ready for execution.

7.5 Rename a recording

How to rename a recording:

1. Select the desired recording in the **Project Overview**.
2. Rename the recording in one of two ways:
 - Press F2, *or*
 - Right-click the recording name and select *Rename Recording* from the list of options
3. The current name is highlighted. Copy & paste or type a new name in its place.

7.6 Export calibration results

For Screen and Scene camera projects, calibration results can be viewed and exported after data collection. This facilitates the evaluation of individual recordings for inclusion or exclusion in the data analysis and specification of calibration quality in reports and publications.

When viewing the recordings list in the Project Overview, you can:

- check whether a calibration has been performed by clicking on the arrow to the left of each recording. If you see a calibration icon, click it to display a results image
- review calibration data by hovering your cursor over the corresponding points in the image
- select one of the view options “Average Results” or “Left and right eyes results” in the bottom left-hand corner of the image
- save the calibration results image as a .png file.

Exporting the calibration data in spreadsheet format

1. Select a participant’s row in Project Overview.
2. Click the down arrow to see more information.
3. If there is a calibration result, you’ll see the estimated accuracy of the result.
4. The calibration results display with the total average results shown in a table at the bottom.
5. To save the image, click the Save image to file button at the bottom right and navigate to where you want to save the .png file. The export file contains the data of all recordings in the recordings list. The values in the export file are average values.



A screen project typically contains a calibration validation. This will be exported in the same way as a calibration without validation.

7.7 Change device settings

Pro Lab can interact with several devices during a recording, for example eye trackers, TTL parallel ports, external scene cameras, GSR units or monitors. Before the recording, you need to configure the devices that are represented as cards above the Timelines on the Record tab.

7.7.1 Add a device

Use the plus (+) button in the upper right on the Record tab. Here you can add a Shimmer GSR, microphone, and/or participant camera.

The color of the device card indicates its current state:

- *Gray* = active
- *Red* = inactive

7.7.2 Remove a device

For some recordings, you might not want to use all devices connected to your Pro Lab computer. For this reason, you can remove them and add them back when you need them.

How to remove a device:

1. Click to expand the card for the desired device.
2. Click the Remove button. The device is now inactive during recordings.

Add the device back again or for the first time using the steps in section [Add a device](#)

7.7.3 Change participant camera and microphone settings



Pro Lab does not support Bluetooth cameras and microphones. References in this manual to cameras and microphones assume a wired connection to the computer.

Participant camera

- Use the plus “+” button on the right-hand side of the Record tab to add a Shimmer GSR, microphone, and/or participant camera.
- If a participant camera is added, click the card to see the drop-down list of available cameras and select one. Click the drop-down list of video resolutions to select one. (See notes below)
- You can drag the participant camera view anywhere within the replay window using the Move button. You can expand the size of the participant camera view by clicking the expand button at the top right of the view window.
- The expander at the top right of the participant camera view window allows you to increase the window size.



Pro Lab displays a list of your camera’s available resolutions. Pro Lab selects the best possible frames per second (fps). Click the **Remove** button on the card to remove the camera.



Choose a lower resolution if you require a smaller file size.

Microphone

Add a microphone by clicking the plus “+” button above the Timelines on the Record tab. To use a microphone, click the card to see the drop-down list of available microphones and select one. Or click the Remove button on the card to remove it.

Replay a recording

In replay there are separate sliders for the sound in the video (stimulus) and the sound from the microphone. Click the sound icon to control recording sound (stimuli sound) or microphone sound. Drag the slider to zero if you don’t want to hear any sound at all.



If you record participant video or audio you can only use x1 (normal) speed to replay it in order for the synchronization between participant video/audio and Stimuli video to function correctly.

7.7.4 Change the eye tracker settings

The eye tracker card informs you about the eye tracker’s type, in what frequency it delivers gaze data, and in which mode it operates.

If the participant is positioned correctly in front of the eye tracker, their eyes are represented in the black square on the right side of the eye tracker's card at the top.

If you click on the eye tracker card, information appears about the eye tracker's serial number and firmware version. If there is a participant in front of the eye tracker, you will also see where in the track box she is positioned.

Change which eye tracker to use during recording

How to select the eye tracker:

1. Click the eye tracker card.
2. Click the drop-down menu on the card. The menu lists which eye trackers are available on the network at the moment (appear as light gray text) or that have been available previously to the Pro Lab computer (appear as dark gray text).
3. Select the desired eye tracker.



You can select an eye tracker which is currently not connected to the Pro Lab computer. If a recording is initiated using such a setting, a warning will be displayed to the moderator. The moderator can then select to abort the recording or to continue without recording gaze data.

Change settings for an eye tracker

Some eye tracker properties can be set and adjusted in the eye tracker card. When you have selected an eye tracker from the list of available ones, its setting options are displayed.



All settings are not available for all eye tracker models. Refer to your eye tracker's manual for more about available settings and how they affect gaze data. This is particularly important for the Eye Tracking Mode settings, as those can impact how gaze data is collected.

How to select gaze data frequency for a recording:

1. Click on the eye tracker card.
2. In the *Frequency* section of the card, select the desired frequency.

How to select the eye tracking mode for a recording:

1. Click the eye tracker card.
2. Open the drop-down menu in the *Eye Tracker Mode* section.
3. Select the desired eye tracker mode.

7.7.5 Change TTL parallel port settings

The TTL parallel port is used for signaling to external listeners about Pro Lab's estimates about how soon a stimulus appears on the presentation screen.



You make the TTL parallel port card visible by enabling the option "Send Stimulus onset markers (TTL)" in the Design module.

When you click on the TTL parallel port card, it provides information about the LPT port number (i.e. LPT1).

How to select the parallel port:

1. Click on the TTL parallel port card.
2. Open the drop-down menu.
3. Select the desired parallel port.

7.8 Stimuli presentation timing data

When analyzing gaze data, you need to know what the participant of the study saw during the data collection. For more information, go to the Tobii Pro Learning Center article [Stimulus presentation timing in Tobii Pro Lab](#).

7.8.1 Modify screen latency settings

Pro Lab timestamps the moment when the application assumes that a stimulus was displayed to the participant. This enables the application to associate the correct stimuli with the data. Its assumption is based on the screen latency, that is, the delay between sending a stimulus to the screen and the screen displaying it. However, a screen usually don't share this information with the connected computer, so if precise timing is required, you must enter this value manually. It is usually listed in the screen's technical specifications.



Pro Lab sends the TTL markers when it estimates that the stimulus has appeared on the screen. This estimation is based, unless otherwise specified, on a zero lag display and a 10ms screen latency. If you need to send TTL markers at stimulus onsets with higher timing precision, read the Tobii Pro Learning Center article, [Stimulus presentation timing in Tobii Pro Lab](#), about how to set up precise stimulus timing.

How to modify the estimated screen latency setting of the presentation screen:

1. Click on the screen in the *Target Screen* section of the desired *Record module* interface.
2. Click on the text ending with "ms*" under the representation of the presentation screen.
3. Enter the estimated latency of the desired screen.

8 Biometrics

8.1 Set up a Shimmer GSR sensor

A Shimmer GSR sensor measures the conductance in the participants' skin during the experiment. Tobii Pro Lab records these values over time.

Before using a Shimmer GSR sensor with Pro Lab. You need to pair the unit to your computer via Bluetooth. This is only necessary the first time you connect the sensor. After the initial pairing, the sensor appears in Pro Lab until you manually unpair it from your PC.

How to pair a Shimmer GSR sensor to your computer:

1. Open the Windows Bluetooth dialogue.
2. Switch on the Shimmer sensor. This enables other Bluetooth units to discover it.
3. When Windows detects the Shimmer sensor, it asks for a four-digit passcode. The default is "1234".
4. Enter the code. After a few seconds, the sensor's indication LED starts flashing slowly in blue to indicate the Bluetooth link is active. You can now use the Shimmer GSR sensor in Pro Lab.

How to add the Shimmer GSR sensor in Pro Lab:

1. Make sure that Shimmer's recording application – "Consensys Software" – is not running when you open Pro Lab.
2. Open the Record tab in Pro Lab.
3. Click on the plus (+) button in the upper-right corner.
4. Click on the Shimmer GSR Add button.
5. A new card for the Shimmer GSR sensor appears among the other recording device cards. It displays basic information about the sensor, such as its serial number, firmware version and frequency settings. The card also displays a GSR live-data representation, which is useful for checking the GSR data quality and see whether you need to adjust the participant's GSR sensor to improve the signal. Pro Lab can only record from one GSR sensor at a time. If you have more than one sensor paired, select the desired one.

8.2 Use the GSR Data Chart



This feature is only available for Screen Projects that use a Shimmer GSR+ sensor for recording GSR data.

The GSR Data Chart plots the galvanic skin response data recorded with the Shimmer GSR+ sensor. The GSR plot shows the participant's skin conductance over time, expressed in microsiemens (μS).

How to enable the GSR Data Chart:

1. Click on the downward pointing triangle next to “Visualizations” to the left of the Timeline. A menu appears.
2. Check the *GSR Data Chart* checkbox.
3. Close the menu by clicking anywhere outside of it.
4. Make all of the GSR Data Chart visible by clicking and dragging on the border between the video display area and the replay controls.

The chart shows the skin conductance value for each sampled GSR data point. If no time zoom is used (read [Replay settings](#)), the data gets down-sampled to fit the resolution of the Gaze Data Chart window and then not all individual data points will appear.

The “Use Autoscaling” toggle switch is located to the left of the GSR Data Chart. When it is off, the chart displays the signal’s complete range from zero to max. When it is on, the chart’s Y axis is autoscaled to the signal’s working range in the selected time window. The autoscaling function is limited to 0.5 μ S. If the GSR signal’s working range in the selected time window is below this value, the GSR signal appears centered on the Y axis with a range of 0.5 μ S.

When you let your mouse cursor hover over the GSR Data Chart, a cross appears with the mouse cursor at its center. You get information to the left of the cursor about where the horizontal line of the cross meets the GSR value axis. If the GSR data is not down-sampled, you will also get precise information to the left of the cursor about the GSR value and the time where the cross is placed.

8.3 GSR data filtering and analysis

Identifying skin conductance responses (SCRs) and determining their main characteristics is a common practice in galvanic skin response (GSR) research. SCRs can be produced in response to a specific event (e.g., stimulus onset), known as event-related skin conductance responses (ER-SCR), or appear spontaneously with varying rates.

In Pro Lab, the GSR analysis is made in the following way:

8.3.1 GSR data filtering

The GSR signal varies slowly over time. Rapid changes in the GSR signal are therefore considered to be external noise. When you analyze your GSR data, you need to remove two common types of noise or artifacts: high-frequency noise and rapid-transient artifacts. Pro Lab removes them by applying a median filter with a time window of 500 ms, followed by a mean filter with a time window of 1000 ms.

8.3.2 SCR detection analysis

After data filtering, Pro Lab applies a SCR detection algorithm to identify SCRs in the GSR data and calculate their main characteristics. Pro Lab does the following:

1. Pro Lab down-samples the GSR by an integer factor (only samples are deleted, not interpolations). Shimmer3 GSR + uses a down-sample factor of 8, resulting in a final sampling rate of approximately 15 Hz
2. Pro Lab applies a standard peak detection method (trough-to-peak) to identify local maxima and minima in the GSR signal. Trough-to-peak pairs get classified as SCRs when their

amplitude is higher than the minimum amplitude threshold of $0.03\mu\text{S}$.

3. Pro Lab calculates the main characteristics for each SCR detected:

- *SCR amplitude*: Amplitude difference between GSR level at SCR onset and GSR level at SCR peak.
- *SCR rise time*: Time difference between SCR peak time and SCR onset time.
- *SCR half recovery time*: Time difference between the time when the GSR level has recovered to 50% of the SCR amplitude and the SCR peak time. In two cases Pro Lab don't calculate this characteristic for a certain SCR(a): (1) when a second SCR (b) starts before the GSR level has recovered to 50% of the SCR(a), or (2) if the recording finishes before the GSR level has recovered to 50% of the SCR(a).

8.3.3 GSR metrics

SCRs can be generated as a response to a specific event (e.g., visual stimulus or unexpected question) known as event related SCR (ER-SCR). ER-SCRs are the most common measure used in research to relate changes in emotional arousal to a specific stimuli. A good stimulus design that allows enough time between stimuli is necessary to avoid uncertainties about which stimulus caused a certain ER-SCR. For more information, read [Understand the metrics](#)

A SCR or ER-SCR is included in an interval when the SCR onset time is within the interval, even when part of the SCR (e.g., SCR peak) is outside of the interval. The reason is that the GSR signal changes very slow, and the closest moment to the external event that triggers the SCR is the SCR onset time.

Here are two examples of when SCR/E-SCR events are valid/counted in a TOI interval:

In Example #1 (below), the SCRs and ER-SCRs with peaks outside a given interval *will be* included in the interval as long as their SCR onset time happened within the interval.




Figure 6. Example #1

In example #2 (below), the SCRs and ER-SCR with peaks within an interval *will not* be included in a given interval if their SCR onset time happened outside the interval.




Figure 7. Example #2

ER-SCR classification

A SCR is considered to be a response to a specific event (e.g., stimulus onset) when its onset falls within a specific time window after the event. In Pro Lab, all identified SCRs with their onset at 1s to 5s after a stimulus onset and custom time-of-interest (TOI) interval start get classified as ER-SCRs.

For each ER-SCR, Pro Lab calculates:

- **Event:** Name of the event that caused the ER-SCR.
- **Latency:** Time difference between the SCR onset time and the event time.

8.4 Export GSR Data

The data from the Shimmer GSR sensor is recorded in microsiemens (μS) and can be exported in sync with other sensor streams. For more information, read [Data export](#). The data can be exported in raw or filtered form. Raw data keeps the original sampling rate of 125Hz, whereas filtered data is down-sampled to 15Hz. When you export GSR data, you select between raw GSR data or filtered GSR data with the Tobii GSR filter on the Galvanic skin response (GSR) column.

9 Analyze data

The Analyzer module provides you with tools that analyze your recordings or export data for further analysis in applications such as MATLAB, Microsoft Excel, or SPSS.

9.1 Replay a recording

Open a replay tab by clicking on a recording in the Project Overview's Recordings list or by selecting a recording from the Analyze drop-down menu.

9.1.1 Replay controls

The replay controls are located under the panel showing the video (i.e. the video display area).

- Start video replay in the Replay tab by clicking the **Play** button at the bottom center.
- Pause the replay of a recording by clicking the **Pause** button at the bottom center. When a replay is in progress, the Play button is transformed into a Pause button.
- Step through a recording by clicking the arrows on either side of the **Play** button. Stepping through fixations or gaze points is determined by your fixation filter setting.

9.1.2 Replay settings

The replay settings are located under the panel showing the video (i.e. the video display area).

How to increase or decrease video replay speed:

1. Click on the *Replay speed control* icon on the bottom left. The *Replay speed* slider appears.
2. Increase (up) or decrease (down) replaying speed by clicking and dragging the slider's handle.

How to increase or decrease replay audio volume:

1. Click on the *Audio volume control* icon on the bottom right. The *Audio Volume* slider appears.
2. Increase (up) or decrease (down) volume by clicking and dragging the slider's handle.

How to change the Recording Timeline Zoom Range:

1. Zoom in to a specific location on the Replay Timeline by clicking and dragging the *Zoom* handles above it.
2. If you need a higher zoom level for the entire Timeline, use the *Time Zoom* slider above to the left of the Timeline.

How to hide the calibration in a recording:

1. Click the Calibration button (looks like a "Z" drawn between calibration points) above the Timeline to hide it.
2. Click the Calibration button again to show it.



The calibration will always be included in the video export since it is part of the recording.

9.1.3 Timeline tools

Timeline tools help you work with your recordings.

Timeline navigation

The Timeline contains a series of thumbnails that facilitates finding your desired recording section. In Glasses projects you can also display a second timeline preview below the first one.

In the upper timeline, open a zoomed-in visualization of the selected part by clicking the **Charts** list box and selecting the *Timeline preview* check box. The Timeline preview is affected by the applied *Time zoom* in the upper timeline.

Fast forward or through or jump to a location in the video replay:

1. Locate the red time needle handle (a red teardrop) on the Timeline below the video display area.
2. Click and drag the time needle to the desired starting point.

Move between event markers:

1. Use CTRL + arrow keys or your mouse to move the time needle to move the time needle between event markers.



When using the yellow selector tool, you can also press the CTRL button to snap to events.

Eye movements visualization

Timeline has a section where you can check whether gaze data is recorded for that time and if so, how the current Gaze Filter classifies it (i.e. fixation, saccade, or unknown). Fixations appear as solid bars, saccades as thin lines, and unclassified data as gray bars.

When you hover the mouse pointer over one of the three elements in this Visualization, a small overlay window appears with information about the eye movement data at that point in the recording:

- **Type:** The gaze point's kind of eye movement: "Fixation," "Saccade," or "Unknown". If a Raw Gaze Filter preset is used, all gaze points become "Fixation." If no gaze data is available for that point on the Timeline, the type is "No gaze data."
- **Gaze samples:** The number of gaze samples indicates how many samples are assigned to each eye movement classification (i.e. how many samples belong to the current "Fixation," "Saccade," or "Unknown" eye movements).
- **Start:** The time here indicates when in the recording, the current eye movement started. It uses the format HH:MM:SS:mmm.
- **Duration:** The time here indicates the duration of the current eye movement. It uses the format HH:MM:SS:mmm.

If you work in a Glasses project where gaze mapping onto a Snapshot has been performed, a gray background is displayed for the recording parts where gaze data has been mapped onto the Snapshot. For more information, read [Map eye tracking data](#).

Snapshot Timeline Visualizations

If you use Snapshots in your Glasses, scene camera, screen or web stimulus project, and you have enabled the **Show snapshot** switch in the Gaze data section of the Analyze tab, each Snapshot gets a row below the Visualization. The Snapshot's file name appears to the left of its row.

If no data is mapped onto the snapshot, the row appears empty except for the Snapshot's file name. If data has been manually mapped to the Snapshot, each fixation is represented by a square whose width corresponds to the fixation's duration on the Timeline. When a Raw Gaze Filter is selected, the mapped data is represented by squares whose width correspond to the time between two samples from the eye tracker (e.g. 20ms for data from Pro Glasses 2 with a 50Hz data rate).

When assisted mapping is used for a snapshot, a diagram is added to the snapshot Visualization for the time when assisted mapping was used. This diagram indicates the similarity between the frames in the recorded video and the snapshot and how well the Real-world mapping algorithm managed to map the gaze point on the snapshot. A high value indicates high similarity and a low value means there is a low similarity score. A low similarity score does not necessarily mean that the data is incorrectly mapped, only that the Real-world mapping algorithm had less information to base the mapping on and the mapping is therefore labeled as less similar.

All points that are mapped by the assisted mapping tool will receive a similarity score. This score can be used to speed up the inspection of the mapping result by highlighting points that received a low similarity score.

In the tools panel to the right you can set the Similarity threshold. All mapped points with a score below this threshold will be marked in orange to make it easier to identify them. The higher you set the threshold the more points, or fixations, will be orange.

You can inspect the points, or fixations, by placing the time needle at the beginning of the section to which you applied assisted mapping, and then step to the next mapped point, or fixation. You can change the behavior of the step to next fixation buttons to step between points, or fixations, that are below or above the threshold. If a mapped point, or fixation, is mapped incorrectly you can map it manually by clicking the snapshot. This overrides the existing mapped point. Manually mapped points receive a similarity score of 100 and are shown in green. If a point, or fixation, has a low similarity score but is correctly mapped on the snapshot, you can confirm it by clicking the confirm button or use the keyboard shortcut "C".

Gaze Data Chart

The Gaze Data Chart plots the angular eye velocity, gaze coordinates (raw or processed), and the classified fixation coordinates of a recording. The chart's main purpose is to help you set the parameters in the Gaze Filter, particularly when using an I-VT fixation filter, but it also provides you with direct access to the raw data when replaying a recording.

How to enable the Gaze Data Chart:

1. Expand the **Chart** menu to the left of the Timeline.
2. Check the *Gaze Data Chart* checkbox.

3. Close the menu by clicking anywhere outside it.
4. Make the whole Gaze Data Chart visible by clicking and dragging on the border between the video display area and the replay controls.

The data in the Gaze Data Chart can be toggled on or off using toggle switches. The axis to the left of the chart shows pixels, while the axis to the right is in degrees per second.

The available data in the chart is:

- **Gaze X:** This is the x-coordinate for each sampled gaze point. If no time zoom is used (read [Timeline tools](#)), the data gets down-sampled to fit the resolution of the Gaze Data Chart window (not all individual data points will appear).
- **Fixation X:** This is the x-coordinate of the fixations determined by the currently used Gaze Filter.
- **Gaze Y:** The y-coordinates for each sampled gaze point. If no time zoom is used (read [Timeline tools](#)), the data gets down-sampled to fit the resolution of the Gaze Data Chart window (not all individual data points will appear).
- **Fixation Y:** This is the y-coordinate of the fixations determined by the currently used Gaze Filter.
- **Gaze velocity w. threshold:** The thin line denotes the velocity calculated by the currently used Gaze Filter, while the thick line indicates the threshold used when determining whether the data gets classified as a fixation or not. (Data with velocities above the thick line is classified as saccades, and data with velocities below is classified as fixations.)
- **Cursor Quick Data:** When you hover your mouse cursor over the Gaze Data Chart, a cross with the mouse cursor at its center appears. If this setting is enabled, you get information to the right of the cursor about where the horizontal line of the cross meets the velocity axis. To the left of the cursor, you get an indication where it meets the pixel axis. Placing the cursor on a point in one of the graphs in the chart displays the value of that point in both pixels and velocity.

Mapped Gaze Data Chart



This feature is only available when the **Snapshot** button in **Replay Settings** is toggled *on*.

The Gaze Data Chart plots the angular eye velocity, gaze coordinates (raw or processed), and the classified fixation coordinates of a Snapshot. The data in the Mapped Gaze Data Chart is based on the fixations that are mapped onto it, either manually or by assisted mapping, for more information, read [Map eye tracking data](#). If the data was mapped using an I-VT Gaze Filter, all individual gaze points in the same fixation get the same x and y-coordinates, and there is no data between fixations. This arrangement always sets the gaze velocity below the velocity threshold.

How to enable the Mapped Gaze Data Chart:

1. Click the Charts drop-down menu to the left of the Timeline.
2. Check the **Mapped Gaze Data Chart** checkbox.
3. Close the menu by clicking anywhere outside it.

4. Make the whole Gaze Data Chart visible by clicking and dragging the border between the video display area and the replay controls.

The data in the Mapped Gaze Data Chart can be switched on or off by toggle switches. The axis to the left of the chart shows pixels while the axis to the right is in degrees per second.

The available data in the chart is:

- **Gaze X:** These are the x-coordinates for each sampled gaze point. If no time zoom is used (read [Timeline tools](#)), the data will be down-sampled to fit the resolution of the Mapped Gaze Data Chart window (not all individual data points will be shown).
- **Fixation X:** This is the x-coordinate of the fixations determined by the currently-used Gaze Filter.
- **Gaze Y:** These are the y-coordinates for each sampled gaze point. If no time zoom is used (read [Timeline tools](#)), the data will be down-sampled to fit the resolution of the Mapped Gaze Data Chart window (not all individual data points will be shown).
- **Fixation Y:** This is the y-coordinate of the fixations determined by the currently-used Gaze Filter.
- **Gaze velocity w. threshold:** The thin line denotes the velocity calculated by the currently-used Gaze Filter, while the thick line shows the threshold used when determining if the data should be classified as a fixation or not. (Data where the velocity is above the thick line is classified as saccades, and data with velocities below is classified as fixations.)
- **Cursor Quick Data:** When you hover your mouse cursor over the Mapped Gaze Data Chart, a cross with the mouse cursor at its center appears. If this setting is enabled, you get information to the right of the cursor about where the horizontal line of the cross meets the velocity axis. To the left of the cursor, you get an indication where it meets the pixel axis. Placing the cursor on a point in one of the graphs in the chart displays the value of that point in both pixels and velocity.

9.1.4 Replay tools

The replay tools are located in the tools panel to the right of the Analyzer module interface.

Recording information

Recording Information is associated with the recording that was generated at the time of its creation. This information cannot be modified.

The Recording Information includes:

- **Recording:** The name that was entered in the **Recording Name** input field in Pro Lab when the recording was created. In Glasses projects, the recording name is generated automatically when the recording is created.
- **Participant:** The name that was entered in the **Participant Name** field in Pro Lab when the recording was created. In Glasses projects, the participant name is entered when the recording is created.
- **Duration:** The duration of the recording displayed as HH:MM:SS,mmm (hours, minutes,

seconds, and milliseconds).

- **Date & Time:** The date and time when the recording was started.

Gaze Data Settings tools

During the replay of a recording, gaze data is first processed and classified by the Gaze Filter and then superimposed on the recorded video. The eye movement classification algorithm (the Gaze Filter) and the appearance of the visualized gaze data can be customized to fit the needs of the researcher. The Gaze Filter processes and classifies the recorded gaze data samples into fixations and other eye movements. The settings of the filter can be saved as presets. In Pro Lab, there are currently three built-in presets. For more details on how they work and what effect it has on your data, read [Gaze Filter functions and effects](#).

When the Tobii I-VT (fixation) filter preset is enabled, the classified eye movements (fixations, saccades, unknown eye movements) are visualized below the Timeline in the replay view. Fixations are visualized as light gray, thick lines, saccades are visualized as light gray, thin lines, and unknown eye movements are visualized as striped gray, thick lines. Due to the short duration of eye movements, it might be necessary to zoom into the Timeline in order to distinguish between individual eye movement instances

How to change the Gaze Filter preset:

1. In the Gaze Data section of the tools panel, select the **Settings** tab, and open the **Gaze Filter** drop-down menu.
2. In the drop-down menu, select the filter of your choice. The available options before any custom presets are created are:
 - **+ Create new Raw filter:** This option creates a custom filter and opens the filter preset configuration panel. See the steps below to create a new Raw filter preset.
 - **+ Create new I-VT filter:** This option creates a custom filter and opens the filter preset configuration panel. See the steps below to create a new I-VT filter preset.
 - **Raw:** If this preset is used, no classification of fixations, saccades, or unknown eye movements is done.
 - **Tobii I-VT (Attention):** This preset is optimized for wearable eye trackers and is used as the default preset for Glasses projects.
 - **Tobii I-VT (Fixation):** This preset is optimized for screen-based eye trackers and is used as the default preset for Screen projects. Read more about the preset settings in the white paper titled "*Determining the Tobii Pro I- VT Fixation Filter's Default Values*" on our website.
3. The new filter preset is now applied during replay.

How to create a Gaze Filter preset:

1. In the Gaze Data section of the tools panel, select the **Settings** tab and open the **Gaze Filter** drop-down menu.
2. In the drop-down menu, select the *+ Create new Raw filter* or *+ Create new I-VT filter* option, depending on which kind of filter preset you want to create.

3. Modify the settings of your choice in the filter. For details about what each setting is and what affect it might have on the gaze data, please refer to the white paper called “*Tobii Pro I-VT Fixation Filter*” on our website.
4. The filter is now automatically saved and set as the currently-used filter. It can be selected again from the *Gaze Filter* drop-down list by referring to the name given to it during creation.
5. Hide the filter preset configuration panel by clicking on the cogwheel next to the *Gaze Filter* drop-down menu.

How to change the Gaze Filter preset settings:

1. In the *Gaze Data* section of the tools panel, select the **Settings** tab and open the **Gaze Filter** drop-down menu.
2. In the drop-down menu, select the filter you want to modify.
3. Click on the cogwheel next to the *Gaze Filter* drop-down. This opens the filter preset configuration panel. Built-in presets in Pro Lab cannot be modified.
4. Modify the settings of your choice in the filter. For details about what each setting and what affect it might have on the gaze data, please refer to the white paper called “*Tobii Pro I-VT Fixation Filter*” on the [Tobii Pro](#) website.
5. The modifications to the filter are automatically saved and applied to the data in the currently-replayed recording. These settings are per recording and not global settings.
6. Hide the filter preset configuration panel by clicking on the cogwheel next to the *Gaze Filter* drop-down menu.

Events tools

Tools for logging interesting and important Events in the recordings are included in the software. These Event Types (e.g. “Participant picks up item 1”) can also be exported together with all the other data collected during a recording. The logged Events, as well as automatically-generated Events, appear on the Timeline and in the Events List in the Events section of the tools panel on the right.

To manually log Events, you must first define the Event Types scheme that contains the Events you want to log and define which keyboard keys you want to use for logging the Events.

It is not possible to create two Event types with the same name within a project. If Event Types are created with the same name but on different computers and coding data from these projects are later merged into the same project, the Event types will also be merged.

How to define an Event Types logging scheme:

1. Locate the **Events** section on the Analyze tab in the tools panel on the right.
2. Click the **New Event Type** button at the top right of the panel. The New Event Type dialogue displays.
3. Enter the name/description of the Event type in the left-hand input field.
4. Click the drop-down menu and select the keyboard shortcut that you would like to use for

logging this Event.

5. Click **OK**.



An Event Types logging scheme can also be created in the Project Overview section of the application. For more information, read [Manage Event Types](#).

How to edit an existing Event Type:

1. Locate the **Events** section on the Analyze tab in the tools panel on the right.
2. Place the mouse cursor over the custom Event you want to edit. Three icons appear to the right.
3. Click the *Pen* edit icon. An *Edit Event Type* dialogue appears.
4. Edit the name/description of the Event in the input field at the left of the dialogue.
5. Click the drop-down menu and select the keyboard shortcut that you would like to use for logging this Event.
6. Click **OK**.



Event Types can also be modified in the Project Overview section of the application. For more information, read [Manage Event Types](#).

How to delete an Event Type from the Event Types scheme:

1. Locate the **Events** section on the Analyze tab in the tools panel on the right.
2. Place the mouse cursor over the custom Event you want to edit. Three icons appear to the right.
3. Click the *Delete* icon that looks like a trash can. A Delete Event Type dialogue opens.
4. Click Delete Event to confirm.



Event Types can also be deleted in the Project Overview section of the application.

How to log an Event manually:

1. Make sure you have defined the Event types that you would like to log.
2. While replaying or pausing a recording, press the shortcut key on your keyboard associated to the Event you want to log. Alternatively, click the Log icon that appears when placing the mouse cursor over an *Event Type* in the *Event Types* list in the *Events* section of the tools panel on the right.

You can filter which kinds of Events are visible in the Events list. The different kinds of Events are grouped into Event groups. The filter also affects which Events are visible on the Timeline.

Available Event groups include:

- **Recording Events:** Events included in this group are general Events relating to the recording (e.g. *Recording Start*, *Recording End*, *Eye Tracker Calibration start*, and *Eye Tracker Calibration end*).
- **Custom Events:** Events included in this group are Events generated from the Event Types created by the user.
- **Sync Events:** Events included in this group are Events generated in order to synchronize the gaze data with data from another data source (e.g. *Sync Port Out High*, *TTL in*, and *TTL out*).
- **Screen project-specific Event groups:**
 - **Participant Events:** Events included in this group are Events generated by the participant (e.g. *Mouse Events* and *Keyboard Events*).
Mouse events can have four possible values from the combinations of Up/Down, Left/Right.
Down means that the mouse button is pressed. *Up* means that the mouse button is released.
Left and *Right* refer to the left and the right mouse buttons). Keyboard events have values of all keys on the keyboard timestamped with the press moment.
 - **Stimuli Events:** Events included in this group are Events relating to the stimulus presentation (e.g. *Image Stimulus Start*, *Image Stimulus End*, *Video Stimulus Start*, and *Video Stimulus End*).
- **Glasses project-specific Event groups:**
 - **Logged live Events:** Events included in this group are user-created Events generated during recording.
 - **Snapshot mapping Events:** Events included in this group are generated when mapping gaze from a recording onto a Snapshot (e.g. *Interval Start* and *Interval End*).

How to use Events to navigate on the Timeline:

1. Locate the Event you want to navigate to in the Events list.
2. Click on the Event. The track slider on the Timeline will now have moved to the time of the Event on the Timeline.

Select a frame and pair it with Times of Interest

When analyzing dynamic content, you may select a background image and assign it to a section of a recording. The same background image can also be assigned to multiple sections of multiple recordings. This workflow is used to produce visualizations and extract AOI statistics from media elements that display dynamic content. This includes screen-based, scene camera, and external presenter recordings.

Example: when analyzing interaction with a phone using the Mobile Device Stand, you may select the starting screen of the phone as the background image and assign the times in each recording when this screen was shown.

How to select and assign a background image:

1. Find a frame in the video that you would like to use as the background image.
2. Click the plus sign list box just above the timeline to the right and select “Save frame as media.”
3. Name the file in the pop-up window that appears and select **OK**.
4. You can assign a time interval to the background image by defining a Custom Time of Interest or by using an auto generated Times of Interest (TOI). The background image can be selected for the TOI when you are editing the TOI.

You are also able to create a TOI and assign a background image (frame) in one step. If you select the plus sign list box just below the window to the right and select “Create TOI with frame.”



Selecting a frame and pairing it with a Time of Interest is only available in Screen and Scene Camera projects.

9.2 Multi-select recordings

You can multi-select recordings on the Project Overview tab using CTRL + click.

Right-click and select *Open in background tab* to open multiple recordings on the Analyze tab simultaneously.

or

Right-click and select *Delete recordings* to delete multiple recordings.

9.3 Export recordings


You can export a video clip or an entire audio and/or video recording.

9.3.1 Export a video clip

You can export video clips with a gaze overlay from Pro Lab by using the recording’s Timeline within a Replay tab.

Exported video clips can also contain participant audio and video.

How to export a video clip with a gaze overlay:

1. Select the interval on the Timeline you want to export by dragging the yellow handles on either side of the red track slider to the desired start and the end.
2. Click the **Export** button  above the timeline of the recording.
3. Select the **Recording media** you want to export.
4. If you export **Sound** from a stimuli *and* participant recording, you need to select which source to use. You can only choose one source which means you can not export both the stimuli and the participant sound recordings.
5. Select **Video resolution** for export. 1920 x 1080 is the default.
6. If Participant video is selected, you can specify the size of it (in % of total screen). You can also choose to display video using *Picture in picture* or *Side by side* mode.

- a. **Picture in picture** enables you to click and drag the participant recording around in the preview window.
 - b. **Side by side** enables you to show participant video on the left or right side using the arrows.
7. Save the .mp4 file to your computer. The default name is “Recording name Participant name” but you can change this to whatever you want.
8. Click the **Export** button to start the export. You can view exporting progress in the activity queue in the top right of Pro Lab. Clicking the activity queue and then clicking a specific export allows you to cancel the export, open the file location when it’s finished exporting, or delete it from the queue.



Figure 8. The Activity queue button lets you see and control the exporting process

9.3.2 Export a whole recording

You can export participant or stimulus audio and video as separate or merged files.


1. Select a recording from the list on the Analyze tab or from the left-hand side menu on the Project Overview tab.
2. Click the **Export** button  above the timeline of the recording.
3. The Export Recording dialog opens.
4. Select the **Recording media** you want to export.
5. If you export **Sound** from a stimuli *and* participant recording, you need to select which source to use. You can only choose one source which means you can not export both the stimuli and the participant sound recordings.
6. Select **Video resolution** for export. 1920 x 1080 is the default.
7. If Participant video is selected, you can specify the size of it (in % of total screen). You can also choose to display video using *Picture in picture* or *Side by side* mode.
 - a. **Picture in picture** enables you to click and drag the participant recording around in the preview window.
 - b. **Side by side** enables you to show participant video on the left or right side using the arrows.
8. Save the .mp4 file to your computer. The default name is “Recording name Participant name” but you can change this to whatever you want.
9. Click the **Export** button to start the export. In the activity queue you can cancel the export, open the file location when its finished exporting, and delete it from the queue.



Figure 9. The Activity queue button lets you see and control the exporting process

9.4 Map eye tracking data

Wearable eye tracking devices, such as Pro Glasses 2 and 3, produce eye gaze data mapped to a coordinate system relative to the wearable eye tracker and the recorded video, not to static objects of interest in the environment.

For most kinds of statistical/numerical analysis to be meaningful, the collected eye tracking data needs to be mapped on to objects, or images, of interest and into a new coordinate system in order to be aggregated.

Remote eye tracking devices produce eye gaze data that is mapped to a coordinate system corresponding to the image/ video stimulus, Screenshot, or frame from the recording. When the screenshots or frames are not corresponding to what the participant has seen the aggregating of data is not possible.

Pro Lab addresses these challenges with a solution called snapshot mapping.

9.4.1 Assisted and manual mapping

The gaze data mapping in Pro Lab lets you map eye gaze data onto still images (snapshots) of environments and objects of interest or external screenshots in two different ways: either semi-automatically using manual mapping or fully automatic using the assisted mapping function. The Snapshot images are typically photos, created by the data collector or the researcher using a screenshot tool or text editor. The assisted mapping function uses Pro Lab's Real World Mapping software technology to automatically map gaze data from a recording onto a Snapshot.

Requirements for Snapshot images:

- Less than 25 megapixels
- PNG, JPG, GIF, or BMP file format

Requirements for Assisted mapping in screen and web recordings:

- Full HD (1920 x 1080) or 1920 x 1200

Requirements for Assisted mapping in scene camera recordings:

- If the camera uses full HD, then (1920 x 1080) or 1920 x 1200



Tobii Pro does not recommend using stitched images with assisted mapping as it may degrade performance. (Stitched images are images that are put together from a series of images covering a greater area than the camera was capable of covering in just one shot.)

If Assisted Mapping is not working as expected, read our online [troubleshooting article](#) on Tobii Pro Connect.

How to import a Snapshot from the Tools panel in the Replay tab:

1. In the Gaze Data section of the Toolspanel, select the **Snapshots** tab.
2. To the right of the text "Snapshot images," click the "+" icon. The file browser opens.
3. Locate the image file (*.bmp, *.gif, *.png, or *.jpg) on your computer that you want to use as a Snapshot.

4. Click **Open**. An import progress bar displays.
5. When the Snapshot image has been imported, click **OK**.

How to import multiple instances of the same Snapshot image:

1. Create and save multiple copies of the image file you want to use as Snapshots.
2. Give each image file a unique file name.
3. For each copy of the image file, follow the instructions in the section above to import the Snapshot images.

To make the manual mapping process efficient, Pro Lab lets you map entire fixations onto a Snapshot with just one click, rather than mapping each gaze point individually. For more information on how this affects the data, read [Effect of mapping data onto Snapshots with Tobii I-VT or Raw Gaze Filter](#). If you have chosen to map data with the Raw Data Gaze Filter preset selected, the data will be mapped gaze point by gaze point. Gaze data from a recording can be mapped on one or several snapshots.

How to map data onto a Snapshot manually:

1. In the Gaze data section of the **Tools** panel, select the **Snapshots** tab.
2. Enable Mapping by toggling on **Show snapshot**.
3. Enable or disable the "Automatically step on next fixation" toggle switch. Enabling this will cause the paused replay to automatically jump to the next fixation/raw data point on the Timeline when a gaze point has been manually mapped. This eliminates the need to use arrow keys to step forward manually on the timeline.
4. In the grid/list of Snapshot images, select the Snapshot onto which you want to map data. You can also select which Snapshot to map data onto from the list of Snapshots located below the replay Timeline. There, each Snapshot is represented by a thumbnail as well as a row on which it will be displayed for which parts of the recording data has been mapped. At any time during the mapping of data, you can switch back and forth between different Snapshots without losing mapped data.
5. While skimming through the recording replay, locate and pause the video at the start of the section that you want to map onto the selected Snapshot.
6. To map data onto the Snapshot, first locate the gaze data point (circle superimposed on the video) in the recorded video. Click once in the corresponding location on the Snapshot image as precisely as possible.
7. Continue this process until all data has been mapped onto the active Snapshot. As data points are mapped onto the Snapshot, the Snapshot timeline will indicate at which times data points have been mapped.
8. Replay or manually step through the recording using the arrow keys once the mapping is completed and compare the mapping on the Snapshot with the gaze locations in the video to verify that data has been mapped correctly. If you want to move a mapped point, right click on it and select "Delete current manually mapped fixation point" from the menu. Then click on the Snapshot to map the gaze point in a new location.

How to map data onto a Snapshot automatically with assisted mapping:

1. In the Gaze Data section of the **Tools** panel, select the Snapshots tab.
2. Enable Mapping by toggling on **Show snapshot**.
3. In the grid/list of Snapshot images, select the Snapshot onto which you want to map data. You can also select which Snapshot to map data onto from the list of Snapshots located below the replay Timeline. There, each Snapshot is represented by a thumbnail as well as a row on which it will be displayed for which parts of the recording data has been mapped. At any time during mapping of data, you can switch back and forth between different Snapshots without losing mapped data.
4. Select the interval on the Timeline you want the gaze points to be mapped automatically in. To select an interval, drag the yellow handles on either side of the red track slider to where you want the start and the end of the interval to be. If needed, you can zoom in on the Timeline to make the interval selection easier. This is most often the part of the recording where the location or object shown on the Snapshot comes into view.
5. Right-click on the selected interval or click the ellipsis (...) located directly over the Timeline, and select "Run assisted mapping." The interval is now placed in the processing queue. The built in Real-world mapping algorithm starts processing the mapping automatically according to the order in the processing queue. If another mapping is already in progress, that mapping will be completed before the next one is initiated. You can check the jobs placed in queue by clicking the number at the top right of the window.
6. You can choose to create another mapping task by repeating steps 4 to 6 and place it in the processing queue, or, if you don't have any more pending tasks, continue to step 8.
7. When the assisted mapping is completed, a diagram is added to the section of the recording for which the mapping has been done on the row representing the Snapshot under the Timeline. The diagram indicates how confident the Real-world mapping algorithm is about the similarity of the gaze point in the recording and the mapped position in the snapshot. A high value indicates high similarity, and a low value, a low similarity level. A low similarity level does not necessarily mean that the data is incorrectly mapped, just that the Real-world mapping algorithm had less information on which to base the mapping. Therefore, it is labeled as less similar.
8. Sections with low similarity should be reviewed and, if necessary, remapped manually (if incorrect mappings are found). Adjust the similarity threshold, in the tool panel on the right, to a level that fits the requirements of your project. Sections above the threshold will be marked in green and sections below the threshold will be marked in orange for easier identification. If you are using a fixation filter, the highlight will be based on the mean similarity of all the mapped gaze points of each fixation. The individual similarity of each gaze point in a fixation can be seen as a green or orange dots.
9. Remap gaze points manually using the following steps:
 - a. For each point or fixation, you can do one of three things:
 - remove it by pressing the **Delete** button on the keyboard or by right-clicking the point on the Snapshot and selecting *Delete current automatically mapped fixation*

point in the dialog

- correct it by manually clicking and moving the gaze point to where it should be
- leave it “as is” by clicking the **Accept** button or clicking the “C” key on your keyboard



You can adjust the behavior of the **Step to next** button to step between mapped gaze points, or fixations, below the similarity threshold or step between all. Buttons for these settings are located to the right of the navigation buttons.



(A) displays when the point is automatically mapped and below the threshold. (M) denotes a new manually mapped point.



Figure 10. An automatically mapped point (A) and a manually mapped point (M)

b. Automatic color-coding:

- In the row representing the Snapshot under the Timeline, manually remapped gaze points and fixations appear in solid green without a graph.
- In the Snapshot image, automatically-generated mappings appear in a green circle, whereas manually mapped points appear in a red circle.
- Overridden (i.e. deleted) points mapped by assisted mapping appear in a gray circle.

c. Repeat this procedure until you are done.

10. Replay or manually step through the recording using the arrow keys once the mapping is completed and compare the mapping on the Snapshot with the gaze locations in the video to verify that data has been mapped correctly.




Figure 11. Step through the recording using the arrow keys to compare the mapping on the Snapshot with the gaze locations

9.5 Visualize eye tracking data

In the Visualizations tool in Pro Lab, data can be visualized using Heat Maps and Gaze Plots.

9.5.1 Create and customize a Heat map

Heat Maps can be of great value when creating reports, papers, or presentations as they help you summarize large quantities of data in an intuitive way. Heat Maps are created on top of stimuli such as Snapshots or images. A Heat Map uses different colors to illustrate the number of fixations participants made within certain areas of the stimulus or for how long they fixated within that area. Red usually indicates the highest number of fixations or the longest time fixating there, and green the least, with varying levels in between. For more detailed information about Heat Maps, read [Calculate heat maps](#).



To be able to create Heat Maps on Snapshots in Glasses projects, at least one Snapshot in the project must have gaze data mapped onto it.

9.5.2 Create a Heat Map

How to create a Heat Map:

1. Click on the *Visualizations* button - either in the *Project Overview* section or select it from the *Analyze* drop-down menu in the top navigation.
2. Go to the *Visualization Type and Settings* tool in the tool panel, to the right of the Visualization interface and select *Heat Map*.
3. Select what Gaze Filter to use when generating your Heat Map. For more information about Gaze Filters, read [Gaze Filter functions and effects](#).
4. Go to the *Times of Interest* tool and select the Time of Interest (TOI) that you want to use as the base of your Visualization. TOIs based on start and end Events for stimuli, such as images or Snapshots, are represented by a thumbnail of that stimulus.
5. Go to the Recordings tool and select or deselect the recordings and intervals that you want to visualize in the Heat Map. Note that all TOIs generated per recording are represented by one interval in that recording. If the same stimulus was displayed multiple times, each time

generates a TOI which then gets represented by an interval in that recording. This means that it is possible to analyze, for example, all data collected when a participant saw a stimulus during a recording or only parts of that data such as the first-time exposure. For TOIs relating to Snapshots in Glasses projects, the segments are generated by Snapshot start and Snapshot end Events. If there is a gap in the data mapped onto a Snapshot of five seconds or more, a Snapshot end Event is generated automatically at the time of the last mapped gaze point and a Snapshot start Event when the next gaze point is mapped after the gap. To see the intervals included, click the arrows to the left of the recording checkboxes to expand the view.

6. Go to the Visualization Type and Settings tool in the tool panel at the right of the Visualization interface, open the Type drop-down menu, and select the calculation basis for the Heat Map. For more information, read [Calculate heat maps](#).

The available options are:

- *Absolute count* – calculated by the number of fixations (or gaze data samples if using a Raw data Gaze Filter).
- *Absolute duration* – calculated by the duration of fixations (or gaze data samples if using a Raw data Gaze Filter).
- *Relative count* – calculated by the number of fixations relative to the total number of fixations made by the participants during the Time of Interest (or gaze data samples instead of fixations if using a Raw data Gaze Filter).
- *Relative duration* – calculated by the duration of the fixations relative to the sum of all fixation durations mapped in the Time of Interest on the Snapshot (or gaze data samples instead of fixations if using a Raw data Gaze Filter)

7. The Heat Map is now generated.



If the selected TOI is not associated with a Snapshot or image stimulus, the background for the generated Heat Map is set to Transparent by default.

9.5.3 Customize the appearance of a Heat Map

The appearance of a Heat Map can be modified to better convey research findings. No underlying data is affected by changing any of the properties listed below.

The following properties of a Heat Map can be modified to change its appearance:

- **Background:** This setting modifies what the Heat Map is generated onto. It provides the following options:
 - **Image:** If this option is selected, the Heat Map will be rendered on top of the stimulus image or Snapshot associated with the selected TOI.
 - **Transparent:** Selecting this option results in a Heat Map generated on a transparent background. This is useful if exporting the Heat Map as an image for further work in (e.g. an image processing tool where it will be combined with some other background image).
 - **Black:** This option produces a Heat Map rendered on a black background.
 - **White:** This option produces a Heat Map rendered on a white background.

- **Radius:** This setting adjusts the size of each fixation’s contribution to the Heat Map. For more information, read [Calculate heat maps](#). The property can be changed by using the slider dragging it to the left will make the area smaller and, to the right, will make it larger. You can also click on the number to the right of the slider. Clicking on the number transforms it into an input field where you can type the desired radius of the fixation contribution.
- **Color:** This setting governs the appearance of the scale used for the Heat Map. To the left is what a low level of fixations or time fixating will look like, while the right side represents what an area with a lot of fixations or time fixating will appear as. The setting allows you to specify the color for the low, mid, and high level on the scale individually by offering a drop-down menu with a color picker for each level. Click on the color you want to use or enter a hexadecimal color code in the number entry field.
- **Opacity:** This setting determines the transparency of the Heat Map data.

Scale max value: This setting adjusts the value for which the max level (color to the right) in the color scale specified by the Color setting is reached. By default, the checkbox for this setting is unchecked, which means that the value is set automatically by Pro Lab. This setting is unavailable for Heat Maps for Relative count or Relative duration where the max value is always set to 100%.

9.5.4 Save a Heat Map as an image

Heat Maps can be saved as image files to be incorporated later into reports or further modified in image-processing software. Available image file formats are *.png and *.jpg.

How to save a Heat Map:

1. Right-click within the Heat Map area. A Context menu appears,
2. Select *Save to File...* A file browser opens.
3. Locate the folder you want to save the file in.
4. In the *File name* field, enter a file name for your Heat Map.
5. Open the *Save as Type* drop-down menu and select the desired file type (*.png or *.jpg). (For Heat Maps with transparent backgrounds, the file type must be *.png to maintain the transparency.)
6. Click *Save*.

9.5.5 Create and customize a Gaze Plot

The Gaze Plot Visualization shows the sequence and position of fixations (dots) on a stimulus, such as an image or Snapshot. The size of the dots indicates the fixation duration and the numbers inside the dots represent the order of the fixations. Gaze Plots can be used to illustrate the gaze pattern of a single test participant throughout an entire recording or of several participants in a short time interval (e.g. when a specific stimulus was shown).



When you create Gaze Plots on Snapshots in Glasses projects, at least one Snapshot in the project must have gaze data mapped onto it.

9.5.6 Create a Gaze Plot

How to create a Gaze Plot:

1. Click the **Visualizations** button, either in **Project Overview** section or select it from the **Analyze** drop-down menu in the top navigation.
2. Use the **Visualization Type and Settings** tool in the tool panel to the right of the Visualization interface and select *Gaze Plot*.
3. Select the desired *Gaze Filter* for generating your Gaze Plot. For more information about Gaze Filters, read [Gaze Filter functions and effects](#).
4. Use the **Times of Interest** tool and select the Time of Interest (TOI) that you want to use as the base for your Visualization. TOIs based on start and end Events for stimuli such as images or Snapshots are represented by a thumbnail of that stimulus.
5. Use the **Recordings** tool to select or deselect the recordings and/or intervals that you want to visualize in the Gaze Plot. All TOIs generated for a recording are represented by one interval in that recording. If the same stimulus was displayed multiple times, each time generates a TOI which then is represented by an interval in that recording. This means that it is possible to analyze, for example, all data collected when a participant saw a stimulus during a recording or only parts of that data, such as the first-time exposure. For TOIs relating to Snapshots in Glasses projects, the segments are generated by the Snapshot start and Snapshot end Events. If there is a gap in the data mapped onto a Snapshot of five seconds or more, a Snapshot end Event is generated automatically at the time of the last mapped gaze point and a Snapshot start Event is made when the next gaze point is mapped after the gap. To see the intervals included, click the arrows to the left of the recording checkboxes to expand the view.
6. The Gaze Plot is now generated.



If the selected TOI is not associated with an image or Snapshot stimulus, the background for the generated Gaze Plot is Transparent by default.

9.5.7 Customize the appearance of a Gaze Plot

The appearance of a Gaze Plot can be modified to better convey research findings. No underlying data is affected by changing any of the properties listed below.

The following properties of a Gaze Plot can be modified to change its appearance:

- **Background:** This setting modifies what the Gaze Plot is generated onto. It provides the following options:
 - **Image:** If this option is selected, the Gaze Plot will be rendered on top of the stimulus image or Snapshot associated with the selected TOI.
 - **Transparent:** Selecting this option results in a Gaze Plot generated on a transparent background. This is useful if you are exporting the Gaze Plot as an image to work further with (e.g. an image-processing tool where it will be combined with some other background image).
 - **Black:** This option produces a Gaze Plot rendered on a black background.
 - **White:** This option produces a Gaze Plot rendered on a white background.
- **Fixation:** This setting specifies what determines the size of each fixation point in the Gaze Plot. The options are the following:

- **Duration:** The size of the point is determined by the duration of the fixation it represents. A long fixation is represented by a large point and a short fixation by a small point.
- **Same size:** If this option is selected, all points will get the same size independent of the duration of the underlying fixation.
- **Scale:** This property can be changed by using the slider dragging it to the left will make the point smaller and moving it to the right makes it larger (or by clicking on the number to the right of the slider). Clicking on the number transforms it into an input field where you can type the desired radius of the fixation point.
- **Border color:** This setting governs the color of the border surrounding the gaze points. There are two options available:
 - **Custom color:** If this option is selected, you can select which border color to use by using a drop-down menu with a color picker which appears to the right of the setting. To select a custom color, click on the drop-down menu to expand it and then on the color you want to use. You can also enter a hexadecimal color code in the number entry field. The same color will be used for all gaze points in the Gaze Plot.
 - **None:** Selecting this option results in the gaze points being shown without any border at all.
- **Fill color:** This setting governs the fill color of the gaze points. There are two options available:
 - **Recording colors:** If this option is selected, the color associated with each recording, as listed in the Data Selection tool, will be used for gaze points representing fixations from that recording.
 - **Same color:** If this option is selected, all gaze points will have the same fill color. To select a fill color, click on the drop-down menu to the right of the Fill color options drop-down menu to expand it. Then, click on the color you want to use. You can also enter a hexadecimal color code in the number entry field.
- **Show gaze order:** This setting either shows or hides the order in which the fixations were made, represented by numbers in the gaze points in the Gaze Plot.
- **Show gaze trail:** This setting either shows or hides the lines connecting the gaze points, which represent fixations that happen sequentially.
- **Opacity:** This setting determines the transparency of the gaze points in the Gaze Plot. This setting only changes the opacity of the gaze point fill. The gaze point border opacity will always be set to 100% unless the “Border color” setting is set to “None.”

9.5.8 Save a Gaze Plot as an image

Gaze Plots can be saved as image files that can be incorporated later into reports or further modified in image-processing software. Available image file formats are *.png and *.jpg.

How to save a Gaze Plot:

1. Right-click within the Gaze Plot area to display the menu.
2. Click **Save to File...** A file browser opens.
3. Locate the folder you want to save the file in.

4. In the **File name** field, enter a file name for your Gaze Plot.
5. Open the **Save as Type** drop-down menu and select the desired file type (*.png or *.jpg). (For Gaze Plots with a transparent background, the file type must be *.png to maintain the transparency.)
6. Click **Save**.

9.5.9 Modify the Gaze Filter setting

The eye movement classification algorithm (the Gaze Filter) and the appearance of the visualized gaze data can be customized to fit the needs of the researcher. The Gaze Filter processes and classifies the recorded gaze data samples into fixations and other eye movements. The settings of the filter can be saved as presets. In Pro Lab, there are currently three built-in presets. For more information on how they work and what effect it has on your data, read [Gaze Filter functions and effects](#).

How to change the Gaze Filter preset:

1. Go to the *Visualization Type and Settings* section of the tools panel and open the *Gaze Filter* drop-down menu.
2. Select the desired filter. The available options (before any custom presets are created) are:
 - **+ Create new Raw filter**: This option creates a custom filter and opens the filter preset configuration panel. See the instructions below on how to create a new Raw filter preset.
 - **+ Create new I-VT filter**: This option creates a custom filter and opens the filter preset configuration panel. See the below instructions on how to create a new I-VT filter preset.
 - **Raw**: If this preset is used, no classification into fixations, saccades, or Unknown Eye Movements is done.
 - **Tobii I-VT (Attention)**: This preset is optimized for wearable eye trackers and is used as the default preset for Glasses projects.
 - **Tobii I-VT**: This preset is optimized for screen-based eye trackers and is used as the default preset for Screen projects. Read more about the preset settings in the white paper called "*Determining the Tobii Pro IVT Fixation Filter's Default Values*" found on our website.
3. The new filter preset is now applied to the data used in the Visualization.

How to create a Gaze Filter preset:

- a. Go to the *Visualization Type and Settings* section of the tools panel and open the *Gaze Filter* drop-down menu.
- b. Select the *+ Create new Raw filter* or *+ Create new I-VT filter* option, depending on which kind of filter preset you want to create.
- c. Modify the desired settings in the filter. For details about what each setting means and what effect it might have on the gaze data, read the white paper "*Tobii Pro I-VT Fixation Filter*" which is available on the Tobii Pro website.

- d. The filter is now automatically saved and set as the currently-used filter. It can be selected again from the *Gaze Filter* drop-down list by referring to the name given to it during creation.
- e. Hide the filter preset configuration panel by clicking on the cogwheel next to the *GazeFilter* drop-down menu.

How to change the Gaze Filter preset settings:

- a. Go to the *Visualization Type and Settings* section of the tools panel and open the *Gaze Filter* drop-down menu.
- b. Select the filter you want to modify.
- c. Click on the cog wheel next to the *Gaze Filter* drop-down menu. The *Filter Preset Configuration* panel appears. Built-in presets in Pro Lab cannot be modified.
- d. Modify the desired settings within the filter. For details about what each setting means and what effect it might have on the gaze data, read the white paper called “*Tobii Pro I-VT Fixation Filter*” which is available on the Tobii Pro website.
- e. The modifications to the filter are automatically saved and applied to the data in the Visualization.
- f. Hide the filter preset configuration panel by clicking on the cogwheel next to the *Gaze Filter* drop-down menu.

9.5.10 Change the Gaze Trail settings

The Gaze Data Settings tab includes the following settings for the Gaze Trail feature:

- *Fill Color* – This is the fill color for the gaze point.
- *Contrast Color* – This is the color of the gaze point’s border.
- *Opacity* – This value decides the transparency of the gaze point. A value of 0 means that the gaze point is completely transparent, just like if no fill color was selected. A value of 100 means that the gaze point is completely solid, filled with the color selected in *Fill Color*.
- *Size* – This is the initial radius size of the gaze point, where 100% represents a tenth of the replay window’s height.
- *Fading duration* – This value decides how long a trail point is shown before it fades out. A longer value creates a longer gaze trail.
- *Maximal size* – This value decides the maximum size the gaze point reaches (where 100% represents a tenth of the replay window’s height, as for *Size* above). If this value is smaller than the *Size* value, no animation will be displayed.



For glasses recordings, please note that the gaze trail only consists of lines. That is, no gaze points will be shown.

9.6 Metrics visualization

Metrics Visualizations allow you to explore a plot of data to quickly discover how the metrics data are distributed, and whether there could be a difference between the study conditions. This tool is


intended to be used for review of a small set of pilot data or to share some initial findings with colleagues before you proceed to collect all your data.

Metrics Visualizations are designed to visualize data from multiple trials or intervals. A set of intervals is typically defined by a Custom Time of Interest, so a Custom Time of Interest is required to generate a metrics visualization. For more information, read [Custom Times of Interest](#).

Metrics Visualizations are designed to compare data that occur in at least two different conditions. There are three types of experimental variables in Pro Lab that can be used to create conditions. One of the following experimental variables is required to generate a metrics visualization: Create Tag Groups, Stimulus Variable, or Participant Variable. To find out more information about these, read [Create Tag Groups](#), [Stimulus/group variables](#), and [Participant Variables](#).

9.6.1 Create a Metrics Visualization

1. In the **Analyze** module, select **Metrics Visualizations** by clicking the arrow next to "Analyze" at the top.
2. The **Settings** panel on the right includes settings for the Times of Interest and Data Selection.
 - a. In the **Times of Interest** section, select the Custom Time of Interest you want to visualize. You can also create new Custom Times of Interest using the plus button (+) and edit an existing one by right-clicking it.
 - b. In **Data selection** all Recordings, Participant Variables, Areas of Interest, Aggregated AOs, and/or Events are selected by default. Deselect any data to exclude it from the visualization.
3. The metrics available for Metrics Visualizations are listed on the left-hand side.
 - Select the metrics you want to visualize.
 - Binned metrics are labeled "Binned." For more information, read [Binned metrics](#).

 To pin a metric to the top of the list, hover over the metric name and select the pin icon.
4. The Condition and Data sections at the bottom of the Metrics Visualizations screen offer further choices for the metrics visualization:

Condition

 - An experimental variable is required to generate a visualization. Select *Tag group*, *Stimulus Variable*, or *Participant Variable*.
 - After you select an experimental variable, then you can select the specific condition to visualize.
 - For binned metrics, adjust the bin size. The default value for each TOI is based on the interval durations in the TOI the first time it is selected.

Data

 - Select the Central tendency measure to visualize: *Median* or *Mean*.

9.6.2 View a Metrics Visualization

The metric visualization is updated automatically.

To view details, hover the mouse cursor over the following parts of the visualization:

Data point: Value, Recording, Participant, Interval, AOI, and Media (when applicable)

Central tendency line: Median or Mean, Sample Size

For binned metrics:

- Central Tendency points - Condition, Bin, Median or Mean, Sample Size
- Bin number or Bin tick - Time interval for bin

9.7 Times of Interest (TOI)

Times of Interest (TOIs) are divided into two general types.

- Automatic TOIs (sometimes called “System-generated TOIs”) are based on automatic events, such as recording, media, web navigation and snapshot events.
- Custom TOIs are created manually by the user and can be based on system generated events or custom events. .

Times of Interest is a highly flexible data analytic tool in Pro Lab. Together with Areas of Interest, they provide useful and fine-grained capabilities for defining not only the spatial extent of your analyses (AOIs) but also the temporal span (TOIs). Used appropriately and with care, researchers can apply these tools to carry out powerful, sophisticated analyses of the most demanding stimulus presentations.

Examples of TOI use:

- selecting data for visualizations such as heat maps and gaze plots
- selecting and aggregating data to calculate metrics associated with a coded task or sub-task periods in your recording (recording task or trial analysis)
- selecting and aggregating data to calculate metrics associated with coded subject behaviors or actions (behavioral coding)
- selecting and aggregating data associated with media exposure and snapshot coverage (media or snapshot restricted task or trial analysis)

Times of Interest are made up of intervals. Intervals are made up of events. Depending on the number of occurrences of Start and End Events, the TOI can consist of one or more intervals in one or more recordings. When you select a TOI for analysis, Pro Lab searches all recordings in the project and uses the gaze data intervals from the appropriate gaze data source to generate the visualization or to export metrics and data.

Intervals

Intervals are spans of time on the recording timeline that have a start and end event. Just as with events, intervals can be generated automatically (using automatic events). For example, there will always be an interval that corresponds to the entire duration of the recording from start to finish. If you are running a screen-based study where images are shown to the participant, then the system will automatically generate intervals that span the period of time when the image appears on the screen. Data on all intervals are available in the metrics output. Statistics on the duration and start/end point time stamps of all intervals can be obtained just by checking a box.

Events

Events are multipurpose markers that are used to identify important timestamps in the recording. Events can be generated automatically by Pro Lab (RecordingStart, ImageStimulusStart, sync events), the participant (key press), or the researcher/analyst (Custom events). Events have an associated time stamp, the exact time the event marker was applied. They can be counted but they do not have any duration since they simply mark a meaningful point in time. In Pro Lab, an Event Type is a class of Events and any instance thereafter is referred to simply as an Event. Once you have a few events, you can do more than just count them and that brings us to the next key concept, intervals.

9.7.1 Automatic Times of Interest

Depending on the type of project you are running, Pro Lab will automatically generate different TOIs. These are Automatic TOIs. The table below shows which types of automatic TOIs are created in recordings of different project types.

Automatic TOI	Intervals	Event types	Screen project	Glasses project	Scene Camera project	External Presenter project
Recording TOI	whole recording	RecordingStart, RecordingEnd	yes	yes	yes	yes
Media TOI	images + videos	ImageStimulusStart, ImageStimulusEnd, VideoStimulusStart, VideoStimulusEnd	yes	no	no	yes
Snapshot TOI	intervals (TOI) mapped to a snapshot	Snapshot (name) Interval Start, Snapshot (name) Interval End	yes	yes	yes	yes
Web navigation TOI	unique URL (webpage visits)	WebStimulusStart, WebStimulusEnd, URLStart, URLEnd	yes	no	no	no

Recording TOI

The Recording TOI exists for every recording in every project type. It includes all data from the start of the recording to the end of the recording.

MediaTOI

Media TOIs are created automatically in external presenter projects and screen-based projects for image, video, and text stimulus.

Snapshot TOI

Snapshot TOIs are created automatically when eye tracking data is mapped onto a snapshot. A snapshot TOI is composed of all the intervals defined by the two automatic events "Snapshot [name] Interval Start" and "Snapshot[name] Interval End". The "Snapshot [name] Interval Start"

event is created by the first gaze point/fixation mapped onto the snapshot. The "Snapshot[name] Interval End" event is created by the last gaze point/fixation mapped onto the snapshot. If there is a gap of 5 s between two mapped gaze points, a "Snapshot[name] Interval End" event is created for the last gaze point/fixation before the gap, and a "Snapshot[name] Interval Start" event is created on the gaze point after the gap.

Web Navigation TOI

Web Navigation TOIs are created automatically for each URL that the participant visits. The intervals are all the different visits to the same URL. A screenshot is captured and stored for each visit. The first screenshot captured will be the one associated with the TOI. The URL (web address) for the web Time of Interest is the default name. You can right-click the Web Navigation TOI and select *Copy URL* or *Open in browser*.

How to change to a different screenshot:

1. Right-click the Web Navigation TOI and select *CloneTOI*.
2. Name the Web Navigation CustomTOI and verify the start and end points.
3. Select a screenshot.



It is very important to match the resolution with the resolution of the snapshot captured by Pro Lab.

4. Click **OK**.

9.7.2 Custom Times of Interest

Custom TOIs are used to specify custom intervals for a visualization, metric calculation, or data export. When you create a custom TOI, you define the interval selection rules and the gaze data source.

A custom TOI's intervals are defined by Start Event types and End Event types with optional offsets. The Event types that can be used to define the start and end of intervals are generated automatically (eg. RecordingStart and RecordingEnd) or are logged manually. For more information, read [Manage Event Types](#).

A Custom TOI's gaze data source will be the recording data, mapped data, or web data.

How to create a Custom Time of Interest:

1. On the **Analyze** tab, make sure that you have defined the events you need. If necessary, create them using Custom Events. For more information, read [Manage Event Types](#).
2. In the **Times of Interest** section of the tool panel, click the plus sign (+).
3. Name your custom Time of Interest with a descriptive name in the **Name** field.
4. Select the gaze data source:
 - Select *Recording data* to create a standard Custom TOI. In a screen-based project, recording data is mapped to the Display Area Coordinate System. In a Glasses project, it is mapped to the Head Unit Coordinate System. For more information, read

Head Unit Coordinate System (HUCS).

- Select *Mapped data* to create a Snapshot Custom TOI. Mapped data is mapped to a snapshot using Assisted or Manual mapping. This creates a new coordinate system, the Media Coordinate System. For more information, read [Scene Camera Projects and the Media Coordinate System \(MCS\)](#).
 - Select *Web data* to create a Web Navigation Custom TOI. Web data is mapped to web screenshots.
5. Select the Event Types that define the start of your interval using the expandable lists under **Start point**.
 6. Select the Event Types that define the end of your interval using the expandable lists under **End point**.
 7. To add an optional offset to an event, toggle the **Offset** switch on. For more information, read [Add offsets to Custom Times of Interest](#).
 8. Select the appropriate media for the custom TOI.
 - Hover over the thumbnail or double click it to see a larger view of the media.
 - Hover over the title to see the full title, and resolution of the media. Double click the title to change it.
 - Web data only - Screenshots associated with the selected Event Types are shown with a link icon. Manually captured screenshots are shown with an M icon. Imported screenshots are shown with a downward facing arrow icon.
 - Web data only - Select the plus button (+) to import a new screenshot. If you import a screenshot, it is very important to match the resolution with the resolution of the snapshot captured by Pro Lab.
 9. Click **OK**.

Edit or delete a custom TOI

To delete or edit a Custom TOI, click the menu (vertical dots) next to the name and select *Delete TOI* or *Edit TOI*. To show or hide a TOI under the timeline, click the eye icon next to the name.

9.7.3 Add offsets to Custom Times of Interest

If you want to analyze data before or after one specific event, you can create a Custom Time of Interest with one event as both start and stop point and add offset time as buffers to extend the time window. You can also use this procedure to tighten the time window. This is useful, for example, if you want to check the closure of a test subject's action without having to watch their initial considerations.

- You remove time after the start point and add time after the end point by setting positive offset values.
- You add time before the start point and remove time before the end point by setting negative offset values.

How to add offsets:

1. In the **Times of Interest** section in the right-hand panel on the **Design** tab, double-click the name of the TOI to which you want to add an offset.
2. Select the desired event types.
3. Toggle on the **Start point** and/or **End point** offset buttons.
4. Set the desired offsets in seconds and milliseconds (limits are +/-100s). Note: The timer displays only when the buttons are set to “on.”
5. Click **OK**. The offsets are applied to all event types that you selected in step 2.
6. Check that the offsets have been properly set by taking a close look at the regular visualizations.

9.8 Areas of Interest (AOI)

Areas of Interest (AOIs) enable numerical and statistical analysis based on regions or objects of interest in your stimuli. AOIs come in different forms:

Manual AOIs refer to the ones the user creates herself by drawing an ellipse, rectangle, or polygon covering an area in the stimulus or Snapshot. Manual AOIs can be either static or dynamic. Static AOIs are the normal kind, and these maintain the same position and shape during the exposure of the stimulus. Dynamic AOIs, however, change in size and shape over the stimulus exposure because the stimulus changes the shape during the exposure. This is normally done in videos, where an object of interest moves around in the video or possibly disappears and reappears again.

Automatic AOIs refer to AOIs that Pro Lab creates by having intricate knowledge of the objects presented on the screen, in this case characters, words, and sentences created in the text.

9.8.1 Draw AOIs on Snapshots, image stimuli, or frames saved as media

(this section applies to Manual AOIs)

How to draw AOIs:

1. Activate the AOI tool by either clicking the **AOI Tool** button at the top on the **Project Overview** tab or select it from the **Analyze** drop-down menu in the top navigation.
2. Use the Media Selection panel on the right to open the snapshot, or media on which you want to draw an AOI.
3. Draw an AOI by selecting one of the following:
 - **Polygons:** Select the **Polygon** tool in the AOI toolbar above the AOI Visualization or press P on the keyboard. Click where you want the AOI to start and then click wherever you want to add vertices. Close the shape by clicking on the first vertex.
 - **Rectangles or ellipses:** Select the **Draw rectangle** or **Draw ellipse** tool in the AOI toolbar above the AOI Visualization, or press E (ellipse) or R (rectangle) on the keyboard. Click and drag the mouse cursor on the image or Snapshot until you have cre-

ated the shape and size you want. Make the shape a perfect square or circle by pressing the **Shift** key while you drag to draw.

How to name an AOI:

Select the AOI in one of the following ways:

- Go to the list of AOIs, double-click the current AOI name and enter the new name.
- Go to the AOI shape, double-click the name label and enter the new name.
- When an AOI is selected in the list of AOIs, press F2 on the keyboard. The “Rename area” dialogue appears.
- Enter a new name in the Input field and click OK.

9.8.2 Dynamic AOIs

The difference between a static AOI and a dynamic AOI is that the dynamic AOI can be moved around to match an object’s location in a video. The shapes and behaviors of dynamic AOIs in Tobii Pro Lab are defined by Keyframes. A Keyframe is created when you add an AOI in one of the frames. Dynamic media containing dynamic AOIs will typically have numerous Keyframes for each AOI. In-between these Keyframes, Pro Lab interpolates the shape and position of the AOI, so that the AOI moves smoothly from one Keyframe to the next. Any changes that are made to an AOI somewhere on the timeline, for example dragging an AOI or moving AOI vertices, will create a new Keyframe unless you are changing an existing Keyframe.

9.8.3 Edit the shape of an AOI

(this section applies to Manual AOIs)

How to change the shape of an existing AOI:

1. Make sure that the desired AOI is visible. (If it’s not, read [Change the name, appearance, or visibility of an AOI.](#))
2. Go to the AOI toolbar above the AOI Visualization and select the *Select/Move Vertices* tool or press V on the keyboard.
3. Drag the vertices to transform the AOI.



Currently, there is no way to add or delete vertices on existing AOIs.

9.8.4 Change the size of an AOI

(this section applies to Manual AOIs)

How to change the size of an AOI:

1. Make sure that the desired AOI is visible. (If it’s not, read [Change the name, appearance, or visibility of an AOI.](#))
2. Go to the AOI toolbar above the AOI Visualization and select the *Select/Move AOI* tool or press 2 on the keyboard.
3. Click on the AOI that you want to resize. A rectangle with sizing handles in the corners appears around the selected AOI.

- Increase or decrease the size by clicking and dragging one of the sizing handles.
- Maintain the object's proportions by pressing the *Shift* while you drag the sizing handle.
- The size of a selected AOI can also be adjusted by modifying the width (W:) and height (H:) values in the property fields on the AOI toolbar.

9.8.5 Move an AOI

(this section applies to Manual AOIs)

How to move an AOI:

1. Make sure that the desired AOI is visible. (If it's not, read [Change the name, appearance, or visibility of an AOI.](#))
2. Go to the AOI toolbar above the AOI Visualization and select the *Select/Move AOI* tool or press S on the keyboard.
3. Click on the AOI that you want to move. A rectangle appears around the selected AOI.
4. Drag the AOI to its new location.
 - You constrain an AOI so that it moves only horizontally or vertically by pressing the Shift key while dragging the object.
 - The position of a selected AOI can also be adjusted by modifying the X: and Y: values in the property fields on the toolbar.

9.8.6 Cut, copy or paste AOIs

(this section applies to Manual AOIs)

An AOI can be cut or copied from one image or Snapshot and pasted onto another Snapshot, as well as into the same image or Snapshot as an identical copy of the original (apart from the name).

How to cut or copy the shape of an AOI:

1. Select the *Select/Move AOI* tool on the AOI toolbar.
2. Right-click on the AOI that you want to cut or copy and select *Copy* or *Cut* from the menu. You can do this by right-clicking either the AOI shape on the image/ Snapshot or the AOI name in the Areas of Interest list in the tool panel to the right.

How to paste the AOI:

1. Paste the AOI in the middle of the image or Snapshot by right-clicking on the Snapshot and selecting *Paste*.
2. Paste the AOI in the same location where the copied AOI was by right-clicking on the Snapshot and selecting *Paste in Place*.

The quickest way to duplicate an AOI is by pressing *Ctrl* on the keyboard while clicking and dragging the desired AOI to a new location. You can also copy and paste AOIs by using the Edit menu on the AOI toolbar or by using the keyboard shortcuts Ctrl+C (to copy), Ctrl+X (to cut), and Ctrl+V (to paste).

9.8.7 Change the name, appearance, or visibility of an AOI

(this section applies to Manual and Automatic AOIs)

How to change the name of an AOI:

- Go to the list of AOIs in the tools panel on the right, double-click the current AOI name and enter the new name.
- Go to the AOI shape, double-click the name label and enter the new name.
- Select an AOI in the list of AOIs in the tools panel on the right and press F2 on the keyboard. The “Rename area” dialog box appears. Go to the input field, enter a new name and click OK.

How to change the appearance of an AOI:

- Change the color of an AOI by clicking the AOI color indicator rectangle next to the AOI name in the *Areas of Interest* list. Pick a new color and click **OK**.
- Change the opacity of the AOIs by adjusting the AOIs *Opacity* slider in the *View Options* section in the tool panel on the right.
- Display or hide the AOI name labels in the AOIs by checking or unchecking *Show AOI Names* in the *View Options* section in the tools panel on the right.
- Display or hide the AOI fill colors by toggling on/off *Show Fill Color* in the *View Options* section in the tools panel on the right.

9.8.8 Undo or redo AOI actions

(this section applies to all AOI types)

How to undo and redo actions:

- Undo your most recent action by clicking the **Undo** button in the top left corner of the AOI tab interface or by pressing Ctrl+Z. Click the arrow (triangle) to see the drop-down list of earlier actions to undo.
- Redo the action you have just undone by clicking the **Redo** button in the top left corner of the AOI tab interface. Click the arrow (triangle) to see the drop-down list of earlier actions you want to redo

9.8.9 Pan and zoom in the AOI editor

(this section applies to all AOI types)

When you work with large images or Snapshots (e.g. Snapshots that cover an entire shopping shelf stitched together by multiple image files), you might find it difficult to see the details. In such cases, it's useful to zoom in to parts of the interface.

How to zoom in to an image or Snapshot:

- Select the **Zoom** tool that looks like a magnifying glass in the Navigation section of the AOI toolbar. The cursor turns into a magnifying glass. Click the image or Snapshot where you want to zoom in.

- Click the **Zoom** tool that looks like a magnifying glass below the AOI Visualization. This opens a zoom slider. Drag the slider to the zoom percentage setting that you want.
- Press Ctrl on the keyboard while scrolling the mouse wheel to zoom in or out. The center of the zoom will be where you have the cursor.
- Click the buttons labeled “Fit to window” or “Actual size” located under the AOI Visualization to switch between the two zoom levels.

When you have zoomed into an image or Snapshot, you cannot always see it in its entirety. It is then useful to pan your field of vision so that you get the desired area in view.

How to pan in the image or Snapshot:

1. Select the *Move Canvas* tool that resemble hand on the *Navigate* section of the AOI toolbar.
2. Drag the Snapshot image to pan it to a new location.

9.8.10 Activate or deactivate an AOI

(this section applies to Dynamic AOIs)

When an AOI is active it collects and records data during the replay. You can choose to deactivate an AOI for any period on the timeline.

How to activate or deactivate an AOI:

1. Select the AOI that you want to activate/deactivate.
2. Right-click on the AOI and select *Activate/Deactivate Selected AOIs* in the context menu. You can also use the *AOI Active* toggle switch in the function bar above the recording.

9.8.11 Export an image of the AOIs on the image or Snapshot

(this section applies to all AOI types)

How to export an image of the AOIs on the image or Snapshot:

1. Right-click on the image or Snapshot in the AOI Visualization.
2. Click *Export Image* in the context menu.
3. In the *Save As* dialogue, enter a file name in the *File name* input field.
4. Navigate to the folder into which you want to save the image and click *Save*.

9.9 AOI tags

(This section applies to all AOI types)

In Pro Lab you can assign Tags to both static and dynamic AOIs. This makes it possible to aggregate data from several AOIs both within and between media for conditions that exist across multiple stimuli elements. This, in turn, creates an easier workflow for many experimental paradigms and research questions. Tags can be used to create Aggregated AOIs which gives you more flexibility in combining different AOI attributes for your Metrics and Data export. Read more about [Aggregated AOIs](#).

Tags can function as a quick selector in case of a large number of AOIs. Read how to [Export metrics data to a file](#).

Tags can also be used to carry information about the AOI or its content into the analysis in the form of metadata. Tags are supported in both Metrics and Data export.

9.9.1 Create AOI Tags

1. Select "AOI Tool" from the drop-down menu at the top next to **Analyze**(the Analyze module).
2. Open AOI Tags Manager by clicking the **Create and edit Tags** button on the **Tags** panel.
3. Click the plus sign to the right of *Ungrouped Tags* to create a tag.
4. Type a descriptive name.
 - ✔ If you want to remove a Tag, hover the mouse cursor over it and click the x symbol.
5. Continue creating tags as needed.




Figure 12. Create tags

Assign Tags to AOIs

When you have created your Tags, then you can assign them to the AOIs in the stimuli.

1. Click the AOI for which you want to assign the Tags. Select several AOIs by Shift-clicking or by clicking a point in the stimuli and then drag a selection box to include the wanted AOIs.
2. Click the Tags you want to associate to the AOI/AOIs on the Tags panel.



Grouped Tags are mutually exclusive, i.e. only one Tag per group can be selected for the selected AOI/ AOIs; but several un-grouped Tags can be selected independently of each other.

3. The AOI now reflects the Tag/Tags that have been selected for it.



Figure 13. In this illustration the *Face1* label is the AOI name, *Mood* is the group label and *Angry* is the tag name.



If you have many tags selected for the AOI, the groups and tag labels may clutter and obscure the AOI. For this reason you can hide the labels by switching off the *Show AOI Tag Names* selector switch in the *View Options* panel. You can still display the tag names by hovering over the “#” symbol on the AOI label.

Create Tag Groups

Groups can be considered the top designator of Tags. For example, you could have the group *Vehicle type* with the Tags *Car*, *Motorcycle* and *Bicycle*, and another group *Brand* with the Tags *Volvo*, *BMW*, *Audi*, *Honda*, *Yamaha*, *Cannondale* and *Specialized*. This way you can export metrics for different types of transportation means, and also drill down by selecting the Tags that denote different vehicle brands.

How to create Tag Groups:

1. Open AOI Tags Manager by clicking the **Create and edit Tags** button on the **Tags** panel. The *AOI Tags Manager* dialog appears.
2. Click the **New Tag Group** button and name the group in the name field.
3. Create as many groups as you need. You can always return and create/edit/delete groups later.
4. Create and add Tags to groups by clicking the Plus (+) sign to the right of the group.
5. Name the Tag when it's added to the group or accept the existing name.



If you want to remove a group, hover the cursor over its name and the **Delete** button (trash can) displays next to the name. If the group has Tags assigned to it, you need to confirm that you want to remove the group and all the associated Tags with it.

9.9.2 Aggregated AOIs

Aggregated AOIs give more flexibility in collecting aggregated data from AOIs labeled by one or multiple Tags both within and between media. By combining the Tags into an Aggregated AOI, you can analyze the gaze data from AOIs that meet specific criteria.


An aggregated AOI is created by selecting Tags using the check boxes or specifying the relationship between them in a logical Expression with a set of available operators.

Suppose you are doing a shopper study and need to categorize your items by type and color. In this case, you have two independent attributes (Tag Groups) applied to your items. The first Group is called “Clothes” and has the following values: Shirt, Trousers, Socks. The second Group is called “Color” and has the following values: White, Blue, and Red.

Now imagine you need to export the gaze data on the Blue and the White Shirts. All the AOIs you have in the project will be labeled by Tags from both Groups at the same time, since they need to be categorized by both the type of clothes and the color. To do this, you need to create an Aggregated AOI “Blue&White Shirts,” select the respective Tags in both Groups and make sure the following logical expression is applied (this is the default):

Clothes:Shirt AND (Color:Blue OR Color:White)

This lets you export the data set that includes all the Shirts that are Blue or White.





Aggregated AOIs are project-wide, which means that creating, editing, or deleting an Aggregated AOI on the Data selection panel in **Metrics** or **Data Export** will respectively create, edit, or delete *the same Aggregated AOI* in the other section.


Create Aggregated AOIs



You must first create tags in the AOI tool in order use them to create aggregated AOIs.

How to create Aggregated AOIs with selected tags:

1. Click the **Metrics** or **Data Export** button either in **Project Overview** or select either one from the **Analyze** drop-down menu in the top navigation.
2. Click the Plus (+) sign next to *Create Aggregated AOIs* in the Aggregated AOIs section on the Data selection panel. The Create aggregated AOI dialog displays.



3. Name your Aggregated AOI with a descriptive name in the Name field.

4. Select tags to aggregate data from the tagged AOIs (tags that represent the attributes you want to use for aggregation). Both Tags in the Tag Groups and Ungrouped Tags are available for selection.
5. If necessary, edit the logical Expression applied to the selected Tags by clicking the Edit button (pen icon) to the right of *Edit expression*.

The following operators can be used:

- **AND** or **&&** - includes all the AOIs tagged by both Tags at the same time
- **OR** or **||** - includes all the AOIs tagged by one of the Tags
- **NOT** or **!** - includes all the AOIs not tagged by the Tag
- **(or)** - brackets are used for more complex expressions (follow rules of mathematical logic)

By default, all the Tags in a Tag Group will be put in brackets with the “OR” relationship between them. All the selected Tag Groups (with the selected Tags in brackets) as well as the Ungrouped Tags will have the “AND” relationship between them as in the formula below:

(Group1:Tag1 OR Group1:Tag2) AND (Group2:Tag3 OR Group2:Tag4) AND UngroupedTag5 AND UngroupedTag6



If the default expression has been changed and saved after editing, the Tag selection tree will not be available to edit (grayed out), which means that you will either have to manually type the Tag names in the expression or start over in order to create the Aggregated AOI.



If there is a syntax error in the expression, Pro Lab shows you the elements that cannot be interpreted and the Create button will be inactive.

6. Click **Create** to save the formula or click the **Close** button (X) to return to the default expression.

Edit an Aggregated AOI

1. Select the **Metrics** or **Data Export** button either in the Project Overview section or select it from the **Analyze** drop-down menu in the top navigation.
2. Hover the mouse cursor over the Aggregated AOI you want to edit in the **Aggregated AOIs** section on the **Data selection** panel.
3. Select the **Edit** button (pen).
4. In the Edit aggregated AOI dialog, you can edit the name, select/unselect the tags to aggregate data from the tagged AOIs, and edit the logical expression applied to the selected Tags.
5. When you are done editing, click **Save** or **Cancel**.



Aggregated AOIs are project-wide, which means that creating, editing, or deleting an Aggregated AOI on the Data selection panel in **Metrics** or **Data Export** will respectively create, edit, or delete *the same Aggregated AOI* in the other section.

Delete an Aggregated AOI

1. Click the **Metrics** or **Data Export** button either in **Project Overview** or select either one from the **Analyze** drop-down menu in the top navigation.
2. Hover the mouse cursor over the Aggregated AOI you want to delete in the *Aggregated AOIs* section on the Data selection panel.
3. Click the **Delete** button (trash can).
4. Click **OK** to confirm the deletion.



Aggregated AOIs are project-wide, which means that creating, editing, or deleting an Aggregated AOI on the Data selection panel in **Metrics** or **Data Export** will respectively create, edit, or delete *the same Aggregated AOI* in the other section.

Export data based on Aggregated AOIs

When you export data, you can select to export only data (or **Metrics**) that are relevant to your Aggregated AOIs .

How to select export data:

1. Click the **Data Export** button either in the **Project Overview** section or select it from the **Analyze** drop-down menu in the top navigation.
2. Expand the *Aggregated AOIs* section on the *Data selection* panel.
3. By default, all the Aggregated AOIs are selected. If you want to select only some of the Aggregated AOIs, deselect everything by clicking the All Aggregated AOIs checkbox. Select the Aggregated AOIs you want to include in your Data Export.
4. Select/deselect any other data related to what you want to export by checking and unchecking the metrics in the "Select data for export list" on the left.
5. Click the blue **Export file(s)** button at the top when ready to export.



The metrics export for Aggregated AOIs will not be affected by the selected Areas of Interest on the Data selection panel. The metrics export will still be based on all AOIs for the selected Aggregated AOIs.

9.9.3 Export and import AOIs

(this section applies to Manual AOIs)

If you want to use the same AOI again with a different image, you can export AOIs and then import them into another project.

Any tags added to the AOIs will be exported/imported along with the AOIs and if you have added them they will be analyzed together (based on the assigned AOI tags). For more information, read [AOI tags](#).



All AOIs in the media will be exported and imported when you select Export AOIs or Import AOIs.

How to export AOIs:

1. Open the AOI Tool from the **Analyze** drop-down list at the top of Pro Lab.
2. Select the media in the Media Selection panel.
3. Right-click the media in the main viewer and select *Export AOIs*.
4. Save the file.



When no AOIs exist, the *Export* option is disabled.

How to import AOIs:

1. Open the AOI Tool from the **Analyze** drop-down list at the top of Pro Lab.
2. Select the media in the **Media Selection** panel.
3. Right-click the media in the main viewer and select *Import AOIs*.
4. Select the file using the file browser.
5. The AOIs will be imported and shown on the media in the main viewer.

9.9.4 Automatic AOIs (Text stimulus)

In addition to drawing Manual AOIs, Automatic AOIs can be generated for a Text stimulus. Select the Text stimulus for which you want to generate Automatic AOIs in Media Selection. The **Select media** section will then be active. Select the desired reading AOI level (Characters, Words or Sentences) to generate Automatic AOI and click **Generate**.




Figure 14. Select the level of AOI and click Generate.




Figure 15. An example of words selected to generate Automatic AOIs

Characters: any type of letter or punctuation mark is one character.

Words: characters divided by [SPACE] are automatically interpreted as a word. Note that punctuation marks like a dot or a comma will count as a part of the Word AOI.

Sentences - A sentence is an array of words ending with a period, question mark, or an exclamation mark (and a space if there are other words after it). You can also create a sentence by using a line break (pressing ENTER).



The Chinese character "。" also marks an array of words combining to form a sentence (and a space if there are other words after it).

The AOIs will be created and visible once they are generated. You can hide them or change the color individually for each AOI the same way you would do for Manual AOIs.

To remove the Automatic AOIs click **Remove**. All Automatic AOIs on the selected level will be removed from this Text stimulus. Automatic AOIs can be regenerated and removed as many times as needed.



Automatic AOIs cannot be imported or exported because they are dependent on the text.



AOI tag groups and tags, as well as Aggregated AOIs work for Automatic AOIs.

For more information, read [Text stimulus \(Reading studies and instruction\)](#).

9.10 Analyze a Screen Recording

You can record a screen by using a Screen Recording stimulus. Adding a Screen Recording stimulus is explained in [Timelines](#) and you can learn more about editing Screen Recording stimulus properties in [Stimulus properties](#).

- When you want to analyze a screen recording, you can log events, just like for any other type of recording. For more information, read [Replay tools](#).
- When you want to aggregate data and create visualizations such as heat maps and gaze plots, you can create TOIs with frames. For more information, read [Replay tools](#).
- When you want to analyze a screen recording, you can create AOIs and dynamic AOIs in the AOI Editor just like for other types of recordings. The name of the recording in the AOI Editor will be "Screen Recording + recording name". For example, if you made a screen recording with the name Version A, the recording name in the Media Selection in the AOI Editor is "Screen Recording Version A". For more information, read [Dynamic AOIs](#).

10 Metrics

To set up a successful eye tracking study you need to define and calculate the appropriate measures for your research question. In addition to choosing the right eye tracking measure, you need to define where and when to calculate this measure, i.e. the Areas of Interest (AOI) that are associated with operationalization of your research question. You also should calculate Times of Interest (TOI), the intervals of the recording when your stimulus or behavior of interest are predicted to occur. Some examples are: the duration of the exposure of a stimulus on the screen, a section of a trial, the time between when a stimulus appears on the screen and when a participant presses a key on the keyboard, the moment someone enters a supermarket aisle and places a product in the shopping basket, etc.

In Pro Lab, the term "metric" is used to define the different measures that are calculated from the recording data. These measures can be exported in different table/file formats that can either be used to get an overview of the data and extract summary statistics, or to organize the data for further processing in statistical software platforms such as R or SPSS.

For best practice, and unless your study is an explorative one, measures should be defined during the planning and design phase of the study.

10.1 Understand the metrics

AOI-based metrics vs. AOI-independent metrics

Some metrics in Pro Lab relate to a particular area of interest (AOI) for example, how much time a participant spends looking at advertisements, or the eyes of a person talking. Other metrics are general for that time of interest interval, for example the average fixation duration during the interval, regardless of what the user specifically looked at. The former says something about the object in the area of interest, but the latter can put this number in context of how the participant behaves in general and when seeing everything else that is visible during this stimulus presentation.

If you select an AOI-based metric, you will get one metric computed per AOI. If it is a saccade metric, then the AOI will determine if a saccade is an "entry saccade" (moving into the AOI) or an "exit saccade" (moving out of the AOI).

However, if you choose a metric that is AOI-independent, then you will get a metric that is computed for the entire interval. For example, an average fixation duration that is AOI-independent will give you the average fixation duration of all fixations occurring in that TOI interval, regardless of what they looked at.

Partial and whole events

Previously Pro Lab offered only fixation-based metrics and these were very inclusive, in that they included all fixations in that Time of Interest (and Area of Interest). All fixations means both whole and partial fixations.

Whole events

An event that starts and stops within the given interval, and is not disrupted by unknown events or data loss at the beginning or end of the event. If any of these two criteria are fulfilled, then it is counted as a whole event.

Partial events

Events that do not fulfill the criteria for whole events are partial events.

Examples of partial events:

- the saccade starts before the interval starts, and stops inside the interval. In this case, there is ambiguity in how to treat this event. It has data outside of the interval you want to analyze, so it is not given that data outside the interval should be considered at all. At the same time, a key part of that event, e.g. the peak velocity, could happen outside of the interval yet be relevant to the entire event.
- A fixation starts before the interval starts. The duration of the fixation is ambiguous as a result. Should it be the time of the fixations that it has spent on the stimulus in question (current interval), or should it be the entire duration of the fixation?
- A fixation is followed by a saccade, but this saccade is interrupted by data loss. In this case, the end of the saccade is uncertain, because we cannot know whether the data samples lost were part of the saccade, or the subsequent event. Considering this ambiguity, this saccade is labeled as a partial event.

For eye movement researchers interested in fixations and saccades as such, it is important to distinguish between whole and partial events in order to exclude partial events from their studies.

For example, a whole saccade had a certain amplitude (e.g. 20 degrees). If it is cut by the TOI interval start (producing a partial saccade) the amplitude may show 8 degrees, it can't be used and should be discarded.

On the other hand, for researchers interested in attention and visual exposure rather than eye movements per se, it is perfectly fine to include partial events.

For example, if an advertisement is suddenly shown (new stimulus) it will be seen and the event can be used, even if the stimulus-interval cuts an ongoing fixation.

Criteria for whole saccades

In order for a saccade to be labeled whole it must fulfill the following criteria:

- A saccade must be preceded and succeeded by a fixation event.

In addition, for the Interval-based export format, partial events pose additional problems as the metrics aggregate all events in an interval. For example, if the amplitude of a partial saccade (cut by the border of the TOI) is included in the calculation of "average amplitude of saccades", the result will be a lower average. This may lead to the wrong conclusion about the average length of saccades in the interval.

To avoid this risk, the following criteria applies for saccades in the interval-based TSV export format:

- A saccade must be wholly contained in the time of interest.

In the Interval- and AOI-based export format partial saccades are excluded from the metric calculation. In the Event-based export format partial saccades are shown as "partial" in the validity column. It is up to the individual researcher to exclude these events if needed.

Criteria for whole fixations

Some researchers are interested in both the saccades and fixations that occur in the same TOI. For example, research suggests that there are phases of visual exploration, with an early orienting phase with high-amplitude saccades and short-duration fixations, followed by an inspection phase with low-amplitude saccades and long-duration fixations. To explore this in a meaningful way, we want fixations to follow the same strict definition as saccades.

Therefore, we also allow the export of “whole fixations”(i.e. "exclude partial fixations"). must fulfill the following criteria:

- A whole fixation must be preceded by a saccade.
- A whole fixation must be succeeded by a saccade.
- A whole fixation must be wholly contained in the time of interest

Finally, for AOI metrics that exclude partial fixations there is an additional criteria.

- The whole fixation duration must be contained within the AOI

These metrics are found in the Interval- and AOI-based export format.

10.2 Export metrics data to a file

You can select the metrics to export in order to create custom spreadsheets or “reports”. Different metrics are available for export depending on which export format you choose. For more information, read [Export metrics data to a file](#) and [Understand the metrics](#).

Exporting metrics



If the Export button is inactive (grayed out), hover the cursor over the **Export** button to see why it is inactive. For example, the tool tip will tell you that at least one recording, one metric, and one TOI must be selected to start the metrics export.

1. Click the **Metrics** button either in the **Project Overview** section or select it from the Analyze drop-down menu in the top navigation.
2. Select the output format for the export file by clicking the **Export Format** drop-down list in the **Settings** panel. Select one of the three formats. For more information, read [Export metrics data to a file](#)
3. If applicable, toggle on **Binning** in the **Settings** panel. Read more about [Binned metrics](#).
4. Select which Gaze Filter you want to apply to the data you will export. For more information, read [Gaze Filter functions and effects](#).
5. Select the data you want to include using the Data Selection section on the right-hand side. The options are:
 - **Recordings**: Choose to export metrics from All Recordings or only the recordings of your choice. Expand the checkbox to see all available recordings.

- **Participant variables:** Choose to export metrics based on a selection of participant variables. Use the checkboxes to select which participant variables the export will be filtered on. All checkboxes are selected as default.
- **Times of Interest:** Choose to export metrics from the complete recordings (selected above) or metrics only from selected TOIs. To simplify the selection, there are two groups of TOI data: Media gaze data, where automatically created TOIs based on when Snapshots or images were active are presented, and Recording gaze data, where your custom TOIs are listed. In the Media gaze data list, each image or Snapshot associated with the project that has AOIs is listed. For each of them, you can also select TOIs which include data from when they were active. Expand the checkbox next to each group to see and select the individual TOIs.
- **Areas of Interest:** Choose to export data relating to all or some of the AOIs created in the project using the checkboxes.
AOIs can also be selected by tags. To do this, select the tags in the drop-down menu "Select by tag" in the Areas of Interest section. All the AOIs associated with the selected tags will be selected automatically in the AOI selection tree. Deselecting one or several of the automatically selected AOIs will remove those tags from the list of selected tags in the Select by tag section.
Note that removing tags from the list will deselect the associated AOIs even if they were selected manually before using the "Selected by tag" option.





Figure 16. Click the drop-down arrow to select by tag in the Areas of Interest side panel

- **Aggregated AOIs:** Choose to export data related to all or some of the Aggregated AOIs created in the project using the checkboxes. Read more about how to [Create Aggregated AOIs](#).
 - **Events:** you can select the recording custom events to include in the Event and interval metrics.
6. Select the metrics you want to include in the export (the selection varies depending on what you chose in steps 2 and 3).

7.  Use the **Select All** toggle to select or deselect all the metrics. Use the up and down arrows in the metric export selection list to collapse (hide) entire sections in the list.
8. Click the **Export** button at the top of the Metrics Export window.
9. In the file browser, locate the folder in which you want to save the file.
10. In the File name field, enter a file name for your exported file.
11. Click **Save**.

10.3 Binned metrics

Sometimes you might want to have more granularity in the metrics in order to see how something unfolds over time in a given interval. This can be achieved by dividing the interval into smaller arbitrary size sections (time bins) and exporting the metric values for each bin falling into a given interval instead of the metric value representatives of the interval as a whole. This process (metrics output type) is referred to as Binning (Binned metrics).

Every bin for every interval for every recording gets its own row.

Binning is available for metrics export in the Interval-based, AOI-based, and Event-based TSV file formats.

To enable Binned metrics, switch on Binning using the toggle in the Settings panel. The Bin size window becomes active. The default bin size is 200 ms (milliseconds). Type the desired value in milliseconds and press ENTER. The available bin size values are between 1 millisecond and 900 000 milliseconds (15 minutes) with the step of 1 millisecond.

The available metrics differ depending on if Binning is toggled on or off. This is reflected in the separate tables for the metrics listed in the [Interval-based TSV file](#).

10.4 Reading metrics

In the example below, you can see a sample gaze path (fixation sequence) and a table showing the AOI word metrics available in Tobii Pro Lab (starting with version 1.162) for this specific sentence example: "The brown fox jumped over the fence despite me shouting."

- *(blank)* - means that the value will be an empty cell when the metrics are exported.
- *Numeric values* (1,0, ...) refer to the actual value in the metrics export file.
- *Fixation bubbles* (① ②) mean that, in order to calculate the metric value, you need to add the duration of the fixations of the respective order number.
- *String metrics* are the actual values in the metrics export file.

The "Total duration of fixations" metric in the AOI reading metrics group (as seen in the table below) differs from the "Total duration of fixations" in the AOI fixation metrics group in the following ways:

- It can be calculated for Automatic AOIs (generated on a Text stimulus) only
- It equals "0" for skipped AOIs and gives blank values for AOIs that were never reached




Figure 17. Sample gaze path (fixation sequence)

Metric name										
Character index	(blank)	(blank)	(blank)	(blank)	(blank)	(blank)	(blank)	(blank)	(blank)	(blank)
Word index	1	2	3	4	5	6	7	8	9	10
Sentence index	1	1	1	1	1	1	1	1	1	1
AOI string	The	brown	fox	jumped	over	the	fence	despite	me	shouting.
Text unit type	Word	Word	Word	Word	Word	Word	Word	Word	Word	Word
Number of units	3	5	3	6	4	3	5	7	2	9
First-pass first fixation duration	1	2	3	5	7	0	8	0	11	(blank)
First-pass duration	1	2	3+4	5	7	0	8	0	11	(blank)
Selective regression-path duration	1	2	3+4	5	7	0	8+9	0	11	(blank)
First-pass regression	0	0	0	1	0	(blank)	1	(blank)	1	(blank)
Total duration of fixations	1	2	3+4+5+9	5	7	0	8+9	12	11	(blank)
Regression-path duration	1	2	3+4	5+6	7	0	8+9+10	0	11+12	(blank)
Re-reading duration	0	0	0	6	0	0	9+10	0	12	(blank)

Figure 18. This table shows all the AOI reading metrics for the sentence example

10.5 Metric export formats

There are four kinds of Metric export formats in Pro Lab:

Interval-based TSV file

Use this format when you want to analyze your data at the TOI interval level. If you want to analyze the metrics data in a statistical analysis software, like R/SPSS/MATLAB etc. This format can be interpreted by the user, but is also especially designed for use in analysis software with meta information and metrics in a column format where the rows contain calculations for the actual Times of Interval.

AOI-based TSV file

This format is preferred for an analysis where you want to have AOI as a grouping factor and thus the AOI name as its own column and gaze metric in their own columns. Every AOI for every interval for every recording gets its own row.

Event-based TSV file

Use this format when you want to analyze individual events during an interval. Each fixation will generate a row in the report and the selected metrics will be shown as columns. Like the Interval-based TSV file it is formatted to be easily to import and analyze in statistical analysis software.

Excel report (.xlsx)

This format is compatible with most spreadsheet software such as Microsoft Excel (2007 and newer), Google Sheets, OpenOffice.org, etc. In this file, each metric is saved in a separate spreadsheet. Each image, Snapshot, or Time of Interest has its own table in the spreadsheet. The data in this export is highly aggregated and is intended to use "as is." It is therefore not the best choice for further analysis in statistical software platforms such as R or SPSS.

General information about export metrics

The metrics available for export in the file formats are shown in the tables below. An interval corresponds to one occurrence of a specific TOI. A TOI can occur multiple times during a recording which means there are multiple intervals. Event-based metrics format where one row corresponds to one eye-movement event, as identified by the Gaze Filter in the top-right corner of the Metrics Export view. Currently, only fixations are supported.

10.5.1 Interval-based TSV file

The metrics available for export in the Interval-based TSV file formats are shown in the table below.

An interval corresponds to one occurrence of a specific time of interest. The interval start is defined as the starting event for the TOI. The interval end is defined as the ending event for the TOI. A specific TOI can occur multiple times during a recording which means there are multiple intervals.



Byte Order Marks (BOM) flags are removed in .tsv files. If you have scripts that rely on this flag, be sure to update them.

General

Metric name	Description	Format
Recording name	Recording name	
Participant	Participant	
Participant variables	Variable value, or values, of the participant. One column for each participant variable.	
Timeline name	Timeline name	
TOI	The name of the current Time of Interest.	
Interval	The interval number of the current TOI interval.	
Media	The name of the media presented to the participant.	

Metric name	Description	Format
Stimulus variables	Stimulus variable value or values of the stimulus. One column for each Stimulus variable.	

General (Binning turned on)

Metric name	Description	Format
Recording name	Recording name	
Participant	Participant	
Participant variables	Variable value, or values, of the participant. One column for each participant variable.	
Timeline name	Timeline name	
TOI	The name of the current Time of Interest.	
Interval	The interval number of the current TOI interval.	
Bin	The index of the current bin in the interval.	
Bin duration	The duration of the current bin.	Milliseconds
Media	The name of the media presented to the participant.	
Stimulus variables	Stimulus variable value or values of the stimulus. One column for each Stimulus variable.	

Interval metrics

Metric name	Description	Format
Duration of interval	The duration of an interval.	Milliseconds
Start of interval	The start time of an interval.	Milliseconds

Interval metrics (Binning turned on)

Metric name	Description	Format
Duration of interval	The duration of an interval.	Milliseconds
Start of interval	The start time of an interval.	Milliseconds

Event metrics

Events can also be used in measures. Event metrics allow you to measure behavior and calculate statistics based on your event coding scheme.

Metric name	Description	Format
Number of Events	The number of Events, including Custom Events and Logged live Events, for an interval.	Count
Time to first Event	The time to the first Event, including Custom Events and Logged live Events, for an interval.	Milliseconds
Last key press	The last registered key press in the interval.	

Event metrics (Binning turned on)

Metric name	Description	Format
Event occurred	Any Event, including Custom Events and Logged live Events, occurs during a bin, indicated by 1/0, for each event type.	Binary
Number of Events	The number of Events, including Custom Events and Logged live Events, for a bin .	Count
Last key press	The last registered key press in the bin.	

AOI Fixation metrics

AOI fixations correspond to fixations that fall within an AOI. The fixations are defined based on the gaze filter you use (e.g. if you use the Raw gaze filter, every valid eye tracking sample is a fixation). AOI fixation metrics allow you to measure statistics based on the fixations within an AOI. They present as an interval (or an occurrence) of the TOI in separate rows in the exported spreadsheet.

Metric name	Description	Format
Total duration of fixations	The total duration of the fixations inside an AOI during an interval.	Milliseconds
Average duration of fixations	The average duration of the fixations inside an AOI during an interval.	Milliseconds
Minimum duration of fixations	The duration of the shortest fixation inside an AOI during an interval.	Milliseconds
Maximum duration of fixations	The duration of the longest fixation inside an AOI during an interval.	Milliseconds

Metric name	Description	Format
Number of fixations	The number of fixations occurring in an AOI during an interval.	Count
Time to first fixation	The time to the first fixation inside an AOI during an interval.	Milliseconds
Duration of first fixation	The duration of the first fixation inside an AOI during an interval.	Milliseconds
Last AOI viewed	The last AOI fixated during an interval.	
AOI at interval end	The AOI fixated at the end of an interval.	
Average pupil diameter	The average pupil diameter of all fixation samples in an AOI in an interval.	

AOI Fixation metrics (Binning turned on)

Metric name	Description	Format
Fixation hit	Any fixation hits inside an AOI, indicated by 1/0, for each bin.	Binary
Total duration of fixations	The total duration of the fixations inside an AOI during a bin.	Milliseconds
Number of fixation starts	The number of fixations inside an AOI that starts in the bin.	Count
Average pupil diameter	The average pupil diameter of all fixation samples in an AOI in a bin.	Millimeters

AOI Fixation metrics (exclude partial fixations)

These metrics exclude fixations that don't fulfill the criteria for whole fixations (see previous section).

Metric name	Description	Format
Total duration of whole fixations	The total duration of the fixations inside an AOI during an interval.	Milliseconds
Average duration of whole fixations	The average duration of the fixations inside an AOI during an interval.	Milliseconds

Metric name	Description	Format
Minimum duration of whole fixations	The duration of the shortest fixation inside an AOI during an interval.	Milliseconds
Maximum duration of whole fixations	The duration of the longest fixation inside an AOI during an interval.	Milliseconds
Number of whole fixations	The number of fixations occurring in an AOI during an interval.	Number
Time to first whole fixation	The time to the first fixation inside an AOI during an interval.	Milliseconds
Duration of first whole fixation	The duration of the first fixation inside this area of interest during an interval.	Milliseconds
Average whole fixation pupil diameter	The average pupil diameter of all whole fixation samples in an AOI in this interval.	Millimeters

AOI Fixation metrics (exclude partial fixations) (Binning turned on)

Metric name	Description	Format
Whole fixation hit	Any fixation hits inside an AOI, indicated by 1/0, for each bin.	Binary
Total duration of whole fixations	The total duration of the fixations inside an AOI during a bin.	Milliseconds
Number of whole fixation starts	The number of fixations inside an AOI that starts in the bin.	Count
Average whole fixation pupil diameter	The average pupil diameter of all whole fixation samples in an AOI in this bin.	Millimeters

AOI Visit metrics

An AOI visit corresponds to all the data between the start of the first fixation inside and AOI to the end of the last fixation in the same AOI. From the first fixation inside the AOI until the last fixation inside the AOI, all data is considered as part of the AOI visit (even saccades, blinks or invalid gaze data).

AOI visit metrics allow you to measure statistics based on visits inside an AOI (e.g. calculating revisiting rate of an AOI).

Metric name	Description	Format
Total duration of Visit	The total duration of the Visits inside an AOI during an interval.	Milliseconds
Average duration of Visit	The average duration of the Visits inside an AOI during an interval.	Milliseconds
Minimum duration of Visit	The duration of the shortest Visit inside an AOI during an interval.	Milliseconds
Maximum duration of Visit	The duration of the longest Visit inside an AOI during an interval.	Milliseconds

Metric name	Description	Format
Number of Visits	The number of Visits occurring in an AOI during an interval.	Count
Time to first Visit	Time in milliseconds to the first Visit inside an AOI during an interval.	Milliseconds
Duration of first Visit	The duration of the first Visit inside an AOI during an interval.	Milliseconds

AOI Visit metrics (Binning turned on)

Metric name	Description	Format
Visit hit	Any Visit hits inside an AOI, indicated by 1/0, for each bin.	Binary
Total duration of Visit	The total duration of the Visits inside an AOI during a bin.	Milliseconds
Number of Visit starts	The number of Visits inside an AOI that starts in the bin.	Count

AOI Glance metrics

All data is considered to be part of the AOI glance (even saccades, blinks or invalid gaze data) from the first saccade leading into the AOI until the last fixation inside the AOI.

Metric name	Description	Format
Total duration of Glances	The total duration of the Glances inside an AOI during an interval.	Milliseconds
Average duration of Glances	The average duration of the Glances inside an AOI during an interval.	Milliseconds
Minimum duration of Glances	The duration of the shortest Glance inside an AOI during an interval.	Milliseconds
Maximum duration of Glances	The duration of the longest Glance inside an AOI during an interval.	Milliseconds
Number of Glances	The number of Glances occurring in an AOI during an interval.	Count
Time to first Glance	Time in milliseconds to the first Glance inside this area of interest during an interval.	Milliseconds
Duration of first Glance	The duration of the first Glance inside this area of interest during an interval.	Milliseconds

AOI Glance metrics (Binning turned on)

Metric name	Description	Format
Glance hit	Any Glance hits inside an AOI, indicated by 1/0, for each bin.	Binary

Metric name	Description	Format
Total duration of Glances	The total duration of the Glances inside an AOI during a bin.	Milliseconds
Number of Glance starts	The number of Glances inside an AOI that starts in the bin.	Count

AOI Click metrics

All AOI Click metrics are based on the primary (left or right, depending on the Windows settings) button clicks only. By default, it is most commonly the left button.

A mouse click itself is a combination of two events - MouseEvent Down & MouseEvent Up where both have their timestamps and position. (Read [Screen project-specific Event groups](#): for more information.) There are two sets of metrics related to mouse clicks:

- **Clicks**: computed based on MouseEvent Down events only
- **Clicks & Releases**: takes both MouseEvent Down & MouseEvent Up events into account

Click & Release in Pro Lab counts only if both events spatially happened inside an AOI or a group of AOIs labeled by the same tag in case of AOI tag selection (even if MouseEvent Up and MouseEvent Down events happened in two different AOIs labeled by the same tag) and MouseEvent Down event occurred in a TOI interval.

All Click & Release temporal (time) metrics are computed (in milliseconds) based on the timestamp of the MouseEvent Down event.



The definition of *Clicks* changed in version 1.162. In the versions 1.152 and earlier, *Clicks* were what we now call *Clicks & Releases*.

Metric name	Description	Format
Number of mouse clicks	The number of times the mouse button is pressed in an AOI during an interval.	Count
Time to first mouse click	The time to when the mouse button is pressed inside an AOI during an interval.	Milliseconds
Time from first fixation to mouse click	The time from the first fixation to when the mouse button is pressed inside an AOI during an interval.	Milliseconds
Number of mouse clicks & releases	The number of times the mouse button is both pressed and released in the same AOI during an interval.	Count
Time to first mouse click & release	The time to the first mouse button is pressed the first time inside an AOI during an interval.	Milliseconds
Time from first fixation to mouse click & release	The time from first fixation to the first time the mouse button is pressed inside an AOI during an interval. This metric requires the mouse button is also released inside the same AOI.	Milliseconds

AOI Click metrics (Binning turned on)

Metric name	Description	Format
Click hit	Any mouse click (button is pressed) hits inside an AOI, indicated by 1/0, for each bin.	Binary
Number of mouse clicks	The number of times the mouse button is pressed in an AOI during a bin.	Count

GSR metrics

SCRs can be generated as a response to a specific event (e.g., visual stimulus or unexpected question) known as event-related SCR (ER-SCR). ER-SCRs are the most common measure used in research to relate changes in emotional arousal to a specific stimuli. A good stimulus design that allows enough time between stimuli is necessary to avoid uncertainties about which stimulus caused a specific ER-SCR.



The SCR is reported in the interval/bin when it starts to rise. It does not reflect the peak. For example, if the onset is in bin 1 and the peak is in bin 3, the value would be "1" for bin 1 and "0" for bins 2 and 3.

Metric name	Description	Format
Average GSR	The average galvanic skin response (GSR) signal, after filtering, for an interval.	Microsiemens
Number of SCR	The number of skin conductance responses (SCRs) for an interval.	Count
Amplitude of event related SCR	The amplitude of each event-related skin conductance response (ER-SCR), for an interval. ER-SCRs are calculated using filtered GSR data.	Microsiemens

GSR metrics (Binning turned on)

Metric name	Description	Format
Average GSR	The average galvanic skin response (GSR) signal, after filtering, for each bin.	Microsiemens
Number of SCR	The number of skin conductance responses (SCRs) in the bin.	Count

Fixation metrics (exclude partial fixations)

Fixation metrics (exclude partial fixations) let you measure statistics based on whole fixations within an interval (occurrence of a TOI) regardless of what the user specifically looked at. The fixations are defined based on the gaze filter you use and exclude fixations that don't fulfill the criteria for whole fixations.

Metric name	Description	Format
Total duration of whole fixations	The total duration of the fixations during an interval.	Milliseconds
Average duration of whole fixations	The average duration of the fixations during an interval.	Milliseconds
Number of whole fixations	The number of whole fixations occurring during an interval.	Count
Duration of first whole fixation	The duration of the first fixation during an interval.	Milliseconds
Average whole fixation pupil diameter	The average pupil diameter of all whole-fixation samples in this interval.	Millimeters

Fixation metrics (exclude partial fixations) (Binning turned on)

Metric name	Description	Format
Number of whole fixation starts	The number of whole fixations that starts during a bin.	Count
Average whole-fixation pupil diameter	The average pupil diameter of all whole-fixation samples in this bin.	Millimeters

Saccade metrics

Saccade metrics let you measure statistics based on saccades within an interval (occurrence of a TOI). You can get general indicators on the velocity, amplitude and direction of saccades.



If you have unrecognizable data, try adjusting the fixation filter settings. This is not a problem that can be fixed in the metrics.

Metric name	Description	Format
Number of saccades	The number of saccades occurring during an interval.	Count
Average peak velocity of saccades	The average peak velocity of all saccades in this interval.	Degrees/second
Minimum peak velocity of saccades	The peak velocity of the saccade with the lowest peak velocity in this interval.	Degrees/second
Maximum peak velocity of saccades	The peak velocity of the saccade with the highest peak velocity in this interval.	Degrees/second
Standard deviation of peak velocity of saccades	The standard deviation of all peak velocities of the saccades in this interval.	Degrees/second
Average amplitude of saccades	The average amplitude of all saccades in this interval.	Degrees
Minimum amplitude of saccades	The amplitude of the saccade with the lowest amplitude in this interval.	Degrees

Metric name	Description	Format
Maximum amplitude of saccades	The amplitude of the saccade with the highest amplitude in this interval.	Degrees
Total amplitude of saccades	The total amplitude of all saccades in this interval.	Degrees
Time to first saccade	The time to the first saccade during an interval.	Milliseconds
Direction of first saccade	The direction of the first saccade in the interval.	Degrees
Peak velocity of first saccade	The peak velocity of the first saccade in the interval.	Degrees/second
Average velocity of first saccade	The average velocity of the first saccade in the interval.	Degrees/second
Amplitude of first saccade	The amplitude of the first saccade in the interval.	Degrees

Saccade metrics (Binning turned on)

Metric name	Description	Format
Number of saccade starts	The number of saccades that starts during a bin.	Count

AOI saccade metrics

AOI saccades are saccades that start, end, or are within an AOI. AOI saccade metrics let you measure statistics based on saccades within an AOI. You can get general indicators on the velocity, amplitude and direction of these saccades.

Metric name	Description	Format
Number of saccades in AOI	The number of saccades occurring in an AOI during an interval.	Count
Time to entry saccade	The duration until the start of the first saccade that ends in an AOI during an interval.	Milliseconds
Time to exit saccade	The duration until the start of the first saccade that exits an AOI during an interval.	Milliseconds
Peak velocity of entry saccade	The peak velocity of the first saccade that ends in an AOI during an interval.	Degrees/second
Peak velocity of exit saccade	The peak velocity of the first saccade that exits an AOI during an interval.	Degrees/second

10.5.2 AOI-based TSV file

General

Metric name	Description	Format
Recording name	Recording name	
Participant	Participant	
Participant variables	Variable value, or values, of the participant. One column for each participant variable.	
Timeline name	Timeline name	
TOI	The name of the current Time of Interest.	
Interval	The interval number of the current TOI interval.	
Media	The name of the media presented to the participant.	
Stimulus variables	Stimulus variable value or values of the stimulus. One column for each Stimulus variable.	
AOI	The Area of Interest name of the current row.	
AOI Tags	The name or names of Tags connected to the AOI. One column for each Tag group and one for Ungrouped tags.	

General (Binning turned on)

Metric name	Description	Format
Recording name	Recording name	
Participant	Participant	
Participant variables	Variable value, or values, of the participant. One column for each participant variable.	
Timeline name	Timeline name	
TOI	The name of the current Time of Interest.	
Interval	The interval number of the current TOI interval.	
Bin	The index of the current bin in the interval.	
Bin duration	The duration of the current bin.	Milliseconds
Media	The name of the media presented to the participant.	
Stimulus variables	Stimulus variable value or values of the stimulus. One column for each Stimulus variable.	
AOI	The Area of Interest name of the current row.	
AOI Tags	The name or names of Tags connected to the AOI. One column for each Tag group and one for Ungrouped tags.	

Events

Metric name	Description	Format
Number of Events	The number of Events, including Custom Events and Logged live Events, for an interval.	Count

Metric name	Description	Format
Time to first Event	The time to the first Event, including Custom Events and Logged live Events, for an interval.	Milliseconds
Last key Press	The last registered key press in the interval.	

Events (Binning turned on)

Metric name	Description	Format
Event occurred	Any Event, including Custom Events and Logged live Events, occurs during a bin, indicated by 1/0, for each event type.	Binary
Number of Events	The number of Events, including Custom Events and Logged live Events, for a bin.	Count
Last key Press	The last registered key press in the bin.	

AOI fixation metrics

AOI fixations correspond to fixations that fall within an AOI. The fixations are defined based on the gaze filter you use (e.g. if you use the Raw gaze filter, every valid eye tracking sample is a fixation). AOI fixation metrics allow you to measure statistics based on the fixations within an AOI. They present as an interval (or an occurrence) of the TOI in separate rows in the exported spreadsheet.

Metric name	Description	Format
Total duration of fixations	The total duration of the fixations inside an AOI during an interval.	Milliseconds
Average duration of fixations	The average duration of the fixations inside an AOI during an interval.	Milliseconds
Minimum duration of fixations	The duration of the shortest fixation inside an AOI during an interval.	Milliseconds
Maximum duration of fixations	The duration of the longest fixation inside an AOI during an interval.	Milliseconds
Number of fixations	The number of fixations occurring in an AOI during an interval.	Count
Time to first fixation	The time to the first fixation inside an AOI during an interval.	Milliseconds
Duration of first fixation	The duration of the first fixation inside an AOI during an interval.	Milliseconds
Last AOI viewed	The last AOI fixated during an interval.	
AOI at interval end	The AOI fixated at the end of an interval.	
Average pupil diameter	The total duration of the Visits inside an AOI during an interval.	Milliseconds

AOI fixation metrics (Binning turned on)

Metric name	Description	Format
Fixation hit	Any fixation hits inside an AOI, indicated by 1/0, for each bin.	Binary
Total duration of fixations	The total duration of the fixations inside an AOI during a bin.	Milliseconds
Number of fixation starts	The number of fixations inside an AOI that starts in the bin.	Count
Average pupil diameter	The average pupil diameter of all fixation samples in an AOI in a bin.	Millimeters

AOI fixation metrics (exclude partial fixations)

These metrics exclude fixations that don't fulfill the criteria for whole fixations (see previous section).

Metric name	Description	Format
Total duration of whole fixations	The total duration of the fixations inside an AOI during an interval.	Milliseconds
Average duration of whole fixations	The average duration of the fixations inside an AOI during an interval.	Milliseconds
Minimum duration of whole fixations	The duration of the shortest fixation inside an AOI during an interval.	Milliseconds
Maximum duration of whole fixations	The duration of the longest fixation inside an AOI during an interval.	Milliseconds
Number of whole fixations	The number of fixations occurring in an AOI during an interval.	Count
Time to first whole fixation	The time to the first fixation inside an AOI during an interval.	Milliseconds
Duration of first whole fixation	The duration of the first fixation inside an AOI during an interval.	Milliseconds
Average whole fixation pupil diameter	The average pupil diameter of all whole fixation samples in an AOI in this interval.	Millimeters

AOI fixation metrics (exclude partial fixations) (Binning turned on)

Metric name	Description	Format
Whole fixation hit	Any fixation hits inside an AOI, indicated by 1/0, for each bin.	Binary
Total duration of whole fixations	The total duration of the fixations inside an AOI during a bin.	Milliseconds
Number of whole fixation starts	The number of fixations inside an AOI that starts in the bin.	Count

Metric name	Description	Format
Average whole-fixation pupil diameter	The average pupil diameter of all whole-fixation samples in an AOI in this bin.	Millimeters

AOI visit metrics

An AOI visit corresponds to all the data between the start of the first fixation inside and AOI to the end of the last fixation in the same AOI. From the first fixation inside the AOI until the last fixation inside the AOI, all data is considered as part of the AOI visit (even saccades, blinks or invalid gaze data).

AOI visit metrics allow you to measure statistics based on visits inside an AOI (e.g. calculating revisiting rate of an AOI).

Metric name	Description	Format
Average whole fixation pupil diameter	The total duration of the Visits inside an AOI during an interval.	Milliseconds
Total duration of Visit	The total duration of the Visits inside an AOI during an interval.	Milliseconds
Average duration of Visit	The average duration of the Visits inside an AOI during an interval.	Milliseconds
Minimum duration of Visit	The duration of the shortest Visit inside an AOI during an interval.	Milliseconds
Maximum duration of Visit	The duration of the longest Visit inside an AOI during an interval.	Milliseconds
Number of Visits	The number of Visits occurring in an AOI during an interval.	Count
Time to first Visit	Time in milliseconds to the first Visit inside an AOI during an interval.	Milliseconds
Duration of first Visit	The duration of the first Visit inside an AOI during an interval.	Milliseconds

AOI visit metrics (Binning turned on)

Metric name	Description	Format
Visit hit	Any Visit hits inside an AOI, indicated by 1/0, for each bin.	Binary
Total duration of Visit	The total duration of the Visits inside an AOI during a bin.	Milliseconds
Number of Visit starts	The number of Visits inside an AOI that starts in the bin.	Count

AOI Glance metrics

All data is considered to be part of the AOI glance (even saccades, blinks or invalid gaze data) from the first saccade leading into the AOI until the last fixation inside the AOI.

Metric name	Description	Format
Total duration of Glances	The total duration of the Glances inside an AOI during an interval.	Milliseconds
Average duration of Glances	The average duration of the Glances inside an AOI during an interval.	Milliseconds
Minimum duration of Glances	The duration of the shortest Glance inside an AOI during an interval.	Milliseconds
Maximum duration of Glances	The duration of the longest Glance inside an AOI during an interval.	Milliseconds
Number of Glances	The number of Glances occurring in an AOI during an interval.	Count
Time to first Glance	Time in milliseconds to the first Glance inside an AOI during an interval.	Milliseconds
Duration of first Glance	The duration of the first Glance inside an AOI during an interval.	Milliseconds

AOI Glance metrics (Binning turned on)

Metric name	Description	Format
Glance hit	Any Glance hits inside an AOI, indicated by 1/0, for each bin.	Binary
Total duration of Glances	The total duration of the Glances inside an AOI during a bin.	Milliseconds
Number of Glance starts	The number of Glances inside an AOI that starts in the bin.	Count

AOI Click metrics

All AOI Click metrics are based on the primary (left or right, depending on the Windows settings) button clicks only. By default, it is most commonly the left button.

A mouse click itself is a combination of two events - MouseEvent Down & MouseEvent Up where both have their timestamps and position. (Read [Screen project-specific Event groups](#): for more information.) There are two sets of metrics related to mouse clicks:

- **Clicks**: computed based on MouseEvent Down events only
- **Clicks & Releases**: takes both MouseEvent Down & MouseEvent Up events into account

Click & Release in Pro Lab counts only if both events spatially happened inside an AOI or a group of AOIs labeled by the same tag in case of AOI tag selection (even if MouseEvent Up and MouseEvent Down events happened in two different AOIs labeled by the same tag) and MouseEvent Down event occurred in a TOI interval.

All Click & Release temporal (time) metrics are computed (in milliseconds) based on the timestamp of the MouseEvent Down event.



The definition of "Clicks" changed in version 1.162. "Clicks & Releases" was called simply "Clicks" in earlier versions.

Metric name	Description	Format
Number of mouse clicks	The number of times the mouse button is pressed in an AOI during an interval.	Count
Time to first mouse click	The time to when the mouse button is pressed inside an AOI during an interval.	Milliseconds
Time from first fixation to mouse click	The time from the first fixation to the first time the mouse button is pressed inside an AOI during an interval.	Milliseconds
Number of mouse clicks & releases	The number of times the mouse button was both pressed and released in the same AOI during an interval.	Count
Time to first mouse click & release	The time to when the mouse button is pressed inside an AOI during an interval. This metric requires the mouse button is also released inside the same AOI.	Milliseconds
Time from first fixation to mouse click & release	The time from the first fixation to the first time the mouse button is pressed inside an AOI during an interval. This metric requires the mouse button is released inside the same AOI.	Milliseconds

AOI Click metrics (Binning turned on)

Metric name	Description	Format
Click hit	Any mouse click (button is pressed) hits inside an AOI, indicated by 1/0, for each bin.	Binary
Number of mouse clicks	The number of times the mouse button is pressed in an AOI during a bin.	Count

AOI saccade metrics

AOI saccades are saccades that start, end, or are within an AOI. AOI saccade metrics let you measure statistics based on saccades within an AOI. You can get general indicators on the velocity, amplitude and direction of these saccades.

Metric name	Description	Format
Number of saccades in AOI	The number of saccades occurring in an AOI during an interval.	Count
Time to entry saccade	The duration until the start of the first saccade that ends in an AOI during an interval.	Milliseconds
Time to exit saccade	The duration until the start of the first saccade that exits an AOI during an interval.	Milliseconds
Peak velocity of entry saccade	The peak velocity of the first saccade that ends in an AOI during an interval.	Degrees/second
Peak velocity of exit saccade	The peak velocity of the first saccade that exits an AOI during an interval.	Degrees/second

AOI reading metrics

AOI reading metrics are only available for text stimuli and are AOI-based metrics on Automatic AOIs, generated for a text stimulus.

These metrics follow the writing system order. AOI1 comes before AOI2 which comes before AOI3, etc. (This determines the occurrence of regressions, progression, and skipping). Metrics for AOIn are computed using other AOIs as well (AOIn-2, AOIn-1, AOIn+1). The AOI selection determines what gets displayed in the data but if the calculations use the data from adjacent AOIs (even if they are not selected), the calculations will still be correct.

AOIs have 3 different levels and all reading metrics are level dependent. Reading metrics on different level AOIs are computed independently from each other. Several Automatic reading-related AOIs can exist at the same time (with their own AOI order):

- Character AOIs
- Word AOIs
- Sentence AOIs

Reading metrics are based only on fixations.

- A regression/progression is determined if there is a fixation afterwards on an area of interest with a respective lower/higher index of the same level (Word, Character, Sentence) regardless of the direction of the saccade itself.
- If the trial/interval ends in the middle of a saccade that would have resulted in a fixation in a regressive position, it will not count as a regression.

Non-AOI data is ignored.

- Fixations landing outside of text AOIs are disregarded, and will not contribute to any metrics calculation, nor terminate any metric calculation.
- A fixation in an AOI, followed by a fixation outside of it, and then followed by a fixation inside the AOI again, will be equivalent to having two fixations in the AOI directly followed by each other. Both fixations in the AOI will count as part of the same pass.

Metric name	Description	Format
Character index	Index of character-level AOI inside its word-level AOI.	Position
Word index	Index of word-level AOI inside its sentence-level AOI.	Position
Sentence index	Index of sentence-level AOI inside this text stimulus.	Position
AOI string	Text string contained in an AOI.	
Text unit type	Type of text unit: character, word, sentence, or custom.	
Number of units	Number of units.	Count

First-pass first fixation duration	The duration of the first fixation during first-pass inside an AOI during an interval.	Milliseconds
First-pass duration	The total duration of the fixations during first-pass inside an AOI during an interval.	Milliseconds
Selective regression-path duration*	The total duration of the fixations from first fixation in this area of interest until a fixation occurs in an area of interest progressive to this one, during an interval. *Previously called "Go-past duration"	Milliseconds
First pass regression	Indicates whether the reader exits the AOI with a regression (1) or reads on progressively (0) during an interval.	Boolean
Total duration of fixations	The total duration of the fixations inside an AOI during an interval.	Milliseconds
Regression-path duration	The total duration of the fixations from first fixation in this area of interest until a fixation occurs in an AOI progressive to this one, including fixations in regressive AOIs, during an interval.	Milliseconds
Re-reading duration	Regression path duration excluding first pass fixations during an interval.	Milliseconds

10.5.3 Event-based TSV file

Use this format when you want to analyze individual events during a trial or interval. Each fixation will generate a row in the report and the selected metrics will be shown as columns. Just like the [Interval-based TSV file](#), it is formatted to be easy to import and analyze in statistical analysis software.

General

Metric name	Description	Format
Recording name	Recording name	
Participant	Participant	
Participant variables	Variable value, or values, of the participant. One column for each participant variable.	
Timeline name	Timeline name	
TOI	The name of the current Time of Interest.	
Interval	The interval number of the current TOI interval.	
Bin duration	The duration of the current bin.	Milliseconds
Media	The name of the media presented to the participant.	
Stimulus variables	Stimulus variable value or values of the stimulus. One column for each Stimulus variable.	
AOI Tags*	The name or names of Tags connected to the AOI. One column for each Tag group and one for Ungrouped tags.	

*AOI Tags is not available when binning is turned on

Event properties

These properties are shared for all events types that are covered by the event-based metrics. Currently, only fixations are supported.

Metric name	Description	Format
Event type	The type of event of the current row.	Fixation; Saccade
Validity	The validity of the event of the row, either whole or partial.	Partial; Whole
EventIndex	Represents the order of the events in the current TOI interval. The index is an auto-increment number starting with 1 for each event type.	Position
Start	The start time counted from current TOI interval start.	Milliseconds
Stop	The stop time counted from current TOI interval start.	Milliseconds
Start bin	The bin the event starts in.	Count
Stop bin	The bin the event stops in.	Count
Duration	The duration of the event.	Milliseconds

Event properties (Binning turned on)

Metric name	Description	Format
Event type	The type of event of the current row.	Fixation; Saccade
Validity	The validity of the event of the row, either whole or partial.	Partial; Whole
EventIndex	Represents the order of the events in the current TOI interval. The index is an auto-increment number starting with 1 for each event type.	Position
Start	The start time counted from current TOI interval start.	Milliseconds
Stop	The stop time counted from current TOI interval start.	Milliseconds
Start bin	The bin the event starts in.	Count
Stop bin	The bin the event ends in.	Count
Duration	The duration of the event.	Milliseconds

Fixation properties

Information specific to each fixation as well as to the general event properties.

Metric name	Description	Format
AOI	The name of the AOI (s) which the current fixation hits.	
AOI proportion	The proportion of the fixation that occurs within the AOI.	

Metric name	Description	Format
Fixation point	The normalized horizontal and vertical coordinate of the fixation point.	Normalized coordinates (DACS)
Average pupil size	The average size of the pupil of the fixation. Note: If the fixation is cut by TOI interval borders, only the gaze samples within the TOI are included.	Millimeters

Fixation properties (Binning turned on)

Metric name	Description	Format
AOI	The name of the AOI (s) which the current fixation hits.	
AOI proportion	The proportion of the fixation that occurs within the AOI.	
Fixation point	The normalized horizontal and vertical coordinate of the fixation point.	Normalized coordinates (DACS)
Average pupil size	The average size of the pupil of the fixation. Note: If the fixation is cut by TOI interval borders, only the gaze samples within the TOI are included.	Millimeters

10.5.4 Excel Report

General

Metric name	Description	Format
Recording name	Recording name	
Participant	Participant	
Participant variables	Variable value or values of the participant. One column for each participant variable	
TOI	The name of the current Time of Interest.	
Interval	The interval number of the current TOI interval.	

Interval metrics

Metric name	Description	Format
Duration of interval	The duration of all time Intervals for each Time of Interest, with averages, medians, sums, counts, variances and standard deviations (n-1).	Seconds
Start of interval	The start time of all time Intervals for each Time of Interest, with averages, medians, counts, variances, and standard deviations (n-1).	Seconds

Event metrics

Events can also be used in measures. Event metrics allow you to measure behavior and calculate statistics based on your event coding scheme.

Metric name	Description	Format
Number of Events	The number of Events, including Custom Events and Logged live Events, for each Time of Interest, with averages, medians, counts, variances, and standard deviations (n-1). Descriptive statistics only include recordings where Events occur.	Count
Number of Events (include zeroes)	The number of Events, including Custom Events and Logged live Events, for each Time of Interest, with averages, medians, counts, variances, and standard deviations (n-1). Descriptive statistics also include recordings where no Events occur.	Count
Time to first Event	The time to first Event, including Custom Events and Logged live Events, for each Time of Interest, with averages, medians, counts, variances, and standard deviations (n-1).	Seconds

AOI fixation metrics

AOI fixations correspond to fixations that fall within an AOI. The fixations are defined based on the gaze filter you use (e.g. if you use the Raw gaze filter, every valid eye tracking sample is a fixation). AOI fixation metrics allow you to measure statistics based on the fixations within an AOI. They present as an interval (or an occurrence) of the TOI in separate rows in the exported spreadsheet.

Metric name	Description	Format
Total duration of fixation in AOI	The total time each participant has fixated each AOI on all Media, with averages, medians, sums, variances, and standard deviations (n-1); the share of total time spent in each AOI out of all AOIs; and the percentage of Participants that fixated within each AOI at least once. Descriptive statistics only based on Recordings with fixations within the AOIs.	Seconds
Total duration of fixation in AOI (include zeroes)	The total time each participant has fixated each AOI on all Media, with averages, medians, sums, variances, and standard deviations (n-1); the share of total time spent in each AOI out of all AOIs; and the percentage of Participants that fixated within each AOI at least once. Descriptive statistics also include Recordings with 0 fixations within the AOIs.	Seconds

Metric name	Description	Format
Average duration of fixation in AOI	The average duration of the fixations within each AOI on all Media, with averages, medians, variances, and standard deviations (n-1); the total Time of Interest and Recording durations.	Seconds
Number of fixations in AOI	The number of fixations within each AOI on all Media, with averages, medians, sums, variances, and standard deviations (n-1); the percentage of Participants that visited each AOI at least once; total number of fixations within the Time of Interest; and the total Time of Interest and Recording Durations. Descriptive statistics only based on Recordings with fixations within the AOIs.	Count
Number of fixations in AOI (include zeroes)	The number of fixations within each AOI on all Media, with averages, medians, sums, variances, and standard deviations (n-1); the percentage of Participants that visited each AOI at least once; total number of fixations within the Time of Interest; and the total Time of Interest and Recording Durations. Descriptive statistics also include Recordings with 0 fixations within the AOIs.	Count
Time to first fixation in AOI	The time to first fixation for each AOI on all Media, with averages, medians, counts, variances, standard deviations (n-1) and Recording durations.	Seconds
Duration of first fixation in AOI	The duration of the first fixation for each AOI on all Media, with averages, medians, counts, variances, standard deviations (n-1) and Recording durations.	Seconds

AOI Visit metrics

An AOI visit corresponds to all the data between the start of the first fixation inside and AOI to the end of the last fixation in the same AOI. From the first fixation inside the AOI until the last fixation inside the AOI, all data is considered as part of the AOI visit (even saccades, blinks or invalid gaze data).

AOI visit metrics allow you to measure statistics based on visits inside an AOI (e.g. calculating revisiting rate of an AOI).

Metric name	Description	Format
Total duration of Visit	The total time each participant has visited each AOI on all Media, with averages, medians, sums, variances, and standard deviations (n-1); the share of total time spent in each AOI out of all AOIs; and the percentage of Participants that visited each AOI at least once. Descriptive statistics are only based on Recordings with fixations within the AOIs.	Seconds

Metric name	Description	Format
Total duration of Visit (include zeroes)	The total time each participant has visited each AOI on all Media, with averages, medians, sums, variances, and standard deviations (n-1); the share of total time spent in each AOI out of all AOIs; and the percentage of Participants that visited each AOI at least once. Descriptive statistics also include Recordings with 0 fixations within the AOIs.	Seconds
Average duration of Visit	The average duration each participant has visited each AOI on all Media, with averages, medians, sums, variances, and standard deviations (n-1).	Seconds
Number of Visits	The number of Visits within each AOI on all Media, with averages, medians, variances, and standard deviations (n-1); and the percentage of Participants that fixated within each AOI at least once. Descriptive statistics only based on Recordings with fixations within the AOIs.	Count
Number of Visits (include zeroes)	The number of Visits within each AOI on all Media, with averages, medians, variances, and standard deviations (n-1); and the percentage of Participants that fixated within each AOI at least once. Descriptive statistics also include Recordings with 0 fixations within the AOIs.	Count

AOI Click metrics

One click is defined as the combination of when the participant presses the primary (left or right) button of the mouse, and when he or she releases it again.

Metric name	Description	Unit
Number of clicks & releases in AOI	The number of times the mouse button is both pressed and released within each AOI on all Media, with averages, medians, variances, and standard deviations (n-1); and the percentage of Participants that clicked within each AOI at least once. Descriptive statistics only based on Recordings with fixations within the AOIs.	Count
Number of clicks & releases in AOI (include zeroes)	The number of times the mouse button is both pressed and released within each AOI on all Media, with averages, medians, variances, and standard deviations (n-1); and the percentage of Participants that clicked within each AOI at least once. Descriptive statistics also include Recordings with 0 clicks within the AOIs.	Count
Time to first click & release in AOI	The time to first mouse button is pressed for each AOI on all Media, with averages, medians, counts, variances, standard deviations (n-1) and Recording durations. This metric requires the mouse button is also released inside the same AOI.	Seconds

Metric name	Description	Unit
Time from first fixation to mouse click in AOI	The time from first fixation to next time the mouse button is pressed for each AOI on all Media, with averages, medians, counts, variances, standard deviations (n-1), Recording durations and the percentage of Participants that fixated and then clicked within each AOI at least once. This metric requires the mouse button is also released inside the same AOI.	Seconds

GSR metrics

SCRs can be generated as a response to a specific event (e.g., visual stimulus or unexpected question) known as event-related SCR (ER-SCR). ER-SCRs are the most common measure used in research to relate changes in emotional arousal to a specific stimuli. A good stimulus design that allows enough time between stimuli is necessary to avoid uncertainties about which stimulus caused a specific ER-SCR.

Metric name	Description	Unit
GSR Average	The average galvanic skin response (GSR) signal, after filtering, for each Time of Interest, with averages, medians, and counts for each participant.	Microsiemens
ER SCR Amplitude	The amplitude of each event related skin conductance response (ER-SCR), for each Interval in Time of Interest, with mean amplitudes, mean magnitudes, response frequencies, and counts for each participant. Time of Interest intervals that does not have an ER-SCR are shown with the symbol "---". ER-SCRs are calculated using filtered GSR data.	Microsiemens
SCR Count	The number of skin conductance responses (SCRs), for each Interval in Time of Interest, with averages, medians, counts, variances, and standard deviations (n-1).	Number

11 Data export

11.1 Export eye tracking, GSR and recording data in a *.tsv file



Gaze filter selection is applicable for eye tracking data only. These settings do not affect the GSR data export.

11.1.1 How to export data

1. Click the **Data Export** button, either in the Project Overview section or select it from the **Analyze** drop-down menu in the top navigation.
2. Select which Gaze Filter you want to apply to the data you will export. For more information about Gaze Filters, read [Gaze Filter functions and effects](#).
 - Select the data you want to include in the *Data Selection* section of the tool on the right of the interface.
 - Under *Recordings*, you can select to export data from All Recordings or only the recordings of your choice.
Click on the triangle next to the checkbox to see all available recordings.
 - Under *Times of Interest*, you can select to export data metrics from the complete recordings (selected above) or data only from selected TOIs. To simplify the selection, you are presented with two groups of TOI data: Media gaze data, where automatically-created TOIs, based on when Snapshots or images were active, are presented, and Recording gaze data, where your custom TOIs are listed. In the *Media gaze* data list, each image or Snapshot associated with the project is listed and, for each of them, you can also select TOIs which include data from when they were active. Click the triangle next to the checkbox by each group to see and select the individual TOIs.
 - Under *Areas of Interest*, you can choose to export data relating to all, or just some, of the AOIs created in the project.
3. Click the **Export** button at the top right of the interface. A file browser appears.
4. Locate the folder you want to save the file in. In the File name field, enter a file name for your export.
5. Click **Save**.

11.2 Data export formats

Exports from Data Export are saved in a tab-separated values file (.tsv) that follows the Unicode standard. The *.tsv output file contains columns. Each column contains data of a type given by the data type name found in the top row for the corresponding column. All data types are described in the tables in [Data export information](#). Images and Snapshots have a set of their own columns with information about the image or Snapshot itself and the gaze data mapped to it. Thus, each added image or Snapshot produces additional columns in the output file. The same is true for Areas of Interest, where each AOI will get its own column in the Data Export.

All rows in a Data Export file have a Recording Timestamp value (except the first row, which contains the column data type name). You can choose whether the timestamp shows milliseconds or microseconds. The timestamp starts at zero at the beginning of each recording.

Since all recorded eye gaze data samples are recorded in a sequence, all eye gaze data points in a recording will have different timestamps. However, some Events may have the same timestamp as eye gaze data points and others may have timestamps between two eye gaze data point timestamps.

Gaze data points and Events have their own rows in the export file so the relationship between the number of rows and time is not linear. Instead, timestamps must be used when plotting/charting eye gaze data from a Data Export file.

In a Data Export file, you will also have a Computer Timestamp value. The tables in [Data export information](#) list the type of information and data types available for export from Pro Lab. Each type has its own column in the Data Export output file.

Read more about timestamps in [Computer timestamp \(screen-based\)](#) and [Computer timestamp \(wearable\)](#).

11.2.1 Computer timestamp (screen-based)

For screen-based recordings (including scene camera and external presenter), the Computer timestamp column contains the value of the win32 clock "QueryPerformanceCounter" (QPC) in microseconds. This means that if other software running on the same computer collects data, and this data is timestamped with QPC, the data can be synced with the data recorded in Pro Lab.

This is the same clock provided by [Tobii Pro SDK](#).

11.2.2 Computer timestamp (wearable)

For Pro Glasses 2 recordings, the Computer timestamp value comes from the internal clock of the recording unit and not the computer running Pro Glasses 2 controller application. This clock starts when the recording unit is booted. It is not possible to use this clock for synchronization of other data sources.

For Pro Glasses 3 recordings, the Computer timestamp value comes from an internal clock that is initialized when the recording is started, so the value will be identical to the Recording Timestamp.

11.3 Data export information

The following tables list the data types and information available for export from Pro Lab. Each type has its own column in the Data Export output file.

11.3.1 General

Data name	Description	Format	Screen project	Glasses project	Scene Camera project	External Presenter project
Project name	Project name	Text	•	•	•	•
Export date	Date when the Data Export is done.	YYYY-MM-DD	•	•	•	•

Data name	Description	Format	Screen project	Glasses project	Scene Camera project	External Presenter project
Participant name	Participant name	Text	•	•	•	•
Participant variables	Variable value or values of the participant.	Text	•	•	•	•
Recording name	Recording name	Text	•	•	•	•
Recording date	Date when the Recording was performed in this time zone.	YYYY-MM- DD	•	•	•	•
Recording date UTC	Date when the Recording was performed in UTC.	YYYY-MM- DD	•	•	•	•
Recording start time	Start time of the Recording in this time zone.	HH:MM: SS:FFF	•	•	•	•
Recording start time UTC	Start time of the Recording in UTC format	HH:MM: SS:FFF	•	•	•	•
Recording duration	Total duration of the recording	Milliseconds	•	•	•	•
Timeline name	Name of the Timeline used during the Recording.	Text	•	•		
Recording Fixation filter name	The name of the Fixation Filter applied to the Recording eye tracking data in the export.	Text	•	•	•	•
Recording software version	The version of the software used to make the Recording.	Text	•		•	•
Recording resolution	Screen resolution used during the Recording.	Pixels	•		•	•
Recording monitor latency	The monitor latency setting for the Recording. Stimulus start and end Event timestamps have been offset by this number to account for the monitor latency.	Milliseconds	•			•
Calibration results	Average accuracy and precision of calibration.	Millimeters, degrees and pixels.	•		•	•
Validation results	Average accuracy and precision of validation.	Millimeters, degrees and pixels.	•			

Data name	Description	Format	Screen project	Glasses project	Scene Camera project	External Presenter project
Eye tracker timestamp	The Recording timestamp in the eye tracker clock.	Microseconds	•		•	•
Event	Name of the Event.	Text	•	•	•	•
Event value	The event value.	Text	•	•	•	•

11.3.2 Eye tracking data

Data name	Description	Format/Units	Screen project	Glasses project	Scene Camera project	External Presenter project
Gaze point 2D	Raw gaze coordinates for each eye individually.	Pixels (DACS)	•			•
Gaze point 2D	Raw gaze coordinates for both eyes combined.	Pixels (MCS)		•	•	
Gaze point 3D	The vergence point of left and right gaze vectors.	Millimeters (HUCS)		•		
Gaze direction	The unit vector for the direction of the gaze, for each eye individually.	Normalized coordinates (DACS)	•		•	•
Gaze direction	The unit vector for the direction of the gaze, for each eye individually.	Normalized coordinates (HUCS)		•		
Pupil position	The 3D coordinates of the pupil position for each eye individually.	Millimeters (HUCS)		•		
Pupil diameter	Estimated size of the pupils.	Millimeters	•	•	•	•
Validity of eye data	Indicates if the eyes have been correctly identified.	Valid/invalid	•	•	•	•
Eye position (DACSmm)	3D position of the eyes.	Millimeters (DACS)	•		•	•
Gaze point 2D (DACSmm)	Raw gaze coordinates for each eye individually.	Millimeters (DACS)	•			•
Gaze point (MCSnorm)	Raw gaze coordinates for each eye individual on the Media.	Normalize coordinates (MCS)	•		•	•
Assisted mapping gaze point	Assisted mapping gaze point coordinates.	Pixels (MCS)	•	•	•	•

Data name	Description	Format/Units	Screen project	Glasses project	Scene Camera project	External Presenter project
Manually mapped gaze point	Manually mapped gaze point coordinates.	Pixels (MCS)	•	•	•	•
Mapped gaze point	The combination of the manually and assisted mapped gaze point coordinates. Manual mapping overrides assisted.	Pixels (MCS)	•	•	•	•
Assisted mapping gaze point score	Similarity score of assisted mapping gaze points.	Normalized	•	•	•	•

11.3.3 Media

Data name	Description	Format/Units	Screen project	Glasses project	Scene Camera project	External presenter project
Presented Stimulus name	The name of the Stimulus being presented to the Participant.	Text	•			•
Presented Media name	The name of the Media presented to the Participant.	Text	•			•
Recording Media name	The name of the Recording Media.	Text		•		
Presented Media dimensions	The dimensions of the Media as presented on the screen to the Participant, including any scaling set in the Stimulus properties.	Pixels	•			•
Recording Media dimensions	The dimensions of the Recording Media.	Pixels		•		
Presented Media position	The position of the Media on the screen. The value represents the positions of the top left corner of the Media in relation to the top left corner of the screen.	Pixels (DACS)	•			•
Original Media dimensions	The original size of the Media presented to the Participant.	Pixels	•			•

Data name	Description	Format/Units	Screen project	Glasses project	Scene Camera project	External presenter project
Media dimensions	The original size of the Snapshot.	Pixels	•	•	•	•

11.3.4 Gaze events

Data name	Description	Format/Units	Screen project	Glasses project	Scene Camera project	External Presenter project
Mapped eye movement type	Type of eye movement event classified by the selected Fixation filter for mapped gaze data.	Fixation Saccade Unclassified EyesNotFound	•	•	•	•
Mapped eye movement type index	Represents the order in which an eye movement was recorded for mapped gaze data. The index is an auto-increment number starting with 1 for each eye movement type.	Number	•	•	•	•
Mapped fixation point	Mapped fixation point. This column is affected by the settings of the Fixation Filter.	Pixels (MCS)	•	•	•	•
Eye movement type	Type of eye movement event classified by the fixation filter settings applied during the gaze data export.	Fixation Saccade Unclassified EyesNotFound	•	•	•	•
Gaze event duration	The duration of the currently active eye movement.	Milliseconds	•	•	•	•
Eye movement type index	Represents the order in which an eye movement was recorded. The index is an auto-increment number starting with 1 for each eye movement type.	Number	•	•	•	•

Data name	Description	Format/Units	Screen project	Glasses project	Scene Camera project	External Presenter project
Fixation point	Coordinates of the fixation point. This column is affected by the settings of the Fixation Filter.	Pixels (DACS)	•	•	•	•
Fixation point (MCSnorm)	Coordinates of the fixation point on the Media.	Normalized coordinates (MCS)	•			•
AOI hit	Reports whether the AOI is active and whether the fixation is located inside of the AOI: -1 = AOI not active; 0 = AOI active, the fixation is not located in the AOI; 1 = AOI active and the fixation is located inside of the AOI; empty cell indicates that the media of the AOI was not visible.	Number	•	•	•	•

11.3.5 Web data

Data name	Description	Format/Units	Screen Project	Glasses project	Scene Camera project	External presenter
Browser client area position	The position of a web browser's client area on the screen. The value represents the position of the top left corner of the client area in relation to the top left corner of the screen.	Pixels (DACS)	•			
Viewport position	The position of the visible area of a web page. The value represents the position of the top left corner of the visible area of a web page in relation to the full web page size.	Pixels	•			
Viewport dimensions	The dimensions of the visible area of a web page.	Pixels	•			

Data name	Description	Format/Units	Screen Project	Glasses project	Scene Camera project	External presenter
Full page size	The full size of the web page. Limited by 5000 px horizontally and 15000 px vertically.	Pixels	•			

11.3.6 Other sensor data

Data name	Description	Format/Units	Screen project	Glasses project	Scene Camera project	External presenter project
Mouse position	The position of the mouse.	Pixels (DACs)	•			
Gyro	Rotation along the X, Y and Z axes.	Degrees/second (HUCS)		•		
Accelerometer	Acceleration along the X, Y and Z axes.	Meters/second ² (HUCS)		•		
Magnetometer	Magnetic field along the X, Y and Z axes.	Microteslas (HUCS)		•		
Galvanic skin response (GSR)	The raw galvanic skin response signal of the Participant.	Microsiemens	•		•	•

11.4 Plof file format

The Pro Lab Output Format (plof) exports all data in a machine-readable format. The main goal for this format is to enable third-party software and researchers an easy and robust import of Tobii Pro eye tracking data. All data stored in Pro Lab can be exported, including raw eye tracking data, eye movement data, manual event coding data, stimulus information data, areas of interest data, raw GSR data and GSR events.

For more detailed information about plof, request the “Tobii Pro Lab Output Format Reference guide” from [Tobii Pro Sales](#).

Appendix A Keyboard shortcuts

A1 Dashboard

Action	Keyboard shortcut
Open selected recording	Enter
Rename selected recording	F2
Delete selected recording	Delete

A2 Design

Action	Keyboard shortcut
Select parent element	Backspace
Expand selected element	Enter
Remove selected element	Delete
Copy	Ctrl + C
Paste	Ctrl + V
Remove cut elements from clipboard	Esc

A3 Recording

Action	Keyboard shortcut
Stop recording	Esc
Calibrate current point (in calibration moderator view / participant screen)	Space

A4 Web recording

Action	Keyboard shortcut
Reload	Ctrl + R or Delete
Reload (in browser)	Refresh button (in browser)
Stop	Esc
Back (in browser)	Back button (in browser)
Forward (in browser)	Forward button (in browser)

A5 Replay

Action	Keyboard shortcut
Play / Pause	Space
Confirm automapped point	C
Clear current selection	Esc
Delete mapped fixation	Delete

A6 Areas of interest (AOI)

Action	Keyboard shortcut
Rename AOI	F2
Go to next frame	F*
Go to previous frame	B
Delete AOI	Delete
Abort tool	Esc
Undo	Ctrl + Z
Redo	Ctrl + Y
Go to next key-frame	Ctrl + F
Go to previous key-frame	Ctrl + B
Select all	Ctrl + A
Copy AOIs	Ctrl + C
Cut AOIs	Ctrl + X
Cut AOIs	Shift + Delete
Paste AOIs in center of viewport	Shift + Insert
Paste AOIs in center of viewport	Ctrl + V
Select AOI width edit box	W**
Select AOI height edit box	H
Zoom tool	Z
Pan tool	Space
Select AOI tool	S
Select AOI vertex tool	V
Draw polygon tool	P
Draw ellipse tool	E
Draw rectangle tool	R
Constraint proportional mode (while drawing)	Shift

* Moving to the next and previous frames in AOI tool only works when hovering mouse over the time bar

** Edit box selection in AOI tool only works when no other edit box is currently being edited

Appendix B Coordinate systems and data mapping



Pro Glasses 2 and Pro Glasses 3 use the intersection of the gaze vectors from the two eyes to calculate the distance at which the participant is looking. If the eye tracker loses tracking for one of the eyes, it will continue using the last known distance until it gets gaze data from both eyes again.

B1 Display Area Coordinate System (DACS)




Figure 19. Display Area Coordinate System (DACS)

When working in a Screen project, the data is mapped onto what is seen on the screen.

DACS is a 3D coordinate system with its origin in the top left corner of the screen (stimuli area) with Y pointing downwards. (A coordinate that is lower on the screen will have a higher Y-value). This may feel unintuitive for a 3D coordinate, but it is common in screen-based coordinate systems. It is oriented so that it is aligned with the plane of the screen, so any coordinate on the screen will have its z-coordinate equal to zero. Gaze points, eye positions, and gaze directions are all available in DACS. The 3D coordinates for gaze point and eye position are expressed in millimeters. Screen positions in pixels or normalized coordinates can easily be translated to DACS 3D coordinates. The User Coordinate System in Tobii eye trackers uses this during the screen setup of the eye tracker. The screen setup defines the three corners of the screen in UCS

coordinates (ULLLUR, Upper Left, Lower Left and Upper Right). This data is then used by Pro Lab to transform UCS coordinates into DACS coordinates.




Figure 20. In a Screen project, the data is mapped onto what is seen on the screen.

When media is presented on the screen, it typically extends over a portion of the screen. Pro Lab keeps track of exactly where on the screen the media is displayed and transforms coordinates (gaze points, mouse cursor positions) from DACS (the entire screen) to media coordinates (MCS). MCS coordinates are either in pixels (media pixels, not necessarily screen pixels) or normalized in order to be comparable between recordings. If the same media was displayed on two different screens, with “fit to screen” scaling, both pixel coordinates and normalized coordinates can still be compared and the same coordinate will always refer to the same position on the media. Origin in MCS is the top left corner of the media. Y is pointing downwards.

B2 Scene Camera Projects and the Media Coordinate System (MCS)

When eye tracking is done in a real-world environment (scene camera projects), what is seen by the study participant is filmed using a scene camera. When viewing the recording, the gaze data is mapped from the coordinate system of the eye tracker to the 2D coordinate system of the scene camera video.




Figure 21. Scene Camera Projects and the Media Coordinate System (MCS)

In scene camera recordings there is not necessarily a computer screen, but there is a surface that the eye tracker is configured on which to map gaze. This could be a computer screen, but it could also be a whiteboard on the wall, a flat surface on a table, or a virtual surface in the scene. If the user is focusing on objects that are not in the plane of the configured surface, the eye tracking data will be incorrect. In order to map gaze from this surface to the scene camera video, it is important that the edges of the surface configured in the eye tracker are marked in the scene camera video image. This allows Pro Lab to calculate a transformation (homography) from the DACS coordinate system to the scene camera video (MCS). Since the scene camera video does not have any physical dimensions and the homography only allows translation in 2D, coordinates in MCS can only be provided in relative coordinates or pixels.

B3 Head Unit Coordinate System (HUCS)




Figure 22. Head Unit Coordinate System (HUCS)

Just like for screen-based trackers, the glasses measures everything in a 3D coordinate system. Eye position and gaze vectors are calculated from the eye images using a 3D eye model that gives positions and angles in a coordinate system with its origin in the center of the scene camera.

The gaze point is calculated as the vergence point between the two gaze vectors. This means that there can be only one 3D gazepoint even though both eyes are tracked. This is because the eye tracker has no way of knowing how far away (or close) you are looking. A screen-based eye tracker uses the intersection of the screen plane and the gaze direction, but in glasses, there is no known screen plane. The user can be looking at his nose or at a mountain miles away. The vergence point will indicate how far away the user is looking (the error in distance is small at short distances and increases rapidly as the distance grows. The distance value can't be used at distances greater than a few meters).

Keep in mind that the coordinate system follows the orientation of the head unit. When the head moves, the coordinate system moves with it. The eyes are typically located in approximately the same position relative to the head unit, but gaze angles and gaze points typically change when you turn your head (and keep your eyes fixed on a stationary object).




Figure 23. The glasses measures everything in a 3D coordinate system

In order to get a coordinate on the scene camera image, the 3D gaze point is projected onto the video surface using the camera parameters. This yields a 2D gaze point in relative coordinates to the video. In Pro Lab, this will be exported as MCS in pixels or relative coordinates. There is no information about physical dimensions in the video.

B4 Map data from Glasses recordings onto Snapshots

The coordinate system for recorded data in a Glasses project is based on the scene camera video in the Glasses, as seen in the picture below. When mapping data onto a Snapshot, a new coordinate system is automatically created. The Media Coordinate System is based on the resolution of the Snapshot image, and its origin is located in the upper-left corner.

All gaze data point coordinates in the Data Export output file relating to gaze points are given in pixels. To understand where a gaze point (x_i, y_i) is located on the scene camera video or on the Snapshot, the size in pixels of the scene camera video or the Snapshot need to be taken into consideration.

For detailed information about the scene camera resolution, please refer to the product description for your wearable eye tracker.

The width of a Snapshot can be derived from the Width [Snapshot Name] column in the Data Export file. The height of a Snapshot can be derived from the Height [Snapshot Name] column in the Data Export file.




Figure 24. The coordinate system for recorded data in a Glasses project is based on the scene camera video in the Glasses

Appendix C Gaze Filter functions and effects

During a recording, Tobii Pro eye trackers collect raw eye movement data points every 1.6 to 33 msec. (depending on the sampling data rate of the eye tracker). Each sample will, along with other data, contain a timestamp and gaze coordinates. These coordinates can then be processed further into fixations, which can be overlaid on a video recording of the stimuli used in the study or used to calculate eye tracking metrics. This process is conducted by applying an eye movement classification algorithm, or Gaze Filter, to the data.

C1 Raw Gaze Filter preset type

The Raw data preset type in the Gaze Filter settings does not classify the gaze data into fixations, saccades, or any other eye movement. This preset type only offers the possibility to apply three data processing functions to the data: Gap fill-in (interpolation), Eye selection, and Noise reduction.

Pro Lab includes one built-in preset for the Raw Gaze Filter preset type. The default parameters for this preset are defined below. If you wish to modify the filter settings for your project, select the option “Create new Raw filter” from the Gaze Filter drop-down menu and change the appropriate variables. Your custom settings will be saved as a new preset available in all parts of the application.

The default parameters in the Raw Gaze Filter built in preset are:

- Gap fill-in (interpolation): Off
- Eye selection: Average
- Noise reduction: Off

C2 Tobii I-VT Gaze Filter definition

The general idea behind an I-VT filter (Velocity-Threshold Identification Gaze Filter) is to classify eye movements based on the velocity of the directional shifts of the eye. The velocity is most commonly given in visual degrees per second ($^{\circ}/s$). If the velocity of the eye movement is below a certain threshold, the samples are classified as part of a fixation. If the velocity is above the threshold, it is classified as a saccade.

The Tobii I-VT filter has been developed to work best in stationary situations, and the default values of the I-VT fixation filter parameters have been set to provide accurate fixation classifications for the most common eye tracking use cases in stationary situations. However, depending on the type of eye tracking study conducted, the settings in the I-VT filter can be adjusted to better suit a particular study, to get finer grained fixation data, or to compensate for high levels of noise. The data processing functions within the I-VT fixation filter are described in the sections below.

Pro Lab includes several presets for the Tobii I-VT filter. The default parameters for these are defined below. If you wish to modify the filter settings for your project, select the option “Create new I-VT filter” from the Gaze Filter drop-down menu and change the appropriate variables. Your custom filter settings will be saved as a new filter available in all parts of the application.



The I-VT fixation filter may fail to correctly classify all eye movement recorded with Pro Glasses 2 and Pro Glasses 3 due to various factors. Read more about I-VT gaze filters in the [Learn & Support](#) section of the Tobii Pro website

If accurate fixation and saccade data is required for your study, it is strongly recommended that you always verify the performance of the fixation filter for your data.

For additional information about the fixation filters, please read the following white papers, which can be found on our website:

- Tobii I-VT Fixation Filter – Algorithm Description
- Determining the Tobii I-VT Fixation Filter’s Default Values

C2.1 Gap fill-in (interpolation)

The purpose of the Gap fill-in function is to fill in data where gaze data points are missing. For example, this prevents a fixation in which a few samples are missing as being interpreted as two separate fixations. This loss of gaze data can occur due to temporary interference of eyelashes or other factors. In these cases, the data loss is usually limited to a short period of time, typically less than 50 msec.

Data can also be lost due to other reasons, such as the participant blinking or looking away. This kind of data loss usually results in data gaps longer than 100 msec. Losing data in these latter cases may be perfectly fine for applications of eye tracking where visual attention is measured, since participants do not see anything while blinking or looking away anyway.

In Pro Lab, you can control the limit of how large the gaps (in msec.) in gaze data are that should be filled in by setting the parameter “Max gap length.”

Data is filled in the data gap through linear interpolation. Data points are added along a straight line between neighboring valid data points. Interpolation is done for each eye separately.

C2.2 Eye selection

The eye selection function enables you to choose to discard the gaze data collected from one of the eyes (left or right) or to average the data from both eyes for fixation classification. The only setting available for Glasses projects is Average, which makes an average of the position data from the left and the right eye. If only one eye is detected, it uses the data from that eye.

C2.3 Noise reduction (Moving Median)

All measurement systems, including eye trackers, experience noise. Noise can come from imperfections in the system setup, as well as from influences and interferences around the environment in which the measurement takes place. In eye tracking research where the fixation is the eye movement of interest, other minor eye movements, such as tremor and microsaccades, can also be seen as noise. In the most basic form, I-VT filters calculate the velocity by multiplying the change in position between two consecutive sample points. If the sampling interval is long, the eye will have time to do quite large shifts of direction between two samples, which makes it easy for the eye tracker to distinguish between real eye movements and noise. The higher the sampling frequency, the smaller the eye movement will be between two consecutive samples at a given eye velocity. Noise will therefore have a greater impact in a high frequency system even though it has the same magnitude as in a low frequency system.

The Moving Median noise reduction function is accessible using the Tobii I-VT fixation filter and the Raw data filter.

The median noise reduction function is a symmetric Moving Median filter. It produces output data by calculating the median value of the number of consecutive data points from the input data series. The number of input data points used to produce each output data point is controlled by the Window size parameter. Each produced output data point is given the timestamp of the input data point that was in the center of the window: the median point of the input window (the window size parameter must be an odd number). In cases where the window is not entirely filled with input data points, no output data will be produced. Most often, this happens when the window is short of data points due to the fact that the window stretches outside a data series with valid gaze or when small fractions of data points are missing within the window. Typically, the window stretches outside a valid data series at the beginning and end of a blink or at the beginning and end of a recording.

The loss of small fractions of data can occur as a result of temporary reflections, occlusions, etc. This can partly be solved by using the Gap fill-in function.

C2.4 Noise reduction (Moving Average)

All measurement systems, including eye trackers, experience noise. Noise can come from imperfections in the system setup, as well as from influences and interferences from the environment in which the measurement takes place. In eye tracking research where the fixation is the eye movement of interest, other minor eye movements, such as tremor and microsaccades, can also be seen as noise. In the most basic form, I-VT filters calculate the velocity by multiplying the change in position between two consecutive sample points. If the sampling interval is long, the eye will have time to do quite large shifts of direction between two samples, which makes it easy for the eye tracker to distinguish between real eye movements and noise. The higher the sampling frequency, the smaller the eye movement will be between two consecutive samples at a given eye velocity. Noise will therefore have a greater impact in a high frequency system even though it has the same magnitude as in a low frequency system.

The Moving Average noise reduction function is accessible using the Tobii I-VT fixation filter and the Raw data filter.

The Moving Average noise reduction function is a symmetric Moving Average filter. It produces output data by creating an arithmetic mean of a number of data points from the input data. The number of input data points used to produce each output point is controlled by the Window size parameter. Each produced output data point is given the timestamp of the input data point that was in the center of the window: the input median point (the window size parameter must be an odd number). The window size is dynamically adjusted so that the window never stretches outside the valid data series.

C2.5 Velocity calculator

The Velocity calculator function assigns an angular velocity (visual degrees/second) to each gaze data point. Angular velocity refers to the angular velocity of the eyes relative to the stimuli.

To calculate the gaze velocity on a stimulus (e.g. a stimulus or a Snapshot), it uses three points: the eye position and two gaze points on the stimulus. For a stationary stimulus, such as an image when using a screen-based eye tracker, this is fairly straight forward, but, for data from a wearable eye tracker, it gets a little more complicated. When it comes to calculating the gaze velocity when using a Snapshot and data mapped from a Pro Glasses 2 or Pro Glasses 3 recording, Pro

Lab uses the Gaze Position 3D coordinates. The Velocity calculator function estimates the eyes' angular velocity for each data point by dividing the angular difference between a preceding and a subsequent data point with the time interval between them. The time interval is set by the parameter window length in the Velocity calculator function. The Velocity calculator will only produce velocity output data if the entire window contains input data. This means that gaps in the input data (like the ones caused by blinks) will result in larger gaps in the output data. The size of the output gap will be equal to the input gap plus the number of data points included by the window length parameter.

C2.6 I-VT fixation classifier

The I-VT fixation classifier is based on the I-VT (Velocity-Threshold identification fixation filter), as described by Salvucci and Goldberg (Salvucci & Goldberg, 2000) and Komogortsev et. al. (Komogortsev, Gobert, Jayarathna, Do Hyong Koh, & Gowda, 2010). It is a threshold function that operates on eye tracking data where each data point has an assigned angular velocity. In Pro Lab, angular velocity is assigned to eye tracking data in the Velocity calculator data processing function, as described in the previous section.

The I-VT fixation classifier applies an angular velocity threshold on each data point. The threshold value is given in degrees/second and is adjusted by setting the parameter Velocity threshold in the Gaze Filter settings. Data points with angular velocity below the threshold value are classified as being part of a fixation, and data points above are classified as being part of a saccade.

According to the I-VT fixation classifier, a fixation is an unbroken chain of raw data points all classified as fixation data points. If the velocity cannot be calculated for a raw data point, it is classified as an Unknown Eye Movement. The fixation coordinate is calculated as the arithmetic mean value of the coordinates of the data points belonging to the specific fixation.

C2.7 Merge adjacent fixations

The purpose of the Merge adjacent fixations function is to merge fixations that have been incorrectly classified as multiple short fixations instead of the same, long fixation. Noise and other disturbances can cause the I-VT classifier to incorrectly classify data points that should belong to a fixation as not being part of it. This will split the fixation into multiple fixations, which are all located close together. The Merge adjacent fixations function can be set to merge these multiple fixations into one fixation.

The Merge adjacent fixations function has two threshold parameters:

- Max time between fixations defines the maximum time interval between separate fixations that should be merged.
- Max angle between fixations defines the maximum visual angle of the eyes between separate fixations that should be merged.

C2.8 Discard short fixations

The purpose of the Discard short fixation function is to remove classified fixations that you believe have too short of a duration to be a meaningful fixation. These short fixations may result from head-movements, post-saccadic oscillations, or as inter-saccadic intervals when one saccade execution overrides another. If fixations are used in order to gauge in-depth processing of some visual information, it may make sense to exclude short fixations. The Discard short fixation function can be used to reclassify an incorrectly classified fixation into an Unknown event. That is,

velocity-wise it is not a saccade, but duration-wise it is not a fixation either, consequently placing it in this third category.

The parameter Minimum fixation duration sets the threshold of how short of a duration a fixation can have to be classified as a “fixation.” All events belonging to fixations shorter than the threshold value will be reclassified as Unknown events. The “Reclassify as saccade” functionality will merge two saccades that are interrupted by an Unknown event. This can typically be useful when exporting Glance or saccade metrics where saccades are terminated by such Unknown events when transitioning into fixations. The data samples belonging to both the preceding saccade, the Unknown event, and the succeeding saccade will all merge into a single saccade and its saccade properties will be recalculated.

C2.9 Tobii I-VT Gaze Filter default parameters

Pro Lab includes a few preset settings for the Tobii I-VT Gaze Filter. The default parameters for these filters are defined below. If you wish to modify the filter settings for your project, select the option “Create new I-VT filter” from the Gaze Filter drop-down menu and change the appropriate variables. Your custom filter settings will be saved as a new filter available in all parts of the application.

Tobii I-VT (Fixation)

This filter is suitable for controlled studies where only fixations and saccades are present in the collected data. As with other velocity-based filters, this filter will not classify smooth pursuit and VOR eye movements correctly. To calculate the gaze velocity for data mapped onto Snapshots, this preset uses the Gaze Position 3D coordinates to calculate the true eye velocity.

- Gap fill-in (interpolation)
 - Default: Disabled
- Eye selection
 - Default: Average
- Noise reduction
 - Default: Moving Median
 - Window size: 3 samples
- Velocity calculator
 - Default: window length 20 msec.
- I-VT fixation classifier
 - Default: Threshold 30 degrees/second
- Merge adjacent fixations
 - Default: Enabled
 - Max time between fixations: 75 msec.
 - Max angle between fixations: 0.5 degrees
- Discard short fixations
 - Default: Enabled
 - Minimum fixation duration: 60 msec.

Tobii I-VT (Attention)

The Attention Filter in Pro Lab is essentially the Tobii Pro IV-T Filter, with the velocity threshold parameter set to 100 degrees/second instead of the default 30 degrees/second.

The Attention Filter was created to handle eye tracking data from glasses recordings, which are conducted under dynamic situations, where either the subject is constantly moving or the objects or targets are moving around the subject. In these situations, we use a large array of eye movements to help us keep our fovea aligned with objects and other visual features in the environment - fixations, saccades, smooth pursuits and VOR. (See our article “Types of Eye Movements” for more information.)

We use the filter to separate the moments when we are trying to stabilize our fovea onto something (fixation, smooth pursuit and vestibular ocular reflex), thus potentially extracting information from the location or object, from the moments the eyes are moving too fast to extract information (saccades).

By setting the threshold to 100 degrees/second, we classify fixation, smooth pursuit and most VOR data as “Attention”. However, we also classify 10-15% of saccades as “Attention” as well, which means we overestimate “Attention” slightly.

Setting the IV-T Filter to the default setting of 30 degrees/second will result in underestimation of the periods of “Attention” or information gathering, since quite a large portion of data belonging to smooth pursuits and VOR periods will be classified as saccades.



When using the “Attention” filter, any fixation based AOI metrics will become “foveal S time to first fixation will become “time to the first moment the fovea is stabilized on an AOI”

C3 Effect of mapping data onto Snapshots with Tobii I-VT or Raw Gaze Filter

When mapping data onto a Snapshot with the I-VT fixation filter turned on, the replay window will only display data points that were classified as fixation data. Taking the previous section into consideration, the conclusion is that the replay will only display fixations if both the participant’s head and the object being researched was still (or in a special case, the participant’s head and the wearable eye tracker follow along with a moving object). To be sure to map all available data in all circumstances, make sure you use a Raw Gaze Filter when mapping data onto Snapshots. The Tobii I-VT (Attention) filter preset is suitable for situations where you only want to study attention without separating any eye movements. This preset has the I-VT fixation classifier set to 100 degrees/second and will, in addition to classifying typical fixations, also classify most VOR (Vestibulo-ocular reflex), smooth pursuit eye movements, and slow saccades as fixations.

Always make sure to double check mapped data by replaying the recording with the Snapshot showing mapped data and the Gaze Filter set to Raw.

Appendix D Calculate heat maps

The Heat Map can be based on either fixation data or raw data, depending on whether you are using one of the I-VT Gaze Filters or a Raw data filter. In both cases, data entries consist of a timestamp, duration, and spatial location (X and Y coordinates). There are four different settings for Heat Maps in Pro Lab:

- **Absolute count:** calculated by the number of fixations (or gaze data samples if using raw data)
- **Absolute duration:** calculated by the duration of fixations (or gaze data samples if using raw data)
- **Relative count:** calculated by the number of fixations relative to the total number of fixations made by the participants in the Time of Interest (or gaze data samples instead of fixations if using raw data)
- **Relative duration:** calculated by the duration of the fixations relative to the sum of all fixation durations mapped in the Time of Interest on the Snapshot (or gaze data samples instead of fixations if using raw data)

Eye tracking raw data is composed of fixation and noise data points. The noise includes eye samples during saccades and other small fixation-related movements. A Heat Map based on raw data will, therefore, contain more noise and show larger areas with high color value. When a Heat Map is based on raw data, each entry corresponds to a raw gaze point from the eye tracker, sampled every 1,6 to 33 msec. (depending on the sampling data rate of the eye tracker). Duration values for each entry are constant as they are defined by the duration of the sampling interval.

When the input is based on fixation data, the eye tracker will group the raw data entries into fixations, and the duration of each fixation depends on the Gaze Filter used to identify the fixations. The input data will typically look like this:



Fixation: 1 Start: 132 msec. Duration: 132 msec. Location: (415,711)

Fixation: 2 Start: 264 msec. Duration: 735 msec. Location: (711,491)


Fixation: 3 Start: 999 msec. Duration: 101 msec. Location: (301,10)

D1 Basic concept

The first step in creating a Heat Map is to map the fixations on the stimulus. This is done by going through all the fixations in all selected recordings one by one and adding their values whenever a fixation shares the same X and Y pixel location as another (see figures below). If Count is selected, Pro Lab adds the number of fixations at the same location. If Absolute duration is selected, Pro Lab instead sums the duration of all fixations in the same location. For Relative duration, the duration of each fixation is divided first by the media viewing time and then added.



Once all of the fixation values have been added together, color values are assigned to all the points with the warmest color usually representing the highest value. However, the colors that are used can be customized in the Heat Map settings.




D2 The Heat Map

The basic concept described in the previous chapter would result in a rather useless image, with small colored dots distributed over the picture. Instead, the application sums up the color values from all the points within a certain distance of the fixation location, and the color values gradually decrease as we move away from the fixation point. This process ensures that the color distribution results in a more “smooth” Heat Map image and gives some biological relevance to the color mapping.

The distribution of values around a fixation point is accomplished by using an approximation to the Gaussian curve that is commonly used in 2D image processing, a cubic Hermite spline polynomial (cspline). The specific function of the polynomial used is $t^2(3-2t)$ - represented by the blue line in the graph below. (The orange line represents the fixation point and maximum color value.)


Heatmap: Colour distribution, Cubic Hermite Spline polynomial



The radius of the function can be adjusted in the Visualization Type and Settings tool, where the option “Heat Map” is selected in a Visualization tab. The default value is set to 50 pixels, corresponding to a total kernel of 100 pixels.

Since the kernel is set in pixels, our Heat Maps are dependent on screen or scene camera video resolution and less representative of the individual’s foveal vision, if the participant varies his distance greatly to the stimulus.

The illustration below shows how the Heat Map calculations using fixation count are done to assign colors to the fixations: For Absolute and Relative duration the calculations are the same. However, instead of using the number of fixations as values, we use the fixation duration or fixation duration/media viewing time, respectively.



Appendix E Pro Lab project and media limitations

Media Type	Maximum resolution recommended system	Maximum resolution minimum system
Image Stimuli	2560 x 1440	1280 x 720
Snapshot Images	50 MP	10 MP
Video Stimuli	1920 x 1080 @ 30 fps	1280 x 720 @ 30 fps
Web and Screen Stimuli	1920 x 1080 @ 30 fps	

Project type	Number of Recordings with recommended system	Number of Recordings with minimum system
Glasses	500 of 20 minutes each	50 of 20 minutes each
Screen based Project with Images	200 with 400 Images	30 with 100 Images
Screen based Project with Videos	200 with 1 x 90 min video or 50 x 2 min videos	30 with 1 x 60 min video or 20 x 2 min videos

Supported file formats	Format	Video codecs
Image Stimuli and Snapshots	jpeg, png, bmp	-
Video	mp4 (recommended), avi	H264 (recommended), DIVX, XVID
Audio	mp3, aac, ac3, pcm	-

Appendix F Analysis of an External Presenter project recorded by E-Prime

Prerequisites:

- If you want to log manual responses, make sure that the relevant E-Object in the experiment is defined in the TPLSetDisplayEventStimulus routine. If this E-Object is defined, you get RESP and ACC values from E-Prime for it as events in Tobii Pro Lab.
- All relevant stimulus variables are encoded in the media name that is sent to Tobii Pro Lab.
- In the example, we assume that you analyze the sample experiment called "TPLFixedPositionAOI" found in your My Experiments folder after installing E-Prime Extensions for Tobii Pro Lab.

What to do:

1. Create a *Time of Interest* (TOI) definition for your trials. In this example, participants have responded manually to the stimulus, either correctly or incorrectly. This accuracy value (0 or 1) has been sent to Tobii Pro Lab and appears in the timeline if you have checked Imported Events (enabled by default). Because of the presence of manual responses, we need to define two Times of Interest: one for correct responses (ACC 1) and one for incorrect responses (ACC 0).
 - a. Create a new custom TOI by clicking on the **Plus (+)** sign on the Time of Interest panel.
 - b. Give the first TOI a name, for example *Correct responses*.
 - c. For its *Start* point, go to *Remote Event* types and pick the event "Stimulus_ACC_1" (it marks all correct responses of the stimulus object in E-Prime).
 - d. For its *End* point, go to *Remote Event* types and select the event "Stimulus_End", which marks the end of all of the objects called "Stimulus" in the E-Prime experiment.
 - e. Create a similar custom TOI for all "Incorrect responses".
2. If you want visualizations of the gaze on your stimuli, go to the *Visualizations* tab and decide what type of visualization you want.
 - a. Select the graphical nature of the visualization: *Gaze Replay* or *Heat Map*, and their parameters.
 - b. Determine what data you want to visualize.
 - You can easily visualize the data from one particular stimulus in the list of Media TOIs in the *Time of Interest* panel to the right.
 - If you want to visualize all data of a similar type, click on the custom TOIs you created previously ("Correct responses" and "Incorrect responses"). That way, you visualize the gaze for those intervals of data.

- However, the background is empty (a gray checkerboard pattern) because there is no stimulus that can be exclusively mapped to a custom TOI (you have probably used several stimuli). If you want a background for a custom TOI, go to the replay of the recording, locate a part of the recording that displays the stimulus you want as a background, and do a *Create TOI with frame* by clicking on the plus sign near the time counter.
 - It is better that you only create this TOI for visualization, because any metrics you export for this TOI will be calculated from AOIs drawn on this separate media frame, and not from the stimuli and AOIs generated from the external presenter client, E-Prime.
3. When you are ready to export metrics for further analysis, select *Analysis - Metrics Export*.
 - a. Select the desired metrics.
 - b. Go to the *Data selection* panel to the right and select all desired recordings.
 - c. Select the desired TOIs. In the example they are "Correct responses" and "Incorrect responses", both mentioned above.
 - d. Finally, select the desired areas of interest. However they are unique per stimulus so they will create many columns in your export. Rather, use AOI tags to group many unique AOIs into a few categories.
 - e. In this experiment, only the ungrouped tags "Non-target" and "Target" are relevant. Deselect all other AOIs, click *Export* and save your file.
 4. Now you have exported your data in a text file with rows and columns. Open the file in Excel or a similar spreadsheet program. (You may need to do a text-to-column operation to properly delimit the file on the tab character.) Every row is a trial (an interval of a TOI), and every column is a property of the trial or the experiment. Before running statistical analysis, you need to recover your experimental variables, which are encoded in the media name of the stimuli, looking like "cow-horse-Target=horse".
 - a. Move this *Media* column to the end after all the other columns.
 - b. Select it and perform a text-to-column operation with "-" as the delimiter, and you will now have three new columns. The left column will be the image to the left, and the center one the image to the right.
 - c. Go ahead and name them as "LeftImage" and "RightImage" now. The final column has content in the form of "Target=horse".
 - d. Do a text-to-column operation with "=" as the delimiter, and delete the first column from this operation.
 - e. You now have the column with just the target name. Name this column "TargetImage". Now you are ready.
 5. For your next experiment, you can use the same delimiter so you have to do only one text-to-column operation. Note that you can use the TOI column, due to the names you gave the TOIs, to separate correct responses from incorrect responses in your analysis.

Appendix G Support, Training, and Warranty

G1 Customer Support

If you need help, please contact [Customer Support](#) at Tobii Pro. In order to receive assistance as quickly as possible, make sure you have access to your Tobii Pro device and, if possible, an Internet connection. You should also be able to supply the serial number of the device, which you will find on a sticker on the back or bottom of the device.

G1.1 Get help online

Many questions can be answered by visiting Tobii Pro Connect. It contains the latest information about contacting Support, helpful articles and FAQs, links to downloads, and much more. Log in or register to see information about your account and to reach Customer Support on [Tobii Pro Connect](#).

G2 Training and Education Services

If you are new to eye tracking, or want to extend your knowledge about eye tracking research, sign up for one of our learning programs and events on our website at [Training and Education Services](#).

G3 Warranty information

Read more online about [Tobii Pro Care and Tobii Pro's eye tracker warranty](#).



Copyright ©Tobii AB (publ). Not all products and services offered in each local market. Specifications are subject to change without prior notice. All trademarks are the property of their respective owners.

Support for your Tobii Pro product

Get Help Online

Visit Tobii Pro Connect for help with your Tobii Pro device. It contains the latest information about contacting Support, links to our Learning Center, and much more. Visit connect.tobii.com.

Contact Your Solution Consultant or Reseller

For questions or problems with your product, contact your Tobii Pro sales representative or authorized reseller for assistance. They are most familiar with your personal setup and can best help you with tips and product training.

Visit tobii.com/contact.